

IP STORAGE TEMPERATURE LOG

Period of Record	Year. Month	2026.6	
Storage	Location	동국대학교일산병원 임상시험센터 약국	9134Room
	Description		실온보관장
	Manufacturer	NA	
	Model Name	NA	
	Serial Number	NA	
	Temperature Range	1~30℃	
	Setting Point	20~23℃	
Thermometer (including alarm system)	Description	Wi-Fi data logger	
	Manufacturer	Testo Korea	
	Model Name	Saveris2-T3	
	Serial Number	0054704396	
	Sensor Channel	Channel 2	
Thermometer Calibration Certification	Min/Max Record	Monthly report 2page / This month Min (24.0)℃, Max(28.6)℃	
	Date of Calibration	2025.11.25	
	Date of Calibration Due	2026.11.25	
Alarm System	Certificate No.	A25-2356-001	
	Alarm threshold	Lower : 10.0℃, Upper : 28.5℃, immediate alarm without delay	
	Alarm Method	Mail, SMS to Pharmacists	

Date	Time	Temperature(℃)	Deviation	Verified by	Comment
2026.06.01	P:30	24.8	☐ Y / ☒ N	ml	
2026.06.02	P:30	24.8	☐ Y / ☒ N	ml	49기 IP 반영
2026.06.03			☐ Y / ☐ N		2026 지방선거
2026.06.04	10:00	26.4	☐ Y / ☒ N	ml	
2026.06.05	P:30	25.6	☐ Y / ☒ N	ml	HL 1113 IP 01.2
2026.06.06			☐ Y / ☐ N		
2026.06.07			☐ Y / ☐ N		
2026.06.08	10:00	25.6	☐ Y / ☒ N	ml	
2026.06.09	P:10	25.8	☐ Y / ☒ N	ml	
2026.06.10	P:20	26.1	☐ Y / ☒ N	ml	
2026.06.11	P:40	25.7	☐ Y / ☒ N	ml	
2026.06.12	P:40	25.5	☐ Y / ☒ N	ml	
2026.06.13			☐ Y / ☐ N		
2026.06.14			☐ Y / ☐ N		
2026.06.15	P:30	26.1	☐ Y / ☒ N	ml	
2026.06.16	P:35	26.8	☐ Y / ☒ N	ml	
2026.06.17	P:40	24.1	☐ Y / ☒ N	ml	
2026.06.18	P:20	24.5	☐ Y / ☒ N	ml	
2026.06.19	P:30	24.2	☐ Y / ☒ N	ml	
2026.06.20			☐ Y / ☐ N		
2026.06.21			☐ Y / ☐ N		
2026.06.22	P:30	24.2	☐ Y / ☒ N	ml	
2026.06.23	P:20	23.7	☐ Y / ☒ N	ml	
2026.06.24	P:15	24.1	☐ Y / ☒ N	ml	
2026.06.25	P:20	23.3	☐ Y / ☒ N	ml	
2026.06.26	P:10	24.1	☐ Y / ☒ N	ml	
2026.06.27			☐ Y / ☐ N		
2026.06.28			☐ Y / ☐ N		
2026.06.29	P:30	24.3	☐ Y / ☒ N	765ml	
2026.06.30	P:20	24.4	☐ Y / ☒ N	ml	
			☐ Y / ☐ N		

date of copying

sign



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 테이블

측정 포인트 요약

동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_CH2

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
Min	23.0
Max	28.6
Average	25.1
26. 6. 1. 오전 12:00:00	24.9
26. 6. 1. 오전 12:30:00	24.9
26. 6. 1. 오전 1:00:00	24.8
26. 6. 1. 오전 1:30:00	24.9
26. 6. 1. 오전 2:00:00	24.9
26. 6. 1. 오전 2:30:00	24.8
26. 6. 1. 오전 3:00:00	24.8
26. 6. 1. 오전 3:30:00	24.7
26. 6. 1. 오전 4:00:00	24.7
26. 6. 1. 오전 4:30:00	24.7
26. 6. 1. 오전 5:00:00	24.6
26. 6. 1. 오전 5:30:00	24.6
26. 6. 1. 오전 6:00:00	24.6
26. 6. 1. 오전 6:30:00	24.6
26. 6. 1. 오전 7:00:00	24.6
26. 6. 1. 오전 7:30:00	24.6
26. 6. 1. 오전 8:00:00	24.7
26. 6. 1. 오전 8:30:00	24.7
26. 6. 1. 오전 9:00:00	24.8
26. 6. 1. 오전 9:30:00	24.8
26. 6. 1. 오전 10:00:00	24.9



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃) CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 1. 오전 10:30:00	25.1
26. 6. 1. 오전 11:00:00	25.2
26. 6. 1. 오전 11:30:00	25.4
26. 6. 1. 오후 12:00:00	25.4
26. 6. 1. 오후 12:30:00	25.5
26. 6. 1. 오후 1:00:00	25.7
26. 6. 1. 오후 1:30:00	25.7
26. 6. 1. 오후 2:00:00	25.7
26. 6. 1. 오후 2:30:00	25.7
26. 6. 1. 오후 3:00:00	25.7
26. 6. 1. 오후 3:30:00	25.7
26. 6. 1. 오후 4:00:00	25.6
26. 6. 1. 오후 4:30:00	25.6
26. 6. 1. 오후 5:00:00	25.6
26. 6. 1. 오후 5:30:00	25.5
26. 6. 1. 오후 6:00:00	25.5
26. 6. 1. 오후 6:30:00	25.3
26. 6. 1. 오후 7:00:00	25.3
26. 6. 1. 오후 7:30:00	25.2
26. 6. 1. 오후 8:00:00	25.1
26. 6. 1. 오후 8:30:00	25.0
26. 6. 1. 오후 9:00:00	25.0
26. 6. 1. 오후 9:30:00	24.9
26. 6. 1. 오후 10:00:00	24.9
26. 6. 1. 오후 10:30:00	24.9
26. 6. 1. 오후 11:00:00	24.8
26. 6. 1. 오후 11:30:00	24.8
26. 6. 2. 오전 12:00:00	24.8
26. 6. 2. 오전 12:30:00	24.8
26. 6. 2. 오전 1:00:00	24.8
26. 6. 2. 오전 1:30:00	24.7
26. 6. 2. 오전 2:00:00	24.7



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 2. 오전 2:30:00	24.7
26. 6. 2. 오전 3:00:00	24.6
26. 6. 2. 오전 3:30:00	24.7
26. 6. 2. 오전 4:00:00	24.6
26. 6. 2. 오전 4:30:00	24.6
26. 6. 2. 오전 5:00:00	24.6
26. 6. 2. 오전 5:30:00	24.5
26. 6. 2. 오전 6:00:00	24.5
26. 6. 2. 오전 6:30:00	24.5
26. 6. 2. 오전 7:00:00	24.5
26. 6. 2. 오전 7:30:00	24.6
26. 6. 2. 오전 8:00:00	24.6
26. 6. 2. 오전 8:30:00	24.7
26. 6. 2. 오전 9:00:00	24.7
26. 6. 2. 오전 9:30:00	24.8
26. 6. 2. 오전 10:00:00	24.8
26. 6. 2. 오전 10:30:00	25.0
26. 6. 2. 오전 11:00:00	25.2
26. 6. 2. 오전 11:30:00	25.4
26. 6. 2. 오후 12:00:00	25.6
26. 6. 2. 오후 12:30:00	25.8
26. 6. 2. 오후 1:00:00	26.0
26. 6. 2. 오후 1:30:00	26.1
26. 6. 2. 오후 2:00:00	26.2
26. 6. 2. 오후 2:30:00	26.3
26. 6. 2. 오후 3:00:00	26.4
26. 6. 2. 오후 3:30:00	26.4
26. 6. 2. 오후 4:00:00	26.4
26. 6. 2. 오후 4:30:00	26.4
26. 6. 2. 오후 5:00:00	26.5
26. 6. 2. 오후 5:30:00	26.5
26. 6. 2. 오후 6:00:00	26.5



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 2. 오후 6:30:00	26.3
26. 6. 2. 오후 7:00:00	26.2
26. 6. 2. 오후 7:30:00	26.1
26. 6. 2. 오후 8:00:00	26.1
26. 6. 2. 오후 8:30:00	26.0
26. 6. 2. 오후 9:00:00	25.9
26. 6. 2. 오후 9:30:00	25.8
26. 6. 2. 오후 10:00:00	25.8
26. 6. 2. 오후 10:30:00	25.8
26. 6. 2. 오후 11:00:00	25.8
26. 6. 2. 오후 11:30:00	25.7
26. 6. 3. 오전 12:00:00	25.7
26. 6. 3. 오전 12:30:00	26.1
26. 6. 3. 오전 1:00:00	25.9
26. 6. 3. 오전 1:30:00	25.7
26. 6. 3. 오전 2:00:00	25.6
26. 6. 3. 오전 2:30:00	25.6
26. 6. 3. 오전 3:00:00	25.5
26. 6. 3. 오전 3:30:00	25.4
26. 6. 3. 오전 4:00:00	25.4
26. 6. 3. 오전 4:30:00	25.3
26. 6. 3. 오전 5:00:00	25.2
26. 6. 3. 오전 5:30:00	25.2
26. 6. 3. 오전 6:00:00	25.1
26. 6. 3. 오전 6:30:00	25.1
26. 6. 3. 오전 7:00:00	25.1
26. 6. 3. 오전 7:30:00	25.1
26. 6. 3. 오전 8:00:00	25.2
26. 6. 3. 오전 8:30:00	25.1
26. 6. 3. 오전 9:00:00	25.2
26. 6. 3. 오전 9:30:00	25.2
26. 6. 3. 오전 10:00:00	25.3



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 3. 오전 10:30:00	25.3
26. 6. 3. 오전 11:00:00	25.4
26. 6. 3. 오전 11:30:00	25.6
26. 6. 3. 오후 12:00:00	25.8
26. 6. 3. 오후 12:30:00	25.9
26. 6. 3. 오후 1:00:00	26.1
26. 6. 3. 오후 1:30:00	26.2
26. 6. 3. 오후 2:00:00	26.3
26. 6. 3. 오후 2:30:00	26.3
26. 6. 3. 오후 3:00:00	26.2
26. 6. 3. 오후 3:30:00	26.2
26. 6. 3. 오후 4:00:00	26.3
26. 6. 3. 오후 4:30:00	26.2
26. 6. 3. 오후 5:00:00	26.2
26. 6. 3. 오후 5:30:00	26.2
26. 6. 3. 오후 6:00:00	26.2
26. 6. 3. 오후 6:30:00	26.1
26. 6. 3. 오후 7:00:00	26.0
26. 6. 3. 오후 7:30:00	26.0
26. 6. 3. 오후 8:00:00	25.9
26. 6. 3. 오후 8:30:00	25.8
26. 6. 3. 오후 9:00:00	25.8
26. 6. 3. 오후 9:30:00	25.7
26. 6. 3. 오후 10:00:00	25.6
26. 6. 3. 오후 10:30:00	25.6
26. 6. 3. 오후 11:00:00	25.5
26. 6. 3. 오후 11:30:00	25.5
26. 6. 4. 오전 12:00:00	25.7
26. 6. 4. 오전 12:30:00	26.2
26. 6. 4. 오전 1:00:00	26.6
26. 6. 4. 오전 1:30:00	27.0
26. 6. 4. 오전 2:00:00	27.2



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 4. 오전 2:30:00	27.3
26. 6. 4. 오전 3:00:00	27.5
26. 6. 4. 오전 3:30:00	27.7
26. 6. 4. 오전 4:00:00	27.8
26. 6. 4. 오전 4:30:00	27.9
26. 6. 4. 오전 5:00:00	28.0
26. 6. 4. 오전 5:30:00	28.1
26. 6. 4. 오전 6:00:00	27.9
26. 6. 4. 오전 6:30:00	27.3
26. 6. 4. 오전 7:00:00	27.0
26. 6. 4. 오전 7:30:00	26.7
26. 6. 4. 오전 8:00:00	26.6
26. 6. 4. 오전 8:30:00	26.4
26. 6. 4. 오전 9:00:00	26.3
26. 6. 4. 오전 9:30:00	26.4
26. 6. 4. 오전 10:00:00	26.5
26. 6. 4. 오전 10:30:00	26.5
26. 6. 4. 오전 11:00:00	26.7
26. 6. 4. 오전 11:30:00	26.8
26. 6. 4. 오후 12:00:00	26.8
26. 6. 4. 오후 12:30:00	26.6
26. 6. 4. 오후 1:00:00	26.5
26. 6. 4. 오후 1:30:00	26.4
26. 6. 4. 오후 2:00:00	26.3
26. 6. 4. 오후 2:30:00	26.3
26. 6. 4. 오후 3:00:00	26.3
26. 6. 4. 오후 3:30:00	26.3
26. 6. 4. 오후 4:00:00	26.3
26. 6. 4. 오후 4:30:00	26.1
26. 6. 4. 오후 5:00:00	26.0
26. 6. 4. 오후 5:30:00	26.0
26. 6. 4. 오후 6:00:00	26.0



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 4. 오후 6:30:00	25.9
26. 6. 4. 오후 7:00:00	25.8
26. 6. 4. 오후 7:30:00	25.7
26. 6. 4. 오후 8:00:00	25.7
26. 6. 4. 오후 8:30:00	25.6
26. 6. 4. 오후 9:00:00	25.5
26. 6. 4. 오후 9:30:00	25.5
26. 6. 4. 오후 10:00:00	25.5
26. 6. 4. 오후 10:30:00	25.4
26. 6. 4. 오후 11:00:00	25.3
26. 6. 4. 오후 11:30:00	25.3
26. 6. 5. 오전 12:00:00	25.3
26. 6. 5. 오전 12:30:00	25.5
26. 6. 5. 오전 1:00:00	25.9
26. 6. 5. 오전 1:30:00	26.2
26. 6. 5. 오전 2:00:00	26.5
26. 6. 5. 오전 2:30:00	26.7
26. 6. 5. 오전 3:00:00	26.9
26. 6. 5. 오전 3:30:00	27.0
26. 6. 5. 오전 4:00:00	27.1
26. 6. 5. 오전 4:30:00	27.1
26. 6. 5. 오전 5:00:00	27.2
26. 6. 5. 오전 5:30:00	27.3
26. 6. 5. 오전 6:00:00	27.3
26. 6. 5. 오전 6:30:00	26.8
26. 6. 5. 오전 7:00:00	26.4
26. 6. 5. 오전 7:30:00	26.1
26. 6. 5. 오전 8:00:00	25.9
26. 6. 5. 오전 8:30:00	25.8
26. 6. 5. 오전 9:00:00	25.6
26. 6. 5. 오전 9:30:00	25.5
26. 6. 5. 오전 10:00:00	25.5



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 5. 오전 10:30:00	25.5
26. 6. 5. 오전 11:00:00	25.5
26. 6. 5. 오전 11:30:00	25.4
26. 6. 5. 오후 12:00:00	25.5
26. 6. 5. 오후 12:30:00	25.6
26. 6. 5. 오후 1:00:00	25.6
26. 6. 5. 오후 1:30:00	25.6
26. 6. 5. 오후 2:00:00	25.7
26. 6. 5. 오후 2:30:00	25.7
26. 6. 5. 오후 3:00:00	25.7
26. 6. 5. 오후 3:30:00	25.8
26. 6. 5. 오후 4:00:00	25.7
26. 6. 5. 오후 4:30:00	25.6
26. 6. 5. 오후 5:00:00	25.6
26. 6. 5. 오후 5:30:00	25.6
26. 6. 5. 오후 6:00:00	25.5
26. 6. 5. 오후 6:30:00	25.5
26. 6. 5. 오후 7:00:00	25.4
26. 6. 5. 오후 7:30:00	25.3
26. 6. 5. 오후 8:00:00	25.3
26. 6. 5. 오후 8:30:00	25.2
26. 6. 5. 오후 9:00:00	25.0
26. 6. 5. 오후 9:30:00	25.0
26. 6. 5. 오후 10:00:00	25.0
26. 6. 5. 오후 10:30:00	24.9
26. 6. 5. 오후 11:00:00	24.9
26. 6. 5. 오후 11:30:00	24.9
26. 6. 6. 오전 12:00:00	24.9
26. 6. 6. 오전 12:30:00	25.0
26. 6. 6. 오전 1:00:00	25.5
26. 6. 6. 오전 1:30:00	25.9
26. 6. 6. 오전 2:00:00	26.2



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 6. 오전 2:30:00	26.4
26. 6. 6. 오전 3:00:00	26.5
26. 6. 6. 오전 3:30:00	26.6
26. 6. 6. 오전 4:00:00	26.7
26. 6. 6. 오전 4:30:00	26.8
26. 6. 6. 오전 5:00:00	26.8
26. 6. 6. 오전 5:30:00	26.9
26. 6. 6. 오전 6:00:00	26.9
26. 6. 6. 오전 6:30:00	26.6
26. 6. 6. 오전 7:00:00	26.1
26. 6. 6. 오전 7:30:00	25.7
26. 6. 6. 오전 8:00:00	25.6
26. 6. 6. 오전 8:30:00	25.5
26. 6. 6. 오전 9:00:00	25.5
26. 6. 6. 오전 9:30:00	25.5
26. 6. 6. 오전 10:00:00	25.5
26. 6. 6. 오전 10:30:00	25.4
26. 6. 6. 오전 11:00:00	25.4
26. 6. 6. 오전 11:30:00	25.5
26. 6. 6. 오후 12:00:00	25.6
26. 6. 6. 오후 12:30:00	25.8
26. 6. 6. 오후 1:00:00	25.8
26. 6. 6. 오후 1:30:00	25.9
26. 6. 6. 오후 2:00:00	26.0
26. 6. 6. 오후 2:30:00	26.0
26. 6. 6. 오후 3:00:00	26.1
26. 6. 6. 오후 3:30:00	26.0
26. 6. 6. 오후 4:00:00	26.1
26. 6. 6. 오후 4:30:00	26.1
26. 6. 6. 오후 5:00:00	26.1
26. 6. 6. 오후 5:30:00	26.1
26. 6. 6. 오후 6:00:00	26.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 6. 오후 6:30:00	26.0
26. 6. 6. 오후 7:00:00	25.9
26. 6. 6. 오후 7:30:00	25.8
26. 6. 6. 오후 8:00:00	25.8
26. 6. 6. 오후 8:30:00	25.8
26. 6. 6. 오후 9:00:00	25.7
26. 6. 6. 오후 9:30:00	25.7
26. 6. 6. 오후 10:00:00	25.6
26. 6. 6. 오후 10:30:00	25.6
26. 6. 6. 오후 11:00:00	25.5
26. 6. 6. 오후 11:30:00	25.4
26. 6. 7. 오전 12:00:00	25.9
26. 6. 7. 오전 12:30:00	26.4
26. 6. 7. 오전 1:00:00	26.8
26. 6. 7. 오전 1:30:00	27.1
26. 6. 7. 오전 2:00:00	27.3
26. 6. 7. 오전 2:30:00	27.5
26. 6. 7. 오전 3:00:00	27.6
26. 6. 7. 오전 3:30:00	27.7
26. 6. 7. 오전 4:00:00	27.8
26. 6. 7. 오전 4:30:00	28.0
26. 6. 7. 오전 5:00:00	28.0
26. 6. 7. 오전 5:30:00	28.1
26. 6. 7. 오전 6:00:00	28.1
26. 6. 7. 오전 6:30:00	27.7
26. 6. 7. 오전 7:00:00	27.2
26. 6. 7. 오전 7:30:00	26.8
26. 6. 7. 오전 8:00:00	26.5
26. 6. 7. 오전 8:30:00	26.3
26. 6. 7. 오전 9:00:00	26.2
26. 6. 7. 오전 9:30:00	26.2
26. 6. 7. 오전 10:00:00	26.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃) CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 7. 오전 10:30:00	26.0
26. 6. 7. 오전 11:00:00	26.0
26. 6. 7. 오전 11:30:00	26.0
26. 6. 7. 오후 12:00:00	26.0
26. 6. 7. 오후 12:30:00	25.9
26. 6. 7. 오후 1:00:00	25.9
26. 6. 7. 오후 1:30:00	25.8
26. 6. 7. 오후 2:00:00	25.8
26. 6. 7. 오후 2:30:00	25.8
26. 6. 7. 오후 3:00:00	25.7
26. 6. 7. 오후 3:30:00	25.7
26. 6. 7. 오후 4:00:00	25.7
26. 6. 7. 오후 4:30:00	25.6
26. 6. 7. 오후 5:00:00	25.6
26. 6. 7. 오후 5:30:00	25.6
26. 6. 7. 오후 6:00:00	25.6
26. 6. 7. 오후 6:30:00	25.6
26. 6. 7. 오후 7:00:00	25.5
26. 6. 7. 오후 7:30:00	25.4
26. 6. 7. 오후 8:00:00	25.4
26. 6. 7. 오후 8:30:00	25.3
26. 6. 7. 오후 9:00:00	25.2
26. 6. 7. 오후 9:30:00	25.1
26. 6. 7. 오후 10:00:00	25.0
26. 6. 7. 오후 10:30:00	25.0
26. 6. 7. 오후 11:00:00	25.0
26. 6. 7. 오후 11:30:00	24.9
26. 6. 8. 오전 12:00:00	25.3
26. 6. 8. 오전 12:30:00	25.8
26. 6. 8. 오전 1:00:00	26.2
26. 6. 8. 오전 1:30:00	26.5
26. 6. 8. 오전 2:00:00	26.7



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 8. 오전 2:30:00	26.8
26. 6. 8. 오전 3:00:00	27.0
26. 6. 8. 오전 3:30:00	27.1
26. 6. 8. 오전 4:00:00	27.2
26. 6. 8. 오전 4:30:00	27.2
26. 6. 8. 오전 5:00:00	27.3
26. 6. 8. 오전 5:30:00	27.4
26. 6. 8. 오전 6:00:00	27.5
26. 6. 8. 오전 6:30:00	27.0
26. 6. 8. 오전 7:00:00	26.5
26. 6. 8. 오전 7:30:00	26.2
26. 6. 8. 오전 8:00:00	25.9
26. 6. 8. 오전 8:30:00	25.8
26. 6. 8. 오전 9:00:00	25.7
26. 6. 8. 오전 9:30:00	25.6
26. 6. 8. 오전 10:00:00	25.7
26. 6. 8. 오전 10:30:00	25.7
26. 6. 8. 오전 11:00:00	25.8
26. 6. 8. 오전 11:30:00	25.9
26. 6. 8. 오후 12:00:00	26.0
26. 6. 8. 오후 12:30:00	26.0
26. 6. 8. 오후 1:00:00	26.1
26. 6. 8. 오후 1:30:00	26.0
26. 6. 8. 오후 2:00:00	26.0
26. 6. 8. 오후 2:30:00	26.0
26. 6. 8. 오후 3:00:00	26.1
26. 6. 8. 오후 3:30:00	26.1
26. 6. 8. 오후 4:00:00	25.9
26. 6. 8. 오후 4:30:00	25.9
26. 6. 8. 오후 5:00:00	25.7
26. 6. 8. 오후 5:30:00	25.7
26. 6. 8. 오후 6:00:00	25.7



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃) CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 8. 오후 6:30:00	25.7
26. 6. 8. 오후 7:00:00	25.8
26. 6. 8. 오후 7:30:00	25.9
26. 6. 8. 오후 8:00:00	25.8
26. 6. 8. 오후 8:30:00	25.7
26. 6. 8. 오후 9:00:00	25.7
26. 6. 8. 오후 9:30:00	25.7
26. 6. 8. 오후 10:00:00	25.7
26. 6. 8. 오후 10:30:00	25.6
26. 6. 8. 오후 11:00:00	25.5
26. 6. 8. 오후 11:30:00	25.5
26. 6. 9. 오전 12:00:00	25.4
26. 6. 9. 오전 12:30:00	25.7
26. 6. 9. 오전 1:00:00	26.2
26. 6. 9. 오전 1:30:00	26.6
26. 6. 9. 오전 2:00:00	26.8
26. 6. 9. 오전 2:30:00	27.0
26. 6. 9. 오전 3:00:00	27.1
26. 6. 9. 오전 3:30:00	27.2
26. 6. 9. 오전 4:00:00	27.3
26. 6. 9. 오전 4:30:00	27.4
26. 6. 9. 오전 5:00:00	27.4
26. 6. 9. 오전 5:30:00	27.5
26. 6. 9. 오전 6:00:00	27.5
26. 6. 9. 오전 6:30:00	27.0
26. 6. 9. 오전 7:00:00	26.5
26. 6. 9. 오전 7:30:00	26.2
26. 6. 9. 오전 8:00:00	26.0
26. 6. 9. 오전 8:30:00	25.8
26. 6. 9. 오전 9:00:00	25.8
26. 6. 9. 오전 9:30:00	25.7
26. 6. 9. 오전 10:00:00	25.7



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 9. 오전 10:30:00	25.6
26. 6. 9. 오전 11:00:00	25.7
26. 6. 9. 오전 11:30:00	25.7
26. 6. 9. 오후 12:00:00	25.8
26. 6. 9. 오후 12:30:00	25.9
26. 6. 9. 오후 1:00:00	25.9
26. 6. 9. 오후 1:30:00	25.9
26. 6. 9. 오후 2:00:00	25.9
26. 6. 9. 오후 2:30:00	25.9
26. 6. 9. 오후 3:00:00	25.9
26. 6. 9. 오후 3:30:00	25.9
26. 6. 9. 오후 4:00:00	25.9
26. 6. 9. 오후 4:30:00	25.9
26. 6. 9. 오후 5:00:00	25.9
26. 6. 9. 오후 5:30:00	25.8
26. 6. 9. 오후 6:00:00	25.7
26. 6. 9. 오후 6:30:00	25.7
26. 6. 9. 오후 7:00:00	25.6
26. 6. 9. 오후 7:30:00	25.6
26. 6. 9. 오후 8:00:00	25.5
26. 6. 9. 오후 8:30:00	25.5
26. 6. 9. 오후 9:00:00	25.4
26. 6. 9. 오후 9:30:00	25.4
26. 6. 9. 오후 10:00:00	25.3
26. 6. 9. 오후 10:30:00	25.3
26. 6. 9. 오후 11:00:00	25.2
26. 6. 9. 오후 11:30:00	25.2
26. 6. 10. 오전 12:00:00	25.4
26. 6. 10. 오전 12:30:00	26.0
26. 6. 10. 오전 1:00:00	26.3
26. 6. 10. 오전 1:30:00	26.6
26. 6. 10. 오전 2:00:00	26.8



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃) CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 10. 오전 2:30:00	27.0
26. 6. 10. 오전 3:00:00	27.1
26. 6. 10. 오전 3:30:00	27.2
26. 6. 10. 오전 4:00:00	27.3
26. 6. 10. 오전 4:30:00	27.4
26. 6. 10. 오전 5:00:00	27.5
26. 6. 10. 오전 5:30:00	27.5
26. 6. 10. 오전 6:00:00	27.6
26. 6. 10. 오전 6:30:00	27.1
26. 6. 10. 오전 7:00:00	26.6
26. 6. 10. 오전 7:30:00	26.4
26. 6. 10. 오전 8:00:00	26.2
26. 6. 10. 오전 8:30:00	26.2
26. 6. 10. 오전 9:00:00	26.1
26. 6. 10. 오전 9:30:00	26.0
26. 6. 10. 오전 10:00:00	26.1
26. 6. 10. 오전 10:30:00	26.1
26. 6. 10. 오전 11:00:00	26.2
26. 6. 10. 오전 11:30:00	26.1
26. 6. 10. 오후 12:00:00	26.2
26. 6. 10. 오후 12:30:00	26.3
26. 6. 10. 오후 1:00:00	26.3
26. 6. 10. 오후 1:30:00	26.4
26. 6. 10. 오후 2:00:00	26.4
26. 6. 10. 오후 2:30:00	26.4
26. 6. 10. 오후 3:00:00	26.3
26. 6. 10. 오후 3:30:00	26.4
26. 6. 10. 오후 4:00:00	26.3
26. 6. 10. 오후 4:30:00	26.2
26. 6. 10. 오후 5:00:00	26.1
26. 6. 10. 오후 5:30:00	26.1
26. 6. 10. 오후 6:00:00	26.0



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 10. 오후 6:30:00	26.0
26. 6. 10. 오후 7:00:00	25.9
26. 6. 10. 오후 7:30:00	25.7
26. 6. 10. 오후 8:00:00	25.6
26. 6. 10. 오후 8:30:00	25.5
26. 6. 10. 오후 9:00:00	25.5
26. 6. 10. 오후 9:30:00	25.4
26. 6. 10. 오후 10:00:00	25.3
26. 6. 10. 오후 10:30:00	25.2
26. 6. 10. 오후 11:00:00	25.1
26. 6. 10. 오후 11:30:00	25.1
26. 6. 11. 오전 12:00:00	25.0
26. 6. 11. 오전 12:30:00	25.1
26. 6. 11. 오전 1:00:00	25.6
26. 6. 11. 오전 1:30:00	25.9
26. 6. 11. 오전 2:00:00	26.2
26. 6. 11. 오전 2:30:00	26.4
26. 6. 11. 오전 3:00:00	26.6
26. 6. 11. 오전 3:30:00	26.7
26. 6. 11. 오전 4:00:00	26.8
26. 6. 11. 오전 4:30:00	26.9
26. 6. 11. 오전 5:00:00	26.9
26. 6. 11. 오전 5:30:00	27.0
26. 6. 11. 오전 6:00:00	27.1
26. 6. 11. 오전 6:30:00	26.9
26. 6. 11. 오전 7:00:00	26.4
26. 6. 11. 오전 7:30:00	26.1
26. 6. 11. 오전 8:00:00	25.9
26. 6. 11. 오전 8:30:00	25.8
26. 6. 11. 오전 9:00:00	25.7
26. 6. 11. 오전 9:30:00	25.7
26. 6. 11. 오전 10:00:00	25.7



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 11. 오전 10:30:00	25.7
26. 6. 11. 오전 11:00:00	25.7
26. 6. 11. 오전 11:30:00	25.7
26. 6. 11. 오후 12:00:00	25.8
26. 6. 11. 오후 12:30:00	25.9
26. 6. 11. 오후 1:00:00	25.9
26. 6. 11. 오후 1:30:00	25.9
26. 6. 11. 오후 2:00:00	26.0
26. 6. 11. 오후 2:30:00	26.0
26. 6. 11. 오후 3:00:00	26.0
26. 6. 11. 오후 3:30:00	26.0
26. 6. 11. 오후 4:00:00	25.9
26. 6. 11. 오후 4:30:00	25.8
26. 6. 11. 오후 5:00:00	25.8
26. 6. 11. 오후 5:30:00	25.8
26. 6. 11. 오후 6:00:00	25.7
26. 6. 11. 오후 6:30:00	25.6
26. 6. 11. 오후 7:00:00	25.5
26. 6. 11. 오후 7:30:00	25.4
26. 6. 11. 오후 8:00:00	25.3
26. 6. 11. 오후 8:30:00	25.3
26. 6. 11. 오후 9:00:00	25.2
26. 6. 11. 오후 9:30:00	25.1
26. 6. 11. 오후 10:00:00	25.1
26. 6. 11. 오후 10:30:00	25.0
26. 6. 11. 오후 11:00:00	24.9
26. 6. 11. 오후 11:30:00	24.8
26. 6. 12. 오전 12:00:00	24.9
26. 6. 12. 오전 12:30:00	25.4
26. 6. 12. 오전 1:00:00	25.8
26. 6. 12. 오전 1:30:00	26.2
26. 6. 12. 오전 2:00:00	26.3



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 12. 오전 2:30:00	26.5
26. 6. 12. 오전 3:00:00	26.6
26. 6. 12. 오전 3:30:00	26.7
26. 6. 12. 오전 4:00:00	26.8
26. 6. 12. 오전 4:30:00	26.9
26. 6. 12. 오전 5:00:00	27.0
26. 6. 12. 오전 5:30:00	27.1
26. 6. 12. 오전 6:00:00	27.0
26. 6. 12. 오전 6:30:00	26.6
26. 6. 12. 오전 7:00:00	26.2
26. 6. 12. 오전 7:30:00	25.9
26. 6. 12. 오전 8:00:00	25.8
26. 6. 12. 오전 8:30:00	25.7
26. 6. 12. 오전 9:00:00	25.6
26. 6. 12. 오전 9:30:00	25.5
26. 6. 12. 오전 10:00:00	25.5
26. 6. 12. 오전 10:30:00	25.5
26. 6. 12. 오전 11:00:00	25.5
26. 6. 12. 오전 11:30:00	25.6
26. 6. 12. 오후 12:00:00	25.6
26. 6. 12. 오후 12:30:00	25.7
26. 6. 12. 오후 1:00:00	25.8
26. 6. 12. 오후 1:30:00	25.8
26. 6. 12. 오후 2:00:00	25.9
26. 6. 12. 오후 2:30:00	25.9
26. 6. 12. 오후 2:50:30	26.1
26. 6. 12. 오후 3:00:00	26.9
26. 6. 12. 오후 3:30:00	26.1
26. 6. 12. 오후 4:00:00	25.8
26. 6. 12. 오후 4:30:00	25.6
26. 6. 12. 오후 5:00:00	25.6
26. 6. 12. 오후 5:30:00	25.6



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃) CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 12. 오후 6:00:00	25.6
26. 6. 12. 오후 6:30:00	25.6
26. 6. 12. 오후 7:00:00	25.5
26. 6. 12. 오후 7:30:00	25.3
26. 6. 12. 오후 8:00:00	25.3
26. 6. 12. 오후 8:30:00	25.2
26. 6. 12. 오후 9:00:00	25.2
26. 6. 12. 오후 9:30:00	25.0
26. 6. 12. 오후 10:00:00	25.0
26. 6. 12. 오후 10:30:00	24.8
26. 6. 12. 오후 11:00:00	24.8
26. 6. 12. 오후 11:30:00	24.9
26. 6. 13. 오전 12:00:00	24.9
26. 6. 13. 오전 12:30:00	25.6
26. 6. 13. 오전 1:00:00	25.9
26. 6. 13. 오전 1:30:00	26.4
26. 6. 13. 오전 2:00:00	26.6
26. 6. 13. 오전 2:30:00	26.6
26. 6. 13. 오전 3:00:00	26.6
26. 6. 13. 오전 3:30:00	26.9
26. 6. 13. 오전 4:00:00	26.9
26. 6. 13. 오전 4:30:00	26.9
26. 6. 13. 오전 5:00:00	26.9
26. 6. 13. 오전 5:30:00	27.0
26. 6. 13. 오전 6:00:00	27.2
26. 6. 13. 오전 6:30:00	26.4
26. 6. 13. 오전 7:00:00	25.9
26. 6. 13. 오전 7:30:00	25.7
26. 6. 13. 오전 8:00:00	25.7
26. 6. 13. 오전 8:30:00	25.5
26. 6. 13. 오전 9:00:00	25.4
26. 6. 13. 오전 9:30:00	25.4



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃) CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 13. 오전 10:00:00	25.5
26. 6. 13. 오전 10:30:00	25.5
26. 6. 13. 오전 11:00:00	25.4
26. 6. 13. 오전 11:30:00	25.5
26. 6. 13. 오후 12:00:00	25.9
26. 6. 13. 오후 12:30:00	26.0
26. 6. 13. 오후 1:00:00	25.9
26. 6. 13. 오후 1:30:00	26.0
26. 6. 13. 오후 2:00:00	26.1
26. 6. 13. 오후 2:30:00	26.1
26. 6. 13. 오후 3:00:00	26.3
26. 6. 13. 오후 3:30:00	26.3
26. 6. 13. 오후 4:00:00	26.0
26. 6. 13. 오후 4:30:00	25.9
26. 6. 13. 오후 5:00:00	25.8
26. 6. 13. 오후 5:30:00	25.9
26. 6. 13. 오후 6:00:00	25.9
26. 6. 13. 오후 6:30:00	25.9
26. 6. 13. 오후 7:00:00	25.7
26. 6. 13. 오후 7:30:00	25.6
26. 6. 13. 오후 8:00:00	25.5
26. 6. 13. 오후 8:30:00	25.3
26. 6. 13. 오후 9:00:00	25.2
26. 6. 13. 오후 9:30:00	25.2
26. 6. 13. 오후 10:00:00	25.2
26. 6. 13. 오후 10:30:00	25.2
26. 6. 13. 오후 11:00:00	25.2
26. 6. 13. 오후 11:30:00	25.2
26. 6. 14. 오전 12:00:00	25.1
26. 6. 14. 오전 12:30:00	26.3
26. 6. 14. 오전 1:00:00	26.8
26. 6. 14. 오전 1:30:00	27.3



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 14. 오전 2:00:00	27.4
26. 6. 14. 오전 2:30:00	27.3
26. 6. 14. 오전 3:00:00	27.6
26. 6. 14. 오전 3:30:00	27.8
26. 6. 14. 오전 4:00:00	27.8
26. 6. 14. 오전 4:30:00	27.8
26. 6. 14. 오전 5:00:00	27.9
26. 6. 14. 오전 5:30:00	28.1
26. 6. 14. 오전 6:00:00	28.1
26. 6. 14. 오전 6:30:00	26.7
26. 6. 14. 오전 7:00:00	26.5
26. 6. 14. 오전 7:30:00	26.3
26. 6. 14. 오전 8:00:00	26.2
26. 6. 14. 오전 8:30:00	26.1
26. 6. 14. 오전 9:00:00	26.0
26. 6. 14. 오전 9:30:00	26.1
26. 6. 14. 오전 10:00:00	26.2
26. 6. 14. 오전 10:30:00	26.2
26. 6. 14. 오전 11:00:00	26.1
26. 6. 14. 오전 11:30:00	26.3
26. 6. 14. 오후 12:00:00	26.4
26. 6. 14. 오후 12:30:00	26.4
26. 6. 14. 오후 1:00:00	26.6
26. 6. 14. 오후 1:30:00	26.6
26. 6. 14. 오후 2:00:00	26.5
26. 6. 14. 오후 2:30:00	26.8
26. 6. 14. 오후 3:00:00	26.5
26. 6. 14. 오후 3:30:00	26.4
26. 6. 14. 오후 4:00:00	26.4
26. 6. 14. 오후 4:30:00	26.3
26. 6. 14. 오후 5:00:00	26.1
26. 6. 14. 오후 5:30:00	26.0



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 14. 오후 6:00:00	26.0
26. 6. 14. 오후 6:30:00	25.8
26. 6. 14. 오후 7:00:00	26.0
26. 6. 14. 오후 7:30:00	25.8
26. 6. 14. 오후 8:00:00	25.7
26. 6. 14. 오후 8:30:00	25.7
26. 6. 14. 오후 9:00:00	25.5
26. 6. 14. 오후 9:30:00	25.4
26. 6. 14. 오후 10:00:00	25.4
26. 6. 14. 오후 10:30:00	25.4
26. 6. 14. 오후 11:00:00	25.4
26. 6. 14. 오후 11:30:00	25.3
26. 6. 15. 오전 12:00:00	25.3
26. 6. 15. 오전 12:30:00	26.3
26. 6. 15. 오전 1:00:00	26.8
26. 6. 15. 오전 1:30:00	27.1
26. 6. 15. 오전 2:00:00	27.3
26. 6. 15. 오전 2:30:00	27.4
26. 6. 15. 오전 3:00:00	27.6
26. 6. 15. 오전 3:30:00	27.7
26. 6. 15. 오전 4:00:00	27.7
26. 6. 15. 오전 4:30:00	27.8
26. 6. 15. 오전 5:00:00	27.9
26. 6. 15. 오전 5:30:00	28.1
26. 6. 15. 오전 6:00:00	28.1
26. 6. 15. 오전 6:30:00	27.2
26. 6. 15. 오전 7:00:00	26.7
26. 6. 15. 오전 7:30:00	26.4
26. 6. 15. 오전 8:00:00	26.3
26. 6. 15. 오전 8:30:00	26.3
26. 6. 15. 오전 9:00:00	26.1
26. 6. 15. 오전 9:30:00	26.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃) CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 15. 오전 10:00:00	26.2
26. 6. 15. 오전 10:30:00	26.2
26. 6. 15. 오전 11:00:00	26.4
26. 6. 15. 오전 11:30:00	26.7
26. 6. 15. 오후 12:00:00	27.0
26. 6. 15. 오후 12:30:00	27.1
26. 6. 15. 오후 1:00:00	27.3
26. 6. 15. 오후 1:30:00	27.1
26. 6. 15. 오후 2:00:00	26.7
26. 6. 15. 오후 2:30:00	26.9
26. 6. 15. 오후 3:00:00	26.7
26. 6. 15. 오후 3:30:00	26.5
26. 6. 15. 오후 4:00:00	26.5
26. 6. 15. 오후 4:30:00	26.3
26. 6. 15. 오후 5:00:00	26.3
26. 6. 15. 오후 5:30:00	26.3
26. 6. 15. 오후 6:00:00	26.2
26. 6. 15. 오후 6:30:00	26.4
26. 6. 15. 오후 7:00:00	26.2
26. 6. 15. 오후 7:30:00	26.3
26. 6. 15. 오후 8:00:00	26.2
26. 6. 15. 오후 8:30:00	26.2
26. 6. 15. 오후 9:00:00	26.2
26. 6. 15. 오후 9:30:00	26.1
26. 6. 15. 오후 10:00:00	26.0
26. 6. 15. 오후 10:30:00	25.8
26. 6. 15. 오후 11:00:00	25.7
26. 6. 15. 오후 11:30:00	25.7
26. 6. 16. 오전 12:00:00	26.1
26. 6. 16. 오전 12:30:00	26.9
26. 6. 16. 오전 1:00:00	27.4
26. 6. 16. 오전 1:30:00	27.6



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 16. 오전 2:00:00	27.8
26. 6. 16. 오전 2:30:00	27.9
26. 6. 16. 오전 3:00:00	28.1
26. 6. 16. 오전 3:30:00	28.2
26. 6. 16. 오전 4:00:00	28.2
26. 6. 16. 오전 4:30:00	28.4
26. 6. 16. 오전 5:00:00	28.6
26. 6. 16. 오전 5:30:00	28.5
26. 6. 16. 오전 6:00:00	28.5
26. 6. 16. 오전 6:30:00	27.5
26. 6. 16. 오전 7:00:00	27.2
26. 6. 16. 오전 7:30:00	27.0
26. 6. 16. 오전 8:00:00	26.7
26. 6. 16. 오전 8:30:00	26.7
26. 6. 16. 오전 9:00:00	26.7
26. 6. 16. 오전 9:30:00	26.8
26. 6. 16. 오전 10:00:00	26.6
26. 6. 16. 오전 10:30:00	26.8
26. 6. 16. 오전 11:00:00	26.9
26. 6. 16. 오전 11:30:00	27.1
26. 6. 16. 오후 12:00:00	25.8
26. 6. 16. 오후 12:30:00	25.1
26. 6. 16. 오후 1:00:00	25.5
26. 6. 16. 오후 1:30:00	24.9
26. 6. 16. 오후 2:00:00	25.0
26. 6. 16. 오후 2:30:00	25.2
26. 6. 16. 오후 3:00:00	24.7
26. 6. 16. 오후 3:30:00	24.8
26. 6. 16. 오후 4:00:00	24.8
26. 6. 16. 오후 4:30:00	25.0
26. 6. 16. 오후 5:00:00	24.9
26. 6. 16. 오후 5:30:00	24.8



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃) CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 16. 오후 6:00:00	24.5
26. 6. 16. 오후 6:30:00	24.6
26. 6. 16. 오후 7:00:00	24.5
26. 6. 16. 오후 7:30:00	24.4
26. 6. 16. 오후 8:00:00	24.2
26. 6. 16. 오후 8:30:00	24.4
26. 6. 16. 오후 9:00:00	24.7
26. 6. 16. 오후 9:30:00	24.5
26. 6. 16. 오후 10:00:00	24.2
26. 6. 16. 오후 10:30:00	24.6
26. 6. 16. 오후 11:00:00	24.4
26. 6. 16. 오후 11:30:00	24.0
26. 6. 17. 오전 12:00:00	24.4
26. 6. 17. 오전 12:30:00	25.2
26. 6. 17. 오전 1:00:00	25.1
26. 6. 17. 오전 1:30:00	24.5
26. 6. 17. 오전 2:00:00	25.1
26. 6. 17. 오전 2:30:00	24.7
26. 6. 17. 오전 3:00:00	25.3
26. 6. 17. 오전 3:30:00	25.0
26. 6. 17. 오전 4:00:00	24.8
26. 6. 17. 오전 4:30:00	25.1
26. 6. 17. 오전 5:00:00	25.0
26. 6. 17. 오전 5:30:00	25.4
26. 6. 17. 오전 6:00:00	24.7
26. 6. 17. 오전 6:30:00	24.5
26. 6. 17. 오전 7:00:00	24.4
26. 6. 17. 오전 7:30:00	24.8
26. 6. 17. 오전 8:00:00	24.7
26. 6. 17. 오전 8:30:00	24.7
26. 6. 17. 오전 9:00:00	24.5
26. 6. 17. 오전 9:30:00	24.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 17. 오전 10:00:00	24.6
26. 6. 17. 오전 10:30:00	24.5
26. 6. 17. 오전 11:00:00	24.5
26. 6. 17. 오전 11:30:00	24.4
26. 6. 17. 오후 12:00:00	24.4
26. 6. 17. 오후 12:30:00	24.6
26. 6. 17. 오후 1:00:00	24.5
26. 6. 17. 오후 1:30:00	24.5
26. 6. 17. 오후 2:00:00	24.7
26. 6. 17. 오후 2:30:00	24.7
26. 6. 17. 오후 3:00:00	24.9
26. 6. 17. 오후 3:30:00	24.8
26. 6. 17. 오후 4:00:00	24.8
26. 6. 17. 오후 4:30:00	24.7
26. 6. 17. 오후 5:00:00	24.5
26. 6. 17. 오후 5:30:00	24.2
26. 6. 17. 오후 6:00:00	24.4
26. 6. 17. 오후 6:30:00	24.5
26. 6. 17. 오후 7:00:00	24.8
26. 6. 17. 오후 7:30:00	24.7
26. 6. 17. 오후 8:00:00	24.4
26. 6. 17. 오후 8:30:00	24.1
26. 6. 17. 오후 9:00:00	24.6
26. 6. 17. 오후 9:30:00	24.5
26. 6. 17. 오후 10:00:00	24.0
26. 6. 17. 오후 10:30:00	24.5
26. 6. 17. 오후 11:00:00	24.2
26. 6. 17. 오후 11:30:00	24.1
26. 6. 18. 오전 12:00:00	24.3
26. 6. 18. 오전 12:30:00	24.7
26. 6. 18. 오전 1:00:00	25.0
26. 6. 18. 오전 1:30:00	24.7



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 18. 오전 2:00:00	24.6
26. 6. 18. 오전 2:30:00	25.1
26. 6. 18. 오전 3:00:00	25.0
26. 6. 18. 오전 3:30:00	24.7
26. 6. 18. 오전 4:00:00	25.2
26. 6. 18. 오전 4:30:00	24.9
26. 6. 18. 오전 5:00:00	24.7
26. 6. 18. 오전 5:30:00	25.3
26. 6. 18. 오전 6:00:00	24.9
26. 6. 18. 오전 6:30:00	24.8
26. 6. 18. 오전 7:00:00	24.8
26. 6. 18. 오전 7:30:00	24.6
26. 6. 18. 오전 8:00:00	24.4
26. 6. 18. 오전 8:30:00	24.1
26. 6. 18. 오전 9:00:00	24.5
26. 6. 18. 오전 9:30:00	24.3
26. 6. 18. 오전 10:00:00	24.5
26. 6. 18. 오전 10:30:00	24.4
26. 6. 18. 오전 11:00:00	24.1
26. 6. 18. 오전 11:30:00	24.4
26. 6. 18. 오후 12:00:00	24.5
26. 6. 18. 오후 12:30:00	24.5
26. 6. 18. 오후 1:00:00	24.5
26. 6. 18. 오후 1:30:00	24.2
26. 6. 18. 오후 2:00:00	24.2
26. 6. 18. 오후 2:30:00	24.3
26. 6. 18. 오후 3:00:00	24.5
26. 6. 18. 오후 3:30:00	24.6
26. 6. 18. 오후 4:00:00	24.4
26. 6. 18. 오후 4:30:00	24.1
26. 6. 18. 오후 5:00:00	23.9
26. 6. 18. 오후 5:30:00	24.3



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃) CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 18. 오후 6:00:00	24.3
26. 6. 18. 오후 6:30:00	23.9
26. 6. 18. 오후 7:00:00	24.1
26. 6. 18. 오후 7:30:00	24.5
26. 6. 18. 오후 8:00:00	24.1
26. 6. 18. 오후 8:30:00	23.7
26. 6. 18. 오후 9:00:00	24.2
26. 6. 18. 오후 9:30:00	23.4
26. 6. 18. 오후 10:00:00	24.1
26. 6. 18. 오후 10:30:00	24.2
26. 6. 18. 오후 11:00:00	24.0
26. 6. 18. 오후 11:30:00	24.2
26. 6. 19. 오전 12:00:00	23.4
26. 6. 19. 오전 12:30:00	24.5
26. 6. 19. 오전 1:00:00	24.7
26. 6. 19. 오전 1:30:00	24.9
26. 6. 19. 오전 2:00:00	24.8
26. 6. 19. 오전 2:30:00	24.2
26. 6. 19. 오전 3:00:00	24.4
26. 6. 19. 오전 3:30:00	25.0
26. 6. 19. 오전 4:00:00	24.7
26. 6. 19. 오전 4:30:00	24.5
26. 6. 19. 오전 5:00:00	24.8
26. 6. 19. 오전 5:30:00	25.1
26. 6. 19. 오전 6:00:00	24.9
26. 6. 19. 오전 6:30:00	24.6
26. 6. 19. 오전 7:00:00	24.2
26. 6. 19. 오전 7:30:00	23.8
26. 6. 19. 오전 8:00:00	24.6
26. 6. 19. 오전 8:30:00	24.1
26. 6. 19. 오전 9:00:00	23.8
26. 6. 19. 오전 9:30:00	24.2



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃) CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 19. 오전 10:00:00	23.8
26. 6. 19. 오전 10:30:00	24.4
26. 6. 19. 오전 11:00:00	24.4
26. 6. 19. 오전 11:30:00	24.3
26. 6. 19. 오후 12:00:00	24.2
26. 6. 19. 오후 12:30:00	24.2
26. 6. 19. 오후 1:00:00	24.1
26. 6. 19. 오후 1:30:00	23.9
26. 6. 19. 오후 2:00:00	23.7
26. 6. 19. 오후 2:30:00	24.1
26. 6. 19. 오후 3:00:00	24.2
26. 6. 19. 오후 3:30:00	24.5
26. 6. 19. 오후 4:00:00	24.3
26. 6. 19. 오후 4:30:00	23.8
26. 6. 19. 오후 5:00:00	24.4
26. 6. 19. 오후 5:30:00	24.1
26. 6. 19. 오후 6:00:00	24.3
26. 6. 19. 오후 6:30:00	24.4
26. 6. 19. 오후 7:00:00	23.7
26. 6. 19. 오후 7:30:00	24.3
26. 6. 19. 오후 8:00:00	24.1
26. 6. 19. 오후 8:30:00	24.0
26. 6. 19. 오후 9:00:00	24.3
26. 6. 19. 오후 9:30:00	23.4
26. 6. 19. 오후 10:00:00	24.0
26. 6. 19. 오후 10:30:00	24.2
26. 6. 19. 오후 11:00:00	23.8
26. 6. 19. 오후 11:30:00	23.9
26. 6. 20. 오전 12:00:00	24.8
26. 6. 20. 오전 12:30:00	24.9
26. 6. 20. 오전 1:00:00	24.3
26. 6. 20. 오전 1:30:00	24.3



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 20. 오전 2:00:00	23.9
26. 6. 20. 오전 2:30:00	24.5
26. 6. 20. 오전 3:00:00	23.6
26. 6. 20. 오전 3:30:00	24.0
26. 6. 20. 오전 4:00:00	24.3
26. 6. 20. 오전 4:30:00	24.1
26. 6. 20. 오전 5:00:00	23.7
26. 6. 20. 오전 5:30:00	24.0
26. 6. 20. 오전 6:00:00	24.2
26. 6. 20. 오전 6:30:00	24.3
26. 6. 20. 오전 7:00:00	23.8
26. 6. 20. 오전 7:30:00	23.7
26. 6. 20. 오전 8:00:00	24.1
26. 6. 20. 오전 8:30:00	23.9
26. 6. 20. 오전 9:00:00	23.7
26. 6. 20. 오전 9:30:00	23.8
26. 6. 20. 오전 10:00:00	23.6
26. 6. 20. 오전 10:30:00	23.5
26. 6. 20. 오전 11:00:00	23.2
26. 6. 20. 오전 11:30:00	23.2
26. 6. 20. 오후 12:00:00	23.5
26. 6. 20. 오후 12:30:00	23.6
26. 6. 20. 오후 1:00:00	23.2
26. 6. 20. 오후 1:30:00	24.0
26. 6. 20. 오후 2:00:00	24.0
26. 6. 20. 오후 2:30:00	24.0
26. 6. 20. 오후 3:00:00	24.0
26. 6. 20. 오후 3:30:00	24.0
26. 6. 20. 오후 4:00:00	23.8
26. 6. 20. 오후 4:30:00	24.0
26. 6. 20. 오후 5:00:00	24.0
26. 6. 20. 오후 5:30:00	23.6



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃) CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 20. 오후 6:00:00	23.5
26. 6. 20. 오후 6:30:00	23.6
26. 6. 20. 오후 7:00:00	23.7
26. 6. 20. 오후 7:30:00	23.7
26. 6. 20. 오후 8:00:00	23.8
26. 6. 20. 오후 8:30:00	23.2
26. 6. 20. 오후 9:00:00	23.3
26. 6. 20. 오후 9:30:00	24.0
26. 6. 20. 오후 10:00:00	24.0
26. 6. 20. 오후 10:30:00	24.0
26. 6. 20. 오후 11:00:00	24.0
26. 6. 20. 오후 11:30:00	24.0
26. 6. 21. 오전 12:00:00	24.2
26. 6. 21. 오전 12:30:00	24.6
26. 6. 21. 오전 1:00:00	24.8
26. 6. 21. 오전 1:30:00	24.6
26. 6. 21. 오전 2:00:00	24.2
26. 6. 21. 오전 2:30:00	23.9
26. 6. 21. 오전 3:00:00	24.1
26. 6. 21. 오전 3:30:00	24.1
26. 6. 21. 오전 4:00:00	24.3
26. 6. 21. 오전 4:30:00	24.5
26. 6. 21. 오전 5:00:00	24.8
26. 6. 21. 오전 5:30:00	24.8
26. 6. 21. 오전 6:00:00	24.6
26. 6. 21. 오전 6:30:00	24.5
26. 6. 21. 오전 7:00:00	23.8
26. 6. 21. 오전 7:30:00	24.1
26. 6. 21. 오전 8:00:00	24.1
26. 6. 21. 오전 8:30:00	23.6
26. 6. 21. 오전 9:00:00	23.7
26. 6. 21. 오전 9:30:00	24.0



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 21. 오전 10:00:00	23.4
26. 6. 21. 오전 10:30:00	24.2
26. 6. 21. 오전 11:00:00	23.5
26. 6. 21. 오전 11:30:00	24.1
26. 6. 21. 오후 12:00:00	23.9
26. 6. 21. 오후 12:30:00	23.7
26. 6. 21. 오후 1:00:00	24.1
26. 6. 21. 오후 1:30:00	24.3
26. 6. 21. 오후 2:00:00	24.0
26. 6. 21. 오후 2:30:00	23.6
26. 6. 21. 오후 3:00:00	24.1
26. 6. 21. 오후 3:30:00	24.3
26. 6. 21. 오후 4:00:00	24.0
26. 6. 21. 오후 4:30:00	24.1
26. 6. 21. 오후 5:00:00	24.0
26. 6. 21. 오후 5:30:00	24.2
26. 6. 21. 오후 6:00:00	23.7
26. 6. 21. 오후 6:30:00	24.1
26. 6. 21. 오후 7:00:00	23.8
26. 6. 21. 오후 7:30:00	23.6
26. 6. 21. 오후 8:00:00	23.9
26. 6. 21. 오후 8:30:00	24.0
26. 6. 21. 오후 9:00:00	24.1
26. 6. 21. 오후 9:30:00	24.0
26. 6. 21. 오후 10:00:00	24.0
26. 6. 21. 오후 10:30:00	24.0
26. 6. 21. 오후 11:00:00	24.0
26. 6. 21. 오후 11:30:00	24.0
26. 6. 22. 오전 12:00:00	23.9
26. 6. 22. 오전 12:30:00	24.9
26. 6. 22. 오전 1:00:00	24.9
26. 6. 22. 오전 1:30:00	24.9



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 22. 오전 2:00:00	24.6
26. 6. 22. 오전 2:30:00	24.3
26. 6. 22. 오전 3:00:00	24.1
26. 6. 22. 오전 3:30:00	24.2
26. 6. 22. 오전 4:00:00	24.4
26. 6. 22. 오전 4:30:00	24.7
26. 6. 22. 오전 5:00:00	24.8
26. 6. 22. 오전 5:30:00	24.8
26. 6. 22. 오전 6:00:00	24.9
26. 6. 22. 오전 6:30:00	24.6
26. 6. 22. 오전 7:00:00	24.2
26. 6. 22. 오전 7:30:00	24.3
26. 6. 22. 오전 8:00:00	24.5
26. 6. 22. 오전 8:30:00	23.5
26. 6. 22. 오전 9:00:00	24.0
26. 6. 22. 오전 9:30:00	24.2
26. 6. 22. 오전 10:00:00	23.5
26. 6. 22. 오전 10:30:00	24.1
26. 6. 22. 오전 11:00:00	23.6
26. 6. 22. 오전 11:30:00	24.1
26. 6. 22. 오후 12:00:00	23.5
26. 6. 22. 오후 12:30:00	24.2
26. 6. 22. 오후 1:00:00	23.6
26. 6. 22. 오후 1:30:00	24.1
26. 6. 22. 오후 2:00:00	23.3
26. 6. 22. 오후 2:30:00	24.0
26. 6. 22. 오후 3:00:00	23.6
26. 6. 22. 오후 3:30:00	24.1
26. 6. 22. 오후 4:00:00	23.2
26. 6. 22. 오후 4:30:00	24.0
26. 6. 22. 오후 5:00:00	24.1
26. 6. 22. 오후 5:30:00	23.6



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 22. 오후 6:00:00	24.2
26. 6. 22. 오후 6:30:00	23.4
26. 6. 22. 오후 7:00:00	24.1
26. 6. 22. 오후 7:30:00	24.0
26. 6. 22. 오후 8:00:00	24.0
26. 6. 22. 오후 8:30:00	23.4
26. 6. 22. 오후 9:00:00	23.6
26. 6. 22. 오후 9:30:00	23.3
26. 6. 22. 오후 10:00:00	23.1
26. 6. 22. 오후 10:30:00	23.3
26. 6. 22. 오후 11:00:00	24.0
26. 6. 22. 오후 11:30:00	23.9
26. 6. 23. 오전 12:00:00	23.9
26. 6. 23. 오전 12:30:00	23.9
26. 6. 23. 오전 1:00:00	24.1
26. 6. 23. 오전 1:30:00	24.0
26. 6. 23. 오전 2:00:00	23.9
26. 6. 23. 오전 2:30:00	24.2
26. 6. 23. 오전 3:00:00	24.2
26. 6. 23. 오전 3:30:00	24.4
26. 6. 23. 오전 4:00:00	24.5
26. 6. 23. 오전 4:30:00	24.7
26. 6. 23. 오전 5:00:00	24.8
26. 6. 23. 오전 5:30:00	24.5
26. 6. 23. 오전 6:00:00	24.7
26. 6. 23. 오전 6:30:00	24.1
26. 6. 23. 오전 7:00:00	24.3
26. 6. 23. 오전 7:30:00	23.4
26. 6. 23. 오전 8:00:00	23.6
26. 6. 23. 오전 8:30:00	23.7
26. 6. 23. 오전 9:00:00	23.7
26. 6. 23. 오전 9:30:00	23.7



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃) CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 23. 오전 10:00:00	23.7
26. 6. 23. 오전 10:30:00	24.0
26. 6. 23. 오전 11:00:00	24.1
26. 6. 23. 오전 11:30:00	23.3
26. 6. 23. 오후 12:00:00	23.9
26. 6. 23. 오후 12:30:00	23.6
26. 6. 23. 오후 1:00:00	23.9
26. 6. 23. 오후 1:30:00	23.3
26. 6. 23. 오후 2:00:00	24.0
26. 6. 23. 오후 2:30:00	23.4
26. 6. 23. 오후 3:00:00	23.9
26. 6. 23. 오후 3:30:00	23.4
26. 6. 23. 오후 4:00:00	23.8
26. 6. 23. 오후 4:30:00	23.9
26. 6. 23. 오후 5:00:00	23.0
26. 6. 23. 오후 5:30:00	23.6
26. 6. 23. 오후 6:00:00	24.1
26. 6. 23. 오후 6:30:00	24.0
26. 6. 23. 오후 7:00:00	23.4
26. 6. 23. 오후 7:30:00	24.2
26. 6. 23. 오후 8:00:00	23.7
26. 6. 23. 오후 8:30:00	24.1
26. 6. 23. 오후 9:00:00	23.9
26. 6. 23. 오후 9:30:00	23.8
26. 6. 23. 오후 10:00:00	24.0
26. 6. 23. 오후 10:30:00	24.0
26. 6. 23. 오후 11:00:00	24.1
26. 6. 23. 오후 11:30:00	23.7
26. 6. 24. 오전 12:00:00	23.5
26. 6. 24. 오전 12:30:00	23.9
26. 6. 24. 오전 1:00:00	23.8
26. 6. 24. 오전 1:30:00	24.0



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 24. 오전 2:00:00	24.5
26. 6. 24. 오전 2:30:00	24.9
26. 6. 24. 오전 3:00:00	24.6
26. 6. 24. 오전 3:30:00	24.1
26. 6. 24. 오전 4:00:00	24.4
26. 6. 24. 오전 4:30:00	24.8
26. 6. 24. 오전 5:00:00	24.1
26. 6. 24. 오전 5:30:00	24.2
26. 6. 24. 오전 6:00:00	23.8
26. 6. 24. 오전 6:30:00	23.3
26. 6. 24. 오전 7:00:00	23.8
26. 6. 24. 오전 7:30:00	23.8
26. 6. 24. 오전 8:00:00	23.9
26. 6. 24. 오전 8:30:00	23.9
26. 6. 24. 오전 9:00:00	24.1
26. 6. 24. 오전 9:30:00	24.2
26. 6. 24. 오전 10:00:00	24.2
26. 6. 24. 오전 10:30:00	24.3
26. 6. 24. 오전 11:00:00	24.5
26. 6. 24. 오전 11:30:00	24.1
26. 6. 24. 오후 12:00:00	23.5
26. 6. 24. 오후 12:30:00	23.5
26. 6. 24. 오후 1:00:00	24.2
26. 6. 24. 오후 1:30:00	24.2
26. 6. 24. 오후 2:00:00	24.0
26. 6. 24. 오후 2:30:00	23.5
26. 6. 24. 오후 3:00:00	24.1
26. 6. 24. 오후 3:30:00	23.9
26. 6. 24. 오후 4:00:00	23.7
26. 6. 24. 오후 4:30:00	23.9
26. 6. 24. 오후 5:00:00	24.0
26. 6. 24. 오후 5:30:00	23.2



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 24. 오후 6:00:00	24.2
26. 6. 24. 오후 6:30:00	23.6
26. 6. 24. 오후 7:00:00	24.0
26. 6. 24. 오후 7:30:00	23.6
26. 6. 24. 오후 8:00:00	24.1
26. 6. 24. 오후 8:30:00	23.3
26. 6. 24. 오후 9:00:00	23.3
26. 6. 24. 오후 9:30:00	23.8
26. 6. 24. 오후 10:00:00	23.9
26. 6. 24. 오후 10:30:00	24.0
26. 6. 24. 오후 11:00:00	24.0
26. 6. 24. 오후 11:30:00	23.9
26. 6. 25. 오전 12:00:00	24.0
26. 6. 25. 오전 12:30:00	23.8
26. 6. 25. 오전 1:00:00	24.0
26. 6. 25. 오전 1:30:00	24.7
26. 6. 25. 오전 2:00:00	24.8
26. 6. 25. 오전 2:30:00	24.3
26. 6. 25. 오전 3:00:00	23.8
26. 6. 25. 오전 3:30:00	24.3
26. 6. 25. 오전 4:00:00	25.0
26. 6. 25. 오전 4:30:00	24.4
26. 6. 25. 오전 5:00:00	24.1
26. 6. 25. 오전 5:30:00	24.6
26. 6. 25. 오전 6:00:00	24.7
26. 6. 25. 오전 6:30:00	24.5
26. 6. 25. 오전 7:00:00	24.0
26. 6. 25. 오전 7:30:00	24.2
26. 6. 25. 오전 8:00:00	24.0
26. 6. 25. 오전 8:30:00	24.1
26. 6. 25. 오전 9:00:00	23.3
26. 6. 25. 오전 9:30:00	23.8



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 25. 오전 10:00:00	24.0
26. 6. 25. 오전 10:30:00	23.6
26. 6. 25. 오전 11:00:00	24.0
26. 6. 25. 오전 11:30:00	23.8
26. 6. 25. 오후 12:00:00	23.8
26. 6. 25. 오후 12:30:00	24.2
26. 6. 25. 오후 1:00:00	23.9
26. 6. 25. 오후 1:30:00	23.4
26. 6. 25. 오후 2:00:00	24.2
26. 6. 25. 오후 2:30:00	23.8
26. 6. 25. 오후 3:00:00	24.1
26. 6. 25. 오후 3:30:00	23.3
26. 6. 25. 오후 4:00:00	23.8
26. 6. 25. 오후 4:30:00	23.9
26. 6. 25. 오후 5:00:00	23.9
26. 6. 25. 오후 5:30:00	23.8
26. 6. 25. 오후 6:00:00	24.1
26. 6. 25. 오후 6:30:00	23.7
26. 6. 25. 오후 7:00:00	23.8
26. 6. 25. 오후 7:30:00	24.0
26. 6. 25. 오후 8:00:00	24.0
26. 6. 25. 오후 8:30:00	24.0
26. 6. 25. 오후 9:00:00	23.4
26. 6. 25. 오후 9:30:00	23.1
26. 6. 25. 오후 10:00:00	23.3
26. 6. 25. 오후 10:30:00	23.2
26. 6. 25. 오후 11:00:00	23.3
26. 6. 25. 오후 11:30:00	23.8
26. 6. 26. 오전 12:00:00	24.2
26. 6. 26. 오전 12:30:00	24.0
26. 6. 26. 오전 1:00:00	23.9
26. 6. 26. 오전 1:30:00	23.8



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 26. 오전 2:00:00	24.2
26. 6. 26. 오전 2:30:00	24.4
26. 6. 26. 오전 3:00:00	24.4
26. 6. 26. 오전 3:30:00	24.6
26. 6. 26. 오전 4:00:00	24.6
26. 6. 26. 오전 4:30:00	24.6
26. 6. 26. 오전 5:00:00	24.6
26. 6. 26. 오전 5:30:00	24.6
26. 6. 26. 오전 6:00:00	24.2
26. 6. 26. 오전 6:30:00	23.4
26. 6. 26. 오전 7:00:00	23.3
26. 6. 26. 오전 7:30:00	24.0
26. 6. 26. 오전 8:00:00	23.9
26. 6. 26. 오전 8:30:00	24.0
26. 6. 26. 오전 9:00:00	24.1
26. 6. 26. 오전 9:30:00	24.2
26. 6. 26. 오전 10:00:00	24.3
26. 6. 26. 오전 10:30:00	24.4
26. 6. 26. 오전 11:00:00	23.9
26. 6. 26. 오전 11:30:00	24.2
26. 6. 26. 오후 12:00:00	23.6
26. 6. 26. 오후 12:30:00	24.2
26. 6. 26. 오후 1:00:00	24.0
26. 6. 26. 오후 1:30:00	23.5
26. 6. 26. 오후 2:00:00	24.2
26. 6. 26. 오후 2:30:00	23.9
26. 6. 26. 오후 3:00:00	23.8
26. 6. 26. 오후 3:30:00	23.8
26. 6. 26. 오후 4:00:00	23.6
26. 6. 26. 오후 4:30:00	23.7
26. 6. 26. 오후 5:00:00	24.0
26. 6. 26. 오후 5:30:00	23.3



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 26. 오후 6:00:00	23.9
26. 6. 26. 오후 6:30:00	24.1
26. 6. 26. 오후 7:00:00	23.8
26. 6. 26. 오후 7:30:00	24.0
26. 6. 26. 오후 8:00:00	24.1
26. 6. 26. 오후 8:30:00	23.3
26. 6. 26. 오후 9:00:00	23.6
26. 6. 26. 오후 9:30:00	23.5
26. 6. 26. 오후 10:00:00	23.7
26. 6. 26. 오후 10:30:00	23.7
26. 6. 26. 오후 11:00:00	23.8
26. 6. 26. 오후 11:30:00	23.8
26. 6. 27. 오전 12:00:00	23.4
26. 6. 27. 오전 12:30:00	24.6
26. 6. 27. 오전 1:00:00	24.5
26. 6. 27. 오전 1:30:00	24.6
26. 6. 27. 오전 2:00:00	24.6
26. 6. 27. 오전 2:30:00	24.6
26. 6. 27. 오전 3:00:00	24.8
26. 6. 27. 오전 3:30:00	23.8
26. 6. 27. 오전 4:00:00	24.0
26. 6. 27. 오전 4:30:00	24.0
26. 6. 27. 오전 5:00:00	23.9
26. 6. 27. 오전 5:30:00	23.9
26. 6. 27. 오전 6:00:00	23.9
26. 6. 27. 오전 6:30:00	23.9
26. 6. 27. 오전 7:00:00	23.8
26. 6. 27. 오전 7:30:00	23.7
26. 6. 27. 오전 8:00:00	23.9
26. 6. 27. 오전 8:30:00	23.9
26. 6. 27. 오전 9:00:00	24.0
26. 6. 27. 오전 9:30:00	23.2



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 27. 오전 10:00:00	23.8
26. 6. 27. 오전 10:30:00	24.0
26. 6. 27. 오전 11:00:00	23.8
26. 6. 27. 오전 11:30:00	24.0
26. 6. 27. 오후 12:00:00	24.1
26. 6. 27. 오후 12:30:00	23.7
26. 6. 27. 오후 1:00:00	23.6
26. 6. 27. 오후 1:30:00	24.1
26. 6. 27. 오후 2:00:00	24.4
26. 6. 27. 오후 2:30:00	24.4
26. 6. 27. 오후 3:00:00	24.5
26. 6. 27. 오후 3:30:00	24.2
26. 6. 27. 오후 4:00:00	24.0
26. 6. 27. 오후 4:30:00	23.7
26. 6. 27. 오후 5:00:00	23.9
26. 6. 27. 오후 5:30:00	24.5
26. 6. 27. 오후 6:00:00	24.4
26. 6. 27. 오후 6:30:00	24.0
26. 6. 27. 오후 7:00:00	23.3
26. 6. 27. 오후 7:30:00	24.3
26. 6. 27. 오후 8:00:00	24.1
26. 6. 27. 오후 8:30:00	23.2
26. 6. 27. 오후 9:00:00	24.2
26. 6. 27. 오후 9:30:00	23.5
26. 6. 27. 오후 10:00:00	24.0
26. 6. 27. 오후 10:30:00	24.0
26. 6. 27. 오후 11:00:00	23.7
26. 6. 27. 오후 11:30:00	24.0
26. 6. 28. 오전 12:00:00	24.0
26. 6. 28. 오전 12:30:00	24.4
26. 6. 28. 오전 1:00:00	24.5
26. 6. 28. 오전 1:30:00	24.5



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 28. 오전 2:00:00	24.7
26. 6. 28. 오전 2:30:00	24.9
26. 6. 28. 오전 3:00:00	24.6
26. 6. 28. 오전 3:30:00	24.2
26. 6. 28. 오전 4:00:00	24.5
26. 6. 28. 오전 4:30:00	24.8
26. 6. 28. 오전 5:00:00	24.9
26. 6. 28. 오전 5:30:00	24.4
26. 6. 28. 오전 6:00:00	24.3
26. 6. 28. 오전 6:30:00	23.9
26. 6. 28. 오전 7:00:00	24.3
26. 6. 28. 오전 7:30:00	23.4
26. 6. 28. 오전 8:00:00	24.2
26. 6. 28. 오전 8:30:00	23.6
26. 6. 28. 오전 9:00:00	24.2
26. 6. 28. 오전 9:30:00	23.6
26. 6. 28. 오전 10:00:00	24.3
26. 6. 28. 오전 10:30:00	24.1
26. 6. 28. 오전 11:00:00	23.8
26. 6. 28. 오전 11:30:00	23.7
26. 6. 28. 오후 12:00:00	24.2
26. 6. 28. 오후 12:30:00	24.6
26. 6. 28. 오후 1:00:00	24.7
26. 6. 28. 오후 1:30:00	24.6
26. 6. 28. 오후 2:00:00	24.8
26. 6. 28. 오후 2:30:00	24.7
26. 6. 28. 오후 3:00:00	24.6
26. 6. 28. 오후 3:30:00	24.5
26. 6. 28. 오후 4:00:00	24.3
26. 6. 28. 오후 4:30:00	24.2
26. 6. 28. 오후 5:00:00	23.9
26. 6. 28. 오후 5:30:00	23.6



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 28. 오후 6:00:00	24.4
26. 6. 28. 오후 6:30:00	24.2
26. 6. 28. 오후 7:00:00	23.8
26. 6. 28. 오후 7:30:00	24.0
26. 6. 28. 오후 8:00:00	24.2
26. 6. 28. 오후 8:30:00	23.7
26. 6. 28. 오후 9:00:00	24.3
26. 6. 28. 오후 9:30:00	23.7
26. 6. 28. 오후 10:00:00	24.1
26. 6. 28. 오후 10:30:00	23.5
26. 6. 28. 오후 11:00:00	24.2
26. 6. 28. 오후 11:30:00	23.4
26. 6. 29. 오전 12:00:00	24.0
26. 6. 29. 오전 12:30:00	24.2
26. 6. 29. 오전 1:00:00	24.1
26. 6. 29. 오전 1:30:00	24.2
26. 6. 29. 오전 2:00:00	24.2
26. 6. 29. 오전 2:30:00	24.6
26. 6. 29. 오전 3:00:00	24.9
26. 6. 29. 오전 3:30:00	24.7
26. 6. 29. 오전 4:00:00	24.4
26. 6. 29. 오전 4:30:00	24.1
26. 6. 29. 오전 5:00:00	24.2
26. 6. 29. 오전 5:30:00	24.4
26. 6. 29. 오전 6:00:00	25.0
26. 6. 29. 오전 6:30:00	25.0
26. 6. 29. 오전 7:00:00	24.6
26. 6. 29. 오전 7:30:00	24.2
26. 6. 29. 오전 8:00:00	23.6
26. 6. 29. 오전 8:30:00	24.3
26. 6. 29. 오전 9:00:00	24.0
26. 6. 29. 오전 9:30:00	24.3



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 29. 오전 10:00:00	24.1
26. 6. 29. 오전 10:30:00	23.8
26. 6. 29. 오전 11:00:00	24.4
26. 6. 29. 오전 11:30:00	24.3
26. 6. 29. 오후 12:00:00	24.1
26. 6. 29. 오후 12:30:00	23.9
26. 6. 29. 오후 1:00:00	23.8
26. 6. 29. 오후 1:30:00	23.9
26. 6. 29. 오후 2:00:00	24.1
26. 6. 29. 오후 2:30:00	24.1
26. 6. 29. 오후 3:00:00	24.2
26. 6. 29. 오후 3:30:00	24.5
26. 6. 29. 오후 4:00:00	24.6
26. 6. 29. 오후 4:30:00	24.5
26. 6. 29. 오후 5:00:00	24.4
26. 6. 29. 오후 5:30:00	24.2
26. 6. 29. 오후 6:00:00	24.1
26. 6. 29. 오후 6:30:00	24.1
26. 6. 29. 오후 7:00:00	23.9
26. 6. 29. 오후 7:30:00	23.9
26. 6. 29. 오후 8:00:00	24.5
26. 6. 29. 오후 8:30:00	24.5
26. 6. 29. 오후 9:00:00	23.8
26. 6. 29. 오후 9:30:00	24.5
26. 6. 29. 오후 10:00:00	24.2
26. 6. 29. 오후 10:30:00	24.3
26. 6. 29. 오후 11:00:00	24.3
26. 6. 29. 오후 11:30:00	23.5
26. 6. 30. 오전 12:00:00	24.8
26. 6. 30. 오전 12:30:00	24.8
26. 6. 30. 오전 1:00:00	24.4
26. 6. 30. 오전 1:30:00	24.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 30. 오전 2:00:00	24.2
26. 6. 30. 오전 2:30:00	24.6
26. 6. 30. 오전 3:00:00	24.9
26. 6. 30. 오전 3:30:00	24.7
26. 6. 30. 오전 4:00:00	24.1
26. 6. 30. 오전 4:30:00	24.2
26. 6. 30. 오전 5:00:00	24.8
26. 6. 30. 오전 5:30:00	25.1
26. 6. 30. 오전 6:00:00	24.2
26. 6. 30. 오전 6:30:00	23.9
26. 6. 30. 오전 7:00:00	24.3
26. 6. 30. 오전 7:30:00	23.7
26. 6. 30. 오전 8:00:00	24.3
26. 6. 30. 오전 8:30:00	23.9
26. 6. 30. 오전 9:00:00	24.4
26. 6. 30. 오전 9:30:00	24.7
26. 6. 30. 오전 10:00:00	24.9
26. 6. 30. 오전 10:30:00	24.9
26. 6. 30. 오전 11:00:00	25.0
26. 6. 30. 오전 11:30:00	25.2
26. 6. 30. 오후 12:00:00	25.4
26. 6. 30. 오후 12:30:00	25.7
26. 6. 30. 오후 1:00:00	25.8
26. 6. 30. 오후 1:30:00	25.9
26. 6. 30. 오후 2:00:00	25.6
26. 6. 30. 오후 2:30:00	24.6
26. 6. 30. 오후 3:00:00	24.4
26. 6. 30. 오후 3:30:00	25.0
26. 6. 30. 오후 4:00:00	24.7
26. 6. 30. 오후 4:30:00	24.8
26. 6. 30. 오후 5:00:00	24.0
26. 6. 30. 오후 5:30:00	24.3



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30℃)_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 30. 오후 6:00:00	24.6
26. 6. 30. 오후 6:30:00	24.5
26. 6. 30. 오후 7:00:00	24.5
26. 6. 30. 오후 7:30:00	24.6
26. 6. 30. 오후 8:00:00	24.5
26. 6. 30. 오후 8:30:00	24.3
26. 6. 30. 오후 9:00:00	24.2
26. 6. 30. 오후 9:30:00	24.0
26. 6. 30. 오후 10:00:00	23.8
26. 6. 30. 오후 10:30:00	23.4
26. 6. 30. 오후 11:00:00	24.2
26. 6. 30. 오후 11:30:00	24.0



IP Storage Temperature Log_Monthly report
2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 55초 KST

알람 테이블

다음에서 발생 [Asia/Seoul]	마감일 [Asia/Seoul]	측정 포인트	알림 유형	상세 정보
26. 6. 16. 오전 5:00:00	26. 6. 16. 오전 5:30:08	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_CH2	온도	28.6 °C > 28.5 °C 상부 경고 값