

IP STORAGE TEMPERATURE LOG

Period of Record	Year. Month	2026.6	
Storage	Location	동국대학교일산병원 임상시험센터 약국	9134Room
	Description	전기냉장냉동기	항온고 1
	Manufacturer	SANYO	
	Model Name	MPR-720R	
	Serial Number	07120004	
	Temperature Range	15~25℃	
	Setting Point	19.0℃	
Thermometer (including alarm system)	Description	Wi-Fi data logger	
	Manufacturer	Testo Korea	
	Model Name	Saveris2-T3	
	Serial Number	0054704396	
	Sensor Channel	Channel 1	
Thermometer Calibration Certification	Min/Max Record	Monthly report 2page / This month Min (19.0)℃, Max(19.3)℃	
	Date of Calibration	2025.11.25	
	Date of Calibration Due	2026.11.25	
Alarm System	Certificate No.	A25-2356-001	
	Alarm threshold	Lower : 15.5℃, Upper : 24.5℃, immediate alarm without delay	
	Alarm Method	Mail, SMS to Pharmacists	

Date	Time	Temperature(℃)	Deviation	Verified by	Comment
2026.06.01	P:30	19.1	☐ Y / ☒ N	ml	
2026.06.02	P:30	19.1	☐ Y / ☒ N	ml	
2026.06.03			☐ Y / ☐ N		2026 지방선거
2026.06.04	10:00	19.1	☐ Y / ☒ N	ml	
2026.06.05	P:30	19.1	☐ Y / ☒ N	ml	
2026.06.06			☐ Y / ☐ N		
2026.06.07			☐ Y / ☐ N		
2026.06.08	10:00	19.2	☐ Y / ☒ N	ml	
2026.06.09	P:10	19.2	☐ Y / ☒ N	ml	
2026.06.10	P:20	19.2	☐ Y / ☒ N	ml	
2026.06.11	P:40	19.2	☐ Y / ☒ N	ml	
2026.06.12	P:40	19.1	☐ Y / ☒ N	ml	
2026.06.13			☐ Y / ☐ N		
2026.06.14			☐ Y / ☐ N		
2026.06.15	P:30	19.2	☐ Y / ☒ N	ml	
2026.06.16	P:35	19.2	☐ Y / ☒ N	ml	
2026.06.17	P:40	19.1	☐ Y / ☒ N	ml	
2026.06.18	P:20	19.2	☐ Y / ☒ N	ml	
2026.06.19	P:30	19.2	☐ Y / ☒ N	ml	
2026.06.20			☐ Y / ☐ N		
2026.06.21			☐ Y / ☐ N		
2026.06.22	P:30	19.1	☐ Y / ☒ N	ml	
2026.06.23	P:20	19.2	☐ Y / ☒ N	ml	
2026.06.24	P:15	19.1	☐ Y / ☒ N	ml	
2026.06.25	P:20	19.1	☐ Y / ☒ N	ml	
2026.06.26	P:10	19.1	☐ Y / ☒ N	ml	
2026.06.27			☐ Y / ☐ N		
2026.06.28			☐ Y / ☐ N		
2026.06.29	P:30	19.2	☐ Y / ☒ N	김성민	
2026.06.30	P:20	19.0	☐ Y / ☒ N	ml	
			☐ Y / ☐ N		

date of copying

sign



측정 테이블

측정 포인트 요약

동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
Min	19.0
Max	19.3
Average	19.1
26. 6. 1. 오전 12:00:00	19.1
26. 6. 1. 오전 12:30:00	19.1
26. 6. 1. 오전 1:00:00	19.2
26. 6. 1. 오전 1:30:00	19.1
26. 6. 1. 오전 2:00:00	19.1
26. 6. 1. 오전 2:30:00	19.2
26. 6. 1. 오전 3:00:00	19.2
26. 6. 1. 오전 3:30:00	19.2
26. 6. 1. 오전 4:00:00	19.1
26. 6. 1. 오전 4:30:00	19.2
26. 6. 1. 오전 5:00:00	19.2
26. 6. 1. 오전 5:30:00	19.1
26. 6. 1. 오전 6:00:00	19.2
26. 6. 1. 오전 6:30:00	19.2
26. 6. 1. 오전 7:00:00	19.2
26. 6. 1. 오전 7:30:00	19.1
26. 6. 1. 오전 8:00:00	19.1
26. 6. 1. 오전 8:30:00	19.2
26. 6. 1. 오전 9:00:00	19.2
26. 6. 1. 오전 9:30:00	19.1
26. 6. 1. 오전 10:00:00	19.2



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 1. 오전 10:30:00	19.2
26. 6. 1. 오전 11:00:00	19.2
26. 6. 1. 오전 11:30:00	19.1
26. 6. 1. 오후 12:00:00	19.2
26. 6. 1. 오후 12:30:00	19.2
26. 6. 1. 오후 1:00:00	19.1
26. 6. 1. 오후 1:30:00	19.1
26. 6. 1. 오후 2:00:00	19.2
26. 6. 1. 오후 2:30:00	19.1
26. 6. 1. 오후 3:00:00	19.1
26. 6. 1. 오후 3:30:00	19.2
26. 6. 1. 오후 4:00:00	19.2
26. 6. 1. 오후 4:30:00	19.1
26. 6. 1. 오후 5:00:00	19.1
26. 6. 1. 오후 5:30:00	19.2
26. 6. 1. 오후 6:00:00	19.1
26. 6. 1. 오후 6:30:00	19.1
26. 6. 1. 오후 7:00:00	19.1
26. 6. 1. 오후 7:30:00	19.2
26. 6. 1. 오후 8:00:00	19.2
26. 6. 1. 오후 8:30:00	19.1
26. 6. 1. 오후 9:00:00	19.2
26. 6. 1. 오후 9:30:00	19.2
26. 6. 1. 오후 10:00:00	19.1
26. 6. 1. 오후 10:30:00	19.1
26. 6. 1. 오후 11:00:00	19.2
26. 6. 1. 오후 11:30:00	19.2
26. 6. 2. 오전 12:00:00	19.1
26. 6. 2. 오전 12:30:00	19.2
26. 6. 2. 오전 1:00:00	19.2
26. 6. 2. 오전 1:30:00	19.2
26. 6. 2. 오전 2:00:00	19.2



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 2. 오전 2:30:00	19.2
26. 6. 2. 오전 3:00:00	19.2
26. 6. 2. 오전 3:30:00	19.2
26. 6. 2. 오전 4:00:00	19.2
26. 6. 2. 오전 4:30:00	19.1
26. 6. 2. 오전 5:00:00	19.2
26. 6. 2. 오전 5:30:00	19.2
26. 6. 2. 오전 6:00:00	19.2
26. 6. 2. 오전 6:30:00	19.1
26. 6. 2. 오전 7:00:00	19.1
26. 6. 2. 오전 7:30:00	19.2
26. 6. 2. 오전 8:00:00	19.2
26. 6. 2. 오전 8:30:00	19.2
26. 6. 2. 오전 9:00:00	19.1
26. 6. 2. 오전 9:30:00	19.1
26. 6. 2. 오전 10:00:00	19.2
26. 6. 2. 오전 10:30:00	19.2
26. 6. 2. 오전 11:00:00	19.1
26. 6. 2. 오전 11:30:00	19.1
26. 6. 2. 오후 12:00:00	19.1
26. 6. 2. 오후 12:30:00	19.2
26. 6. 2. 오후 1:00:00	19.2
26. 6. 2. 오후 1:30:00	19.1
26. 6. 2. 오후 2:00:00	19.1
26. 6. 2. 오후 2:30:00	19.2
26. 6. 2. 오후 3:00:00	19.2
26. 6. 2. 오후 3:30:00	19.1
26. 6. 2. 오후 4:00:00	19.1
26. 6. 2. 오후 4:30:00	19.1
26. 6. 2. 오후 5:00:00	19.2
26. 6. 2. 오후 5:30:00	19.3
26. 6. 2. 오후 6:00:00	19.2



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 2. 오후 6:30:00	19.2
26. 6. 2. 오후 7:00:00	19.1
26. 6. 2. 오후 7:30:00	19.1
26. 6. 2. 오후 8:00:00	19.2
26. 6. 2. 오후 8:30:00	19.2
26. 6. 2. 오후 9:00:00	19.2
26. 6. 2. 오후 9:30:00	19.2
26. 6. 2. 오후 10:00:00	19.2
26. 6. 2. 오후 10:30:00	19.1
26. 6. 2. 오후 11:00:00	19.1
26. 6. 2. 오후 11:30:00	19.1
26. 6. 3. 오전 12:00:00	19.1
26. 6. 3. 오전 12:30:00	19.1
26. 6. 3. 오전 1:00:00	19.2
26. 6. 3. 오전 1:30:00	19.2
26. 6. 3. 오전 2:00:00	19.2
26. 6. 3. 오전 2:30:00	19.2
26. 6. 3. 오전 3:00:00	19.2
26. 6. 3. 오전 3:30:00	19.1
26. 6. 3. 오전 4:00:00	19.1
26. 6. 3. 오전 4:30:00	19.1
26. 6. 3. 오전 5:00:00	19.1
26. 6. 3. 오전 5:30:00	19.1
26. 6. 3. 오전 6:00:00	19.1
26. 6. 3. 오전 6:30:00	19.1
26. 6. 3. 오전 7:00:00	19.1
26. 6. 3. 오전 7:30:00	19.2
26. 6. 3. 오전 8:00:00	19.2
26. 6. 3. 오전 8:30:00	19.2
26. 6. 3. 오전 9:00:00	19.1
26. 6. 3. 오전 9:30:00	19.2
26. 6. 3. 오전 10:00:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 3. 오전 10:30:00	19.2
26. 6. 3. 오전 11:00:00	19.2
26. 6. 3. 오전 11:30:00	19.2
26. 6. 3. 오후 12:00:00	19.1
26. 6. 3. 오후 12:30:00	19.1
26. 6. 3. 오후 1:00:00	19.1
26. 6. 3. 오후 1:30:00	19.1
26. 6. 3. 오후 2:00:00	19.1
26. 6. 3. 오후 2:30:00	19.1
26. 6. 3. 오후 3:00:00	19.2
26. 6. 3. 오후 3:30:00	19.2
26. 6. 3. 오후 4:00:00	19.2
26. 6. 3. 오후 4:30:00	19.2
26. 6. 3. 오후 5:00:00	19.2
26. 6. 3. 오후 5:30:00	19.1
26. 6. 3. 오후 6:00:00	19.1
26. 6. 3. 오후 6:30:00	19.1
26. 6. 3. 오후 7:00:00	19.1
26. 6. 3. 오후 7:30:00	19.2
26. 6. 3. 오후 8:00:00	19.2
26. 6. 3. 오후 8:30:00	19.2
26. 6. 3. 오후 9:00:00	19.3
26. 6. 3. 오후 9:30:00	19.2
26. 6. 3. 오후 10:00:00	19.1
26. 6. 3. 오후 10:30:00	19.1
26. 6. 3. 오후 11:00:00	19.1
26. 6. 3. 오후 11:30:00	19.1
26. 6. 4. 오전 12:00:00	19.1
26. 6. 4. 오전 12:30:00	19.1
26. 6. 4. 오전 1:00:00	19.1
26. 6. 4. 오전 1:30:00	19.1
26. 6. 4. 오전 2:00:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 4. 오전 2:30:00	19.1
26. 6. 4. 오전 3:00:00	19.1
26. 6. 4. 오전 3:30:00	19.1
26. 6. 4. 오전 4:00:00	19.2
26. 6. 4. 오전 4:30:00	19.2
26. 6. 4. 오전 5:00:00	19.1
26. 6. 4. 오전 5:30:00	19.1
26. 6. 4. 오전 6:00:00	19.2
26. 6. 4. 오전 6:30:00	19.3
26. 6. 4. 오전 7:00:00	19.3
26. 6. 4. 오전 7:30:00	19.2
26. 6. 4. 오전 8:00:00	19.2
26. 6. 4. 오전 8:30:00	19.1
26. 6. 4. 오전 9:00:00	19.1
26. 6. 4. 오전 9:30:00	19.1
26. 6. 4. 오전 10:00:00	19.2
26. 6. 4. 오전 10:30:00	19.2
26. 6. 4. 오전 11:00:00	19.2
26. 6. 4. 오전 11:30:00	19.2
26. 6. 4. 오후 12:00:00	19.1
26. 6. 4. 오후 12:30:00	19.1
26. 6. 4. 오후 1:00:00	19.2
26. 6. 4. 오후 1:30:00	19.2
26. 6. 4. 오후 2:00:00	19.2
26. 6. 4. 오후 2:30:00	19.2
26. 6. 4. 오후 3:00:00	19.2
26. 6. 4. 오후 3:30:00	19.2
26. 6. 4. 오후 4:00:00	19.2
26. 6. 4. 오후 4:30:00	19.1
26. 6. 4. 오후 5:00:00	19.2
26. 6. 4. 오후 5:30:00	19.1
26. 6. 4. 오후 6:00:00	19.2



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 4. 오후 6:30:00	19.2
26. 6. 4. 오후 7:00:00	19.2
26. 6. 4. 오후 7:30:00	19.2
26. 6. 4. 오후 8:00:00	19.2
26. 6. 4. 오후 8:30:00	19.2
26. 6. 4. 오후 9:00:00	19.1
26. 6. 4. 오후 9:30:00	19.2
26. 6. 4. 오후 10:00:00	19.1
26. 6. 4. 오후 10:30:00	19.2
26. 6. 4. 오후 11:00:00	19.1
26. 6. 4. 오후 11:30:00	19.1
26. 6. 5. 오전 12:00:00	19.1
26. 6. 5. 오전 12:30:00	19.1
26. 6. 5. 오전 1:00:00	19.1
26. 6. 5. 오전 1:30:00	19.1
26. 6. 5. 오전 2:00:00	19.1
26. 6. 5. 오전 2:30:00	19.1
26. 6. 5. 오전 3:00:00	19.1
26. 6. 5. 오전 3:30:00	19.2
26. 6. 5. 오전 4:00:00	19.2
26. 6. 5. 오전 4:30:00	19.2
26. 6. 5. 오전 5:00:00	19.1
26. 6. 5. 오전 5:30:00	19.2
26. 6. 5. 오전 6:00:00	19.1
26. 6. 5. 오전 6:30:00	19.2
26. 6. 5. 오전 7:00:00	19.3
26. 6. 5. 오전 7:30:00	19.3
26. 6. 5. 오전 8:00:00	19.3
26. 6. 5. 오전 8:30:00	19.2
26. 6. 5. 오전 9:00:00	19.2
26. 6. 5. 오전 9:30:00	19.2
26. 6. 5. 오전 10:00:00	19.2



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 5. 오전 10:30:00	19.1
26. 6. 5. 오전 11:00:00	19.1
26. 6. 5. 오전 11:30:00	19.2
26. 6. 5. 오후 12:00:00	19.2
26. 6. 5. 오후 12:30:00	19.2
26. 6. 5. 오후 1:00:00	19.2
26. 6. 5. 오후 1:30:00	19.1
26. 6. 5. 오후 2:00:00	19.1
26. 6. 5. 오후 2:30:00	19.1
26. 6. 5. 오후 3:00:00	19.2
26. 6. 5. 오후 3:30:00	19.1
26. 6. 5. 오후 4:00:00	19.1
26. 6. 5. 오후 4:30:00	19.1
26. 6. 5. 오후 5:00:00	19.1
26. 6. 5. 오후 5:30:00	19.1
26. 6. 5. 오후 6:00:00	19.1
26. 6. 5. 오후 6:30:00	19.1
26. 6. 5. 오후 7:00:00	19.1
26. 6. 5. 오후 7:30:00	19.1
26. 6. 5. 오후 8:00:00	19.2
26. 6. 5. 오후 8:30:00	19.2
26. 6. 5. 오후 9:00:00	19.2
26. 6. 5. 오후 9:30:00	19.2
26. 6. 5. 오후 10:00:00	19.2
26. 6. 5. 오후 10:30:00	19.1
26. 6. 5. 오후 11:00:00	19.1
26. 6. 5. 오후 11:30:00	19.1
26. 6. 6. 오전 12:00:00	19.1
26. 6. 6. 오전 12:30:00	19.1
26. 6. 6. 오전 1:00:00	19.0
26. 6. 6. 오전 1:30:00	19.0
26. 6. 6. 오전 2:00:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 6. 오전 2:30:00	19.1
26. 6. 6. 오전 3:00:00	19.1
26. 6. 6. 오전 3:30:00	19.1
26. 6. 6. 오전 4:00:00	19.1
26. 6. 6. 오전 4:30:00	19.1
26. 6. 6. 오전 5:00:00	19.2
26. 6. 6. 오전 5:30:00	19.2
26. 6. 6. 오전 6:00:00	19.2
26. 6. 6. 오전 6:30:00	19.2
26. 6. 6. 오전 7:00:00	19.2
26. 6. 6. 오전 7:30:00	19.2
26. 6. 6. 오전 8:00:00	19.2
26. 6. 6. 오전 8:30:00	19.2
26. 6. 6. 오전 9:00:00	19.2
26. 6. 6. 오전 9:30:00	19.2
26. 6. 6. 오전 10:00:00	19.2
26. 6. 6. 오전 10:30:00	19.2
26. 6. 6. 오전 11:00:00	19.1
26. 6. 6. 오전 11:30:00	19.1
26. 6. 6. 오후 12:00:00	19.1
26. 6. 6. 오후 12:30:00	19.1
26. 6. 6. 오후 1:00:00	19.2
26. 6. 6. 오후 1:30:00	19.2
26. 6. 6. 오후 2:00:00	19.2
26. 6. 6. 오후 2:30:00	19.1
26. 6. 6. 오후 3:00:00	19.2
26. 6. 6. 오후 3:30:00	19.1
26. 6. 6. 오후 4:00:00	19.2
26. 6. 6. 오후 4:30:00	19.2
26. 6. 6. 오후 5:00:00	19.2
26. 6. 6. 오후 5:30:00	19.2
26. 6. 6. 오후 6:00:00	19.2



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 6. 오후 6:30:00	19.2
26. 6. 6. 오후 7:00:00	19.2
26. 6. 6. 오후 7:30:00	19.2
26. 6. 6. 오후 8:00:00	19.2
26. 6. 6. 오후 8:30:00	19.2
26. 6. 6. 오후 9:00:00	19.2
26. 6. 6. 오후 9:30:00	19.1
26. 6. 6. 오후 10:00:00	19.1
26. 6. 6. 오후 10:30:00	19.2
26. 6. 6. 오후 11:00:00	19.2
26. 6. 6. 오후 11:30:00	19.1
26. 6. 7. 오전 12:00:00	19.1
26. 6. 7. 오전 12:30:00	19.0
26. 6. 7. 오전 1:00:00	19.1
26. 6. 7. 오전 1:30:00	19.1
26. 6. 7. 오전 2:00:00	19.0
26. 6. 7. 오전 2:30:00	19.1
26. 6. 7. 오전 3:00:00	19.1
26. 6. 7. 오전 3:30:00	19.1
26. 6. 7. 오전 4:00:00	19.1
26. 6. 7. 오전 4:30:00	19.1
26. 6. 7. 오전 5:00:00	19.1
26. 6. 7. 오전 5:30:00	19.2
26. 6. 7. 오전 6:00:00	19.1
26. 6. 7. 오전 6:30:00	19.1
26. 6. 7. 오전 7:00:00	19.3
26. 6. 7. 오전 7:30:00	19.2
26. 6. 7. 오전 8:00:00	19.2
26. 6. 7. 오전 8:30:00	19.2
26. 6. 7. 오전 9:00:00	19.2
26. 6. 7. 오전 9:30:00	19.1
26. 6. 7. 오전 10:00:00	19.2



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 7. 오전 10:30:00	19.2
26. 6. 7. 오전 11:00:00	19.1
26. 6. 7. 오전 11:30:00	19.1
26. 6. 7. 오후 12:00:00	19.1
26. 6. 7. 오후 12:30:00	19.2
26. 6. 7. 오후 1:00:00	19.1
26. 6. 7. 오후 1:30:00	19.1
26. 6. 7. 오후 2:00:00	19.1
26. 6. 7. 오후 2:30:00	19.2
26. 6. 7. 오후 3:00:00	19.1
26. 6. 7. 오후 3:30:00	19.1
26. 6. 7. 오후 4:00:00	19.1
26. 6. 7. 오후 4:30:00	19.1
26. 6. 7. 오후 5:00:00	19.2
26. 6. 7. 오후 5:30:00	19.1
26. 6. 7. 오후 6:00:00	19.1
26. 6. 7. 오후 6:30:00	19.1
26. 6. 7. 오후 7:00:00	19.2
26. 6. 7. 오후 7:30:00	19.2
26. 6. 7. 오후 8:00:00	19.2
26. 6. 7. 오후 8:30:00	19.1
26. 6. 7. 오후 9:00:00	19.1
26. 6. 7. 오후 9:30:00	19.2
26. 6. 7. 오후 10:00:00	19.2
26. 6. 7. 오후 10:30:00	19.2
26. 6. 7. 오후 11:00:00	19.2
26. 6. 7. 오후 11:30:00	19.2
26. 6. 8. 오전 12:00:00	19.1
26. 6. 8. 오전 12:30:00	19.0
26. 6. 8. 오전 1:00:00	19.1
26. 6. 8. 오전 1:30:00	19.1
26. 6. 8. 오전 2:00:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 8. 오전 2:30:00	19.1
26. 6. 8. 오전 3:00:00	19.1
26. 6. 8. 오전 3:30:00	19.1
26. 6. 8. 오전 4:00:00	19.1
26. 6. 8. 오전 4:30:00	19.2
26. 6. 8. 오전 5:00:00	19.2
26. 6. 8. 오전 5:30:00	19.1
26. 6. 8. 오전 6:00:00	19.1
26. 6. 8. 오전 6:30:00	19.2
26. 6. 8. 오전 7:00:00	19.3
26. 6. 8. 오전 7:30:00	19.3
26. 6. 8. 오전 8:00:00	19.2
26. 6. 8. 오전 8:30:00	19.2
26. 6. 8. 오전 9:00:00	19.1
26. 6. 8. 오전 9:30:00	19.1
26. 6. 8. 오전 10:00:00	19.2
26. 6. 8. 오전 10:30:00	19.2
26. 6. 8. 오전 11:00:00	19.2
26. 6. 8. 오전 11:30:00	19.2
26. 6. 8. 오후 12:00:00	19.1
26. 6. 8. 오후 12:30:00	19.1
26. 6. 8. 오후 1:00:00	19.2
26. 6. 8. 오후 1:30:00	19.2
26. 6. 8. 오후 2:00:00	19.2
26. 6. 8. 오후 2:30:00	19.1
26. 6. 8. 오후 3:00:00	19.1
26. 6. 8. 오후 3:30:00	19.2
26. 6. 8. 오후 4:00:00	19.2
26. 6. 8. 오후 4:30:00	19.2
26. 6. 8. 오후 5:00:00	19.1
26. 6. 8. 오후 5:30:00	19.2
26. 6. 8. 오후 6:00:00	19.2



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 8. 오후 6:30:00	19.2
26. 6. 8. 오후 7:00:00	19.2
26. 6. 8. 오후 7:30:00	19.1
26. 6. 8. 오후 8:00:00	19.2
26. 6. 8. 오후 8:30:00	19.2
26. 6. 8. 오후 9:00:00	19.2
26. 6. 8. 오후 9:30:00	19.1
26. 6. 8. 오후 10:00:00	19.1
26. 6. 8. 오후 10:30:00	19.1
26. 6. 8. 오후 11:00:00	19.2
26. 6. 8. 오후 11:30:00	19.2
26. 6. 9. 오전 12:00:00	19.1
26. 6. 9. 오전 12:30:00	19.1
26. 6. 9. 오전 1:00:00	19.1
26. 6. 9. 오전 1:30:00	19.1
26. 6. 9. 오전 2:00:00	19.1
26. 6. 9. 오전 2:30:00	19.1
26. 6. 9. 오전 3:00:00	19.1
26. 6. 9. 오전 3:30:00	19.2
26. 6. 9. 오전 4:00:00	19.2
26. 6. 9. 오전 4:30:00	19.1
26. 6. 9. 오전 5:00:00	19.1
26. 6. 9. 오전 5:30:00	19.2
26. 6. 9. 오전 6:00:00	19.1
26. 6. 9. 오전 6:30:00	19.2
26. 6. 9. 오전 7:00:00	19.2
26. 6. 9. 오전 7:30:00	19.2
26. 6. 9. 오전 8:00:00	19.2
26. 6. 9. 오전 8:30:00	19.1
26. 6. 9. 오전 9:00:00	19.1
26. 6. 9. 오전 9:30:00	19.2
26. 6. 9. 오전 10:00:00	19.2



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 9. 오전 10:30:00	19.1
26. 6. 9. 오전 11:00:00	19.1
26. 6. 9. 오전 11:30:00	19.2
26. 6. 9. 오후 12:00:00	19.1
26. 6. 9. 오후 12:30:00	19.1
26. 6. 9. 오후 1:00:00	19.1
26. 6. 9. 오후 1:30:00	19.2
26. 6. 9. 오후 2:00:00	19.1
26. 6. 9. 오후 2:30:00	19.2
26. 6. 9. 오후 3:00:00	19.2
26. 6. 9. 오후 3:30:00	19.1
26. 6. 9. 오후 4:00:00	19.1
26. 6. 9. 오후 4:30:00	19.2
26. 6. 9. 오후 5:00:00	19.2
26. 6. 9. 오후 5:30:00	19.1
26. 6. 9. 오후 6:00:00	19.2
26. 6. 9. 오후 6:30:00	19.2
26. 6. 9. 오후 7:00:00	19.2
26. 6. 9. 오후 7:30:00	19.1
26. 6. 9. 오후 8:00:00	19.2
26. 6. 9. 오후 8:30:00	19.2
26. 6. 9. 오후 9:00:00	19.2
26. 6. 9. 오후 9:30:00	19.2
26. 6. 9. 오후 10:00:00	19.2
26. 6. 9. 오후 10:30:00	19.2
26. 6. 9. 오후 11:00:00	19.1
26. 6. 9. 오후 11:30:00	19.2
26. 6. 10. 오전 12:00:00	19.2
26. 6. 10. 오전 12:30:00	19.0
26. 6. 10. 오전 1:00:00	19.1
26. 6. 10. 오전 1:30:00	19.1
26. 6. 10. 오전 2:00:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 10. 오전 2:30:00	19.1
26. 6. 10. 오전 3:00:00	19.1
26. 6. 10. 오전 3:30:00	19.2
26. 6. 10. 오전 4:00:00	19.1
26. 6. 10. 오전 4:30:00	19.1
26. 6. 10. 오전 5:00:00	19.2
26. 6. 10. 오전 5:30:00	19.2
26. 6. 10. 오전 6:00:00	19.2
26. 6. 10. 오전 6:30:00	19.3
26. 6. 10. 오전 7:00:00	19.3
26. 6. 10. 오전 7:30:00	19.2
26. 6. 10. 오전 8:00:00	19.2
26. 6. 10. 오전 8:30:00	19.2
26. 6. 10. 오전 9:00:00	19.2
26. 6. 10. 오전 9:30:00	19.2
26. 6. 10. 오전 10:00:00	19.2
26. 6. 10. 오전 10:30:00	19.1
26. 6. 10. 오전 11:00:00	19.2
26. 6. 10. 오전 11:30:00	19.3
26. 6. 10. 오후 12:00:00	19.2
26. 6. 10. 오후 12:30:00	19.2
26. 6. 10. 오후 1:00:00	19.2
26. 6. 10. 오후 1:30:00	19.2
26. 6. 10. 오후 2:00:00	19.2
26. 6. 10. 오후 2:30:00	19.2
26. 6. 10. 오후 3:00:00	19.2
26. 6. 10. 오후 3:30:00	19.1
26. 6. 10. 오후 4:00:00	19.1
26. 6. 10. 오후 4:30:00	19.1
26. 6. 10. 오후 5:00:00	19.2
26. 6. 10. 오후 5:30:00	19.2
26. 6. 10. 오후 6:00:00	19.2



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 10. 오후 6:30:00	19.1
26. 6. 10. 오후 7:00:00	19.1
26. 6. 10. 오후 7:30:00	19.1
26. 6. 10. 오후 8:00:00	19.1
26. 6. 10. 오후 8:30:00	19.1
26. 6. 10. 오후 9:00:00	19.2
26. 6. 10. 오후 9:30:00	19.2
26. 6. 10. 오후 10:00:00	19.1
26. 6. 10. 오후 10:30:00	19.1
26. 6. 10. 오후 11:00:00	19.1
26. 6. 10. 오후 11:30:00	19.1
26. 6. 11. 오전 12:00:00	19.1
26. 6. 11. 오전 12:30:00	19.1
26. 6. 11. 오전 1:00:00	19.1
26. 6. 11. 오전 1:30:00	19.1
26. 6. 11. 오전 2:00:00	19.1
26. 6. 11. 오전 2:30:00	19.1
26. 6. 11. 오전 3:00:00	19.0
26. 6. 11. 오전 3:30:00	19.1
26. 6. 11. 오전 4:00:00	19.1
26. 6. 11. 오전 4:30:00	19.2
26. 6. 11. 오전 5:00:00	19.1
26. 6. 11. 오전 5:30:00	19.1
26. 6. 11. 오전 6:00:00	19.1
26. 6. 11. 오전 6:30:00	19.2
26. 6. 11. 오전 7:00:00	19.2
26. 6. 11. 오전 7:30:00	19.2
26. 6. 11. 오전 8:00:00	19.2
26. 6. 11. 오전 8:30:00	19.2
26. 6. 11. 오전 9:00:00	19.1
26. 6. 11. 오전 9:30:00	19.1
26. 6. 11. 오전 10:00:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25℃) CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 11. 오전 10:30:00	19.1
26. 6. 11. 오전 11:00:00	19.1
26. 6. 11. 오전 11:30:00	19.1
26. 6. 11. 오후 12:00:00	19.1
26. 6. 11. 오후 12:30:00	19.1
26. 6. 11. 오후 1:00:00	19.2
26. 6. 11. 오후 1:30:00	19.1
26. 6. 11. 오후 2:00:00	19.1
26. 6. 11. 오후 2:30:00	19.1
26. 6. 11. 오후 3:00:00	19.1
26. 6. 11. 오후 3:30:00	19.1
26. 6. 11. 오후 4:00:00	19.2
26. 6. 11. 오후 4:30:00	19.2
26. 6. 11. 오후 5:00:00	19.1
26. 6. 11. 오후 5:30:00	19.1
26. 6. 11. 오후 6:00:00	19.1
26. 6. 11. 오후 6:30:00	19.1
26. 6. 11. 오후 7:00:00	19.2
26. 6. 11. 오후 7:30:00	19.1
26. 6. 11. 오후 8:00:00	19.2
26. 6. 11. 오후 8:30:00	19.2
26. 6. 11. 오후 9:00:00	19.2
26. 6. 11. 오후 9:30:00	19.1
26. 6. 11. 오후 10:00:00	19.2
26. 6. 11. 오후 10:30:00	19.1
26. 6. 11. 오후 11:00:00	19.2
26. 6. 11. 오후 11:30:00	19.1
26. 6. 12. 오전 12:00:00	19.1
26. 6. 12. 오전 12:30:00	19.0
26. 6. 12. 오전 1:00:00	19.0
26. 6. 12. 오전 1:30:00	19.0
26. 6. 12. 오전 2:00:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 12. 오전 2:30:00	19.2
26. 6. 12. 오전 3:00:00	19.2
26. 6. 12. 오전 3:30:00	19.2
26. 6. 12. 오전 4:00:00	19.1
26. 6. 12. 오전 4:30:00	19.1
26. 6. 12. 오전 5:00:00	19.1
26. 6. 12. 오전 5:30:00	19.1
26. 6. 12. 오전 6:00:00	19.2
26. 6. 12. 오전 6:30:00	19.2
26. 6. 12. 오전 7:00:00	19.2
26. 6. 12. 오전 7:30:00	19.2
26. 6. 12. 오전 8:00:00	19.1
26. 6. 12. 오전 8:30:00	19.1
26. 6. 12. 오전 9:00:00	19.1
26. 6. 12. 오전 9:30:00	19.1
26. 6. 12. 오전 10:00:00	19.1
26. 6. 12. 오전 10:30:00	19.1
26. 6. 12. 오전 11:00:00	19.1
26. 6. 12. 오전 11:30:00	19.1
26. 6. 12. 오후 12:00:00	19.1
26. 6. 12. 오후 12:30:00	19.1
26. 6. 12. 오후 1:00:00	19.1
26. 6. 12. 오후 1:30:00	19.1
26. 6. 12. 오후 2:00:00	19.1
26. 6. 12. 오후 2:30:00	19.1
26. 6. 12. 오후 2:50:30	19.3
26. 6. 12. 오후 3:00:00	19.1
26. 6. 12. 오후 3:30:00	19.1
26. 6. 12. 오후 4:00:00	19.1
26. 6. 12. 오후 4:30:00	19.1
26. 6. 12. 오후 5:00:00	19.1
26. 6. 12. 오후 5:30:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 12. 오후 6:00:00	19.1
26. 6. 12. 오후 6:30:00	19.1
26. 6. 12. 오후 7:00:00	19.2
26. 6. 12. 오후 7:30:00	19.1
26. 6. 12. 오후 8:00:00	19.1
26. 6. 12. 오후 8:30:00	19.1
26. 6. 12. 오후 9:00:00	19.1
26. 6. 12. 오후 9:30:00	19.1
26. 6. 12. 오후 10:00:00	19.2
26. 6. 12. 오후 10:30:00	19.2
26. 6. 12. 오후 11:00:00	19.1
26. 6. 12. 오후 11:30:00	19.1
26. 6. 13. 오전 12:00:00	19.0
26. 6. 13. 오전 12:30:00	19.0
26. 6. 13. 오전 1:00:00	19.0
26. 6. 13. 오전 1:30:00	19.0
26. 6. 13. 오전 2:00:00	19.1
26. 6. 13. 오전 2:30:00	19.1
26. 6. 13. 오전 3:00:00	19.1
26. 6. 13. 오전 3:30:00	19.1
26. 6. 13. 오전 4:00:00	19.1
26. 6. 13. 오전 4:30:00	19.1
26. 6. 13. 오전 5:00:00	19.1
26. 6. 13. 오전 5:30:00	19.1
26. 6. 13. 오전 6:00:00	19.1
26. 6. 13. 오전 6:30:00	19.1
26. 6. 13. 오전 7:00:00	19.2
26. 6. 13. 오전 7:30:00	19.2
26. 6. 13. 오전 8:00:00	19.2
26. 6. 13. 오전 8:30:00	19.2
26. 6. 13. 오전 9:00:00	19.1
26. 6. 13. 오전 9:30:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25℃)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 13. 오전 10:00:00	19.1
26. 6. 13. 오전 10:30:00	19.1
26. 6. 13. 오전 11:00:00	19.1
26. 6. 13. 오전 11:30:00	19.1
26. 6. 13. 오후 12:00:00	19.1
26. 6. 13. 오후 12:30:00	19.1
26. 6. 13. 오후 1:00:00	19.1
26. 6. 13. 오후 1:30:00	19.1
26. 6. 13. 오후 2:00:00	19.1
26. 6. 13. 오후 2:30:00	19.2
26. 6. 13. 오후 3:00:00	19.1
26. 6. 13. 오후 3:30:00	19.1
26. 6. 13. 오후 4:00:00	19.1
26. 6. 13. 오후 4:30:00	19.1
26. 6. 13. 오후 5:00:00	19.2
26. 6. 13. 오후 5:30:00	19.1
26. 6. 13. 오후 6:00:00	19.1
26. 6. 13. 오후 6:30:00	19.1
26. 6. 13. 오후 7:00:00	19.1
26. 6. 13. 오후 7:30:00	19.1
26. 6. 13. 오후 8:00:00	19.2
26. 6. 13. 오후 8:30:00	19.1
26. 6. 13. 오후 9:00:00	19.1
26. 6. 13. 오후 9:30:00	19.1
26. 6. 13. 오후 10:00:00	19.1
26. 6. 13. 오후 10:30:00	19.2
26. 6. 13. 오후 11:00:00	19.2
26. 6. 13. 오후 11:30:00	19.2
26. 6. 14. 오전 12:00:00	19.1
26. 6. 14. 오전 12:30:00	19.1
26. 6. 14. 오전 1:00:00	19.1
26. 6. 14. 오전 1:30:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 14. 오전 2:00:00	19.1
26. 6. 14. 오전 2:30:00	19.1
26. 6. 14. 오전 3:00:00	19.1
26. 6. 14. 오전 3:30:00	19.2
26. 6. 14. 오전 4:00:00	19.2
26. 6. 14. 오전 4:30:00	19.2
26. 6. 14. 오전 5:00:00	19.1
26. 6. 14. 오전 5:30:00	19.1
26. 6. 14. 오전 6:00:00	19.2
26. 6. 14. 오전 6:30:00	19.3
26. 6. 14. 오전 7:00:00	19.2
26. 6. 14. 오전 7:30:00	19.2
26. 6. 14. 오전 8:00:00	19.2
26. 6. 14. 오전 8:30:00	19.2
26. 6. 14. 오전 9:00:00	19.2
26. 6. 14. 오전 9:30:00	19.3
26. 6. 14. 오전 10:00:00	19.2
26. 6. 14. 오전 10:30:00	19.2
26. 6. 14. 오전 11:00:00	19.1
26. 6. 14. 오전 11:30:00	19.2
26. 6. 14. 오후 12:00:00	19.2
26. 6. 14. 오후 12:30:00	19.2
26. 6. 14. 오후 1:00:00	19.2
26. 6. 14. 오후 1:30:00	19.2
26. 6. 14. 오후 2:00:00	19.2
26. 6. 14. 오후 2:30:00	19.2
26. 6. 14. 오후 3:00:00	19.2
26. 6. 14. 오후 3:30:00	19.2
26. 6. 14. 오후 4:00:00	19.2
26. 6. 14. 오후 4:30:00	19.1
26. 6. 14. 오후 5:00:00	19.2
26. 6. 14. 오후 5:30:00	19.2



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 14. 오후 6:00:00	19.2
26. 6. 14. 오후 6:30:00	19.2
26. 6. 14. 오후 7:00:00	19.2
26. 6. 14. 오후 7:30:00	19.2
26. 6. 14. 오후 8:00:00	19.2
26. 6. 14. 오후 8:30:00	19.2
26. 6. 14. 오후 9:00:00	19.2
26. 6. 14. 오후 9:30:00	19.2
26. 6. 14. 오후 10:00:00	19.1
26. 6. 14. 오후 10:30:00	19.1
26. 6. 14. 오후 11:00:00	19.1
26. 6. 14. 오후 11:30:00	19.2
26. 6. 15. 오전 12:00:00	19.2
26. 6. 15. 오전 12:30:00	19.1
26. 6. 15. 오전 1:00:00	19.1
26. 6. 15. 오전 1:30:00	19.1
26. 6. 15. 오전 2:00:00	19.1
26. 6. 15. 오전 2:30:00	19.1
26. 6. 15. 오전 3:00:00	19.1
26. 6. 15. 오전 3:30:00	19.1
26. 6. 15. 오전 4:00:00	19.2
26. 6. 15. 오전 4:30:00	19.2
26. 6. 15. 오전 5:00:00	19.2
26. 6. 15. 오전 5:30:00	19.1
26. 6. 15. 오전 6:00:00	19.1
26. 6. 15. 오전 6:30:00	19.2
26. 6. 15. 오전 7:00:00	19.3
26. 6. 15. 오전 7:30:00	19.3
26. 6. 15. 오전 8:00:00	19.2
26. 6. 15. 오전 8:30:00	19.2
26. 6. 15. 오전 9:00:00	19.2
26. 6. 15. 오전 9:30:00	19.2



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 15. 오전 10:00:00	19.2
26. 6. 15. 오전 10:30:00	19.2
26. 6. 15. 오전 11:00:00	19.2
26. 6. 15. 오전 11:30:00	19.2
26. 6. 15. 오후 12:00:00	19.2
26. 6. 15. 오후 12:30:00	19.1
26. 6. 15. 오후 1:00:00	19.2
26. 6. 15. 오후 1:30:00	19.1
26. 6. 15. 오후 2:00:00	19.2
26. 6. 15. 오후 2:30:00	19.3
26. 6. 15. 오후 3:00:00	19.3
26. 6. 15. 오후 3:30:00	19.2
26. 6. 15. 오후 4:00:00	19.2
26. 6. 15. 오후 4:30:00	19.2
26. 6. 15. 오후 5:00:00	19.2
26. 6. 15. 오후 5:30:00	19.2
26. 6. 15. 오후 6:00:00	19.3
26. 6. 15. 오후 6:30:00	19.2
26. 6. 15. 오후 7:00:00	19.2
26. 6. 15. 오후 7:30:00	19.2
26. 6. 15. 오후 8:00:00	19.2
26. 6. 15. 오후 8:30:00	19.2
26. 6. 15. 오후 9:00:00	19.2
26. 6. 15. 오후 9:30:00	19.2
26. 6. 15. 오후 10:00:00	19.1
26. 6. 15. 오후 10:30:00	19.1
26. 6. 15. 오후 11:00:00	19.1
26. 6. 15. 오후 11:30:00	19.2
26. 6. 16. 오전 12:00:00	19.1
26. 6. 16. 오전 12:30:00	19.1
26. 6. 16. 오전 1:00:00	19.1
26. 6. 16. 오전 1:30:00	19.0



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 16. 오전 2:00:00	19.1
26. 6. 16. 오전 2:30:00	19.1
26. 6. 16. 오전 3:00:00	19.1
26. 6. 16. 오전 3:30:00	19.1
26. 6. 16. 오전 4:00:00	19.1
26. 6. 16. 오전 4:30:00	19.1
26. 6. 16. 오전 5:00:00	19.2
26. 6. 16. 오전 5:30:00	19.1
26. 6. 16. 오전 6:00:00	19.1
26. 6. 16. 오전 6:30:00	19.3
26. 6. 16. 오전 7:00:00	19.2
26. 6. 16. 오전 7:30:00	19.2
26. 6. 16. 오전 8:00:00	19.2
26. 6. 16. 오전 8:30:00	19.2
26. 6. 16. 오전 9:00:00	19.2
26. 6. 16. 오전 9:30:00	19.2
26. 6. 16. 오전 10:00:00	19.2
26. 6. 16. 오전 10:30:00	19.1
26. 6. 16. 오전 11:00:00	19.2
26. 6. 16. 오전 11:30:00	19.2
26. 6. 16. 오후 12:00:00	19.2
26. 6. 16. 오후 12:30:00	19.2
26. 6. 16. 오후 1:00:00	19.2
26. 6. 16. 오후 1:30:00	19.1
26. 6. 16. 오후 2:00:00	19.2
26. 6. 16. 오후 2:30:00	19.2
26. 6. 16. 오후 3:00:00	19.2
26. 6. 16. 오후 3:30:00	19.2
26. 6. 16. 오후 4:00:00	19.2
26. 6. 16. 오후 4:30:00	19.2
26. 6. 16. 오후 5:00:00	19.1
26. 6. 16. 오후 5:30:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 16. 오후 6:00:00	19.1
26. 6. 16. 오후 6:30:00	19.1
26. 6. 16. 오후 7:00:00	19.1
26. 6. 16. 오후 7:30:00	19.1
26. 6. 16. 오후 8:00:00	19.1
26. 6. 16. 오후 8:30:00	19.1
26. 6. 16. 오후 9:00:00	19.1
26. 6. 16. 오후 9:30:00	19.2
26. 6. 16. 오후 10:00:00	19.2
26. 6. 16. 오후 10:30:00	19.1
26. 6. 16. 오후 11:00:00	19.1
26. 6. 16. 오후 11:30:00	19.2
26. 6. 17. 오전 12:00:00	19.1
26. 6. 17. 오전 12:30:00	19.1
26. 6. 17. 오전 1:00:00	19.1
26. 6. 17. 오전 1:30:00	19.1
26. 6. 17. 오전 2:00:00	19.1
26. 6. 17. 오전 2:30:00	19.1
26. 6. 17. 오전 3:00:00	19.1
26. 6. 17. 오전 3:30:00	19.1
26. 6. 17. 오전 4:00:00	19.1
26. 6. 17. 오전 4:30:00	19.1
26. 6. 17. 오전 5:00:00	19.2
26. 6. 17. 오전 5:30:00	19.2
26. 6. 17. 오전 6:00:00	19.2
26. 6. 17. 오전 6:30:00	19.2
26. 6. 17. 오전 7:00:00	19.1
26. 6. 17. 오전 7:30:00	19.1
26. 6. 17. 오전 8:00:00	19.1
26. 6. 17. 오전 8:30:00	19.1
26. 6. 17. 오전 9:00:00	19.1
26. 6. 17. 오전 9:30:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 17. 오전 10:00:00	19.1
26. 6. 17. 오전 10:30:00	19.1
26. 6. 17. 오전 11:00:00	19.2
26. 6. 17. 오전 11:30:00	19.1
26. 6. 17. 오후 12:00:00	19.1
26. 6. 17. 오후 12:30:00	19.2
26. 6. 17. 오후 1:00:00	19.1
26. 6. 17. 오후 1:30:00	19.1
26. 6. 17. 오후 2:00:00	19.1
26. 6. 17. 오후 2:30:00	19.1
26. 6. 17. 오후 3:00:00	19.1
26. 6. 17. 오후 3:30:00	19.1
26. 6. 17. 오후 4:00:00	19.1
26. 6. 17. 오후 4:30:00	19.1
26. 6. 17. 오후 5:00:00	19.2
26. 6. 17. 오후 5:30:00	19.1
26. 6. 17. 오후 6:00:00	19.1
26. 6. 17. 오후 6:30:00	19.2
26. 6. 17. 오후 7:00:00	19.1
26. 6. 17. 오후 7:30:00	19.2
26. 6. 17. 오후 8:00:00	19.2
26. 6. 17. 오후 8:30:00	19.2
26. 6. 17. 오후 9:00:00	19.1
26. 6. 17. 오후 9:30:00	19.1
26. 6. 17. 오후 10:00:00	19.1
26. 6. 17. 오후 10:30:00	19.1
26. 6. 17. 오후 11:00:00	19.1
26. 6. 17. 오후 11:30:00	19.1
26. 6. 18. 오전 12:00:00	19.1
26. 6. 18. 오전 12:30:00	19.1
26. 6. 18. 오전 1:00:00	19.1
26. 6. 18. 오전 1:30:00	19.2



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 18. 오전 2:00:00	19.2
26. 6. 18. 오전 2:30:00	19.2
26. 6. 18. 오전 3:00:00	19.1
26. 6. 18. 오전 3:30:00	19.1
26. 6. 18. 오전 4:00:00	19.1
26. 6. 18. 오전 4:30:00	19.1
26. 6. 18. 오전 5:00:00	19.1
26. 6. 18. 오전 5:30:00	19.1
26. 6. 18. 오전 6:00:00	19.1
26. 6. 18. 오전 6:30:00	19.2
26. 6. 18. 오전 7:00:00	19.2
26. 6. 18. 오전 7:30:00	19.2
26. 6. 18. 오전 8:00:00	19.2
26. 6. 18. 오전 8:30:00	19.2
26. 6. 18. 오전 9:00:00	19.2
26. 6. 18. 오전 9:30:00	19.2
26. 6. 18. 오전 10:00:00	19.1
26. 6. 18. 오전 10:30:00	19.1
26. 6. 18. 오전 11:00:00	19.1
26. 6. 18. 오전 11:30:00	19.1
26. 6. 18. 오후 12:00:00	19.1
26. 6. 18. 오후 12:30:00	19.1
26. 6. 18. 오후 1:00:00	19.1
26. 6. 18. 오후 1:30:00	19.1
26. 6. 18. 오후 2:00:00	19.1
26. 6. 18. 오후 2:30:00	19.1
26. 6. 18. 오후 3:00:00	19.1
26. 6. 18. 오후 3:30:00	19.1
26. 6. 18. 오후 4:00:00	19.2
26. 6. 18. 오후 4:30:00	19.2
26. 6. 18. 오후 5:00:00	19.1
26. 6. 18. 오후 5:30:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 18. 오후 6:00:00	19.1
26. 6. 18. 오후 6:30:00	19.1
26. 6. 18. 오후 7:00:00	19.1
26. 6. 18. 오후 7:30:00	19.1
26. 6. 18. 오후 8:00:00	19.1
26. 6. 18. 오후 8:30:00	19.1
26. 6. 18. 오후 9:00:00	19.2
26. 6. 18. 오후 9:30:00	19.2
26. 6. 18. 오후 10:00:00	19.2
26. 6. 18. 오후 10:30:00	19.2
26. 6. 18. 오후 11:00:00	19.2
26. 6. 18. 오후 11:30:00	19.1
26. 6. 19. 오전 12:00:00	19.1
26. 6. 19. 오전 12:30:00	19.1
26. 6. 19. 오전 1:00:00	19.1
26. 6. 19. 오전 1:30:00	19.1
26. 6. 19. 오전 2:00:00	19.1
26. 6. 19. 오전 2:30:00	19.1
26. 6. 19. 오전 3:00:00	19.1
26. 6. 19. 오전 3:30:00	19.1
26. 6. 19. 오전 4:00:00	19.1
26. 6. 19. 오전 4:30:00	19.2
26. 6. 19. 오전 5:00:00	19.2
26. 6. 19. 오전 5:30:00	19.1
26. 6. 19. 오전 6:00:00	19.2
26. 6. 19. 오전 6:30:00	19.1
26. 6. 19. 오전 7:00:00	19.2
26. 6. 19. 오전 7:30:00	19.2
26. 6. 19. 오전 8:00:00	19.2
26. 6. 19. 오전 8:30:00	19.2
26. 6. 19. 오전 9:00:00	19.2
26. 6. 19. 오전 9:30:00	19.2



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 19. 오전 10:00:00	19.2
26. 6. 19. 오전 10:30:00	19.2
26. 6. 19. 오전 11:00:00	19.1
26. 6. 19. 오전 11:30:00	19.1
26. 6. 19. 오후 12:00:00	19.1
26. 6. 19. 오후 12:30:00	19.1
26. 6. 19. 오후 1:00:00	19.1
26. 6. 19. 오후 1:30:00	19.1
26. 6. 19. 오후 2:00:00	19.1
26. 6. 19. 오후 2:30:00	19.2
26. 6. 19. 오후 3:00:00	19.2
26. 6. 19. 오후 3:30:00	19.2
26. 6. 19. 오후 4:00:00	19.2
26. 6. 19. 오후 4:30:00	19.2
26. 6. 19. 오후 5:00:00	19.1
26. 6. 19. 오후 5:30:00	19.1
26. 6. 19. 오후 6:00:00	19.1
26. 6. 19. 오후 6:30:00	19.1
26. 6. 19. 오후 7:00:00	19.1
26. 6. 19. 오후 7:30:00	19.1
26. 6. 19. 오후 8:00:00	19.1
26. 6. 19. 오후 8:30:00	19.1
26. 6. 19. 오후 9:00:00	19.1
26. 6. 19. 오후 9:30:00	19.2
26. 6. 19. 오후 10:00:00	19.2
26. 6. 19. 오후 10:30:00	19.1
26. 6. 19. 오후 11:00:00	19.1
26. 6. 19. 오후 11:30:00	19.1
26. 6. 20. 오전 12:00:00	19.1
26. 6. 20. 오전 12:30:00	19.1
26. 6. 20. 오전 1:00:00	19.1
26. 6. 20. 오전 1:30:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 20. 오전 2:00:00	19.2
26. 6. 20. 오전 2:30:00	19.2
26. 6. 20. 오전 3:00:00	19.2
26. 6. 20. 오전 3:30:00	19.2
26. 6. 20. 오전 4:00:00	19.1
26. 6. 20. 오전 4:30:00	19.1
26. 6. 20. 오전 5:00:00	19.1
26. 6. 20. 오전 5:30:00	19.1
26. 6. 20. 오전 6:00:00	19.1
26. 6. 20. 오전 6:30:00	19.1
26. 6. 20. 오전 7:00:00	19.1
26. 6. 20. 오전 7:30:00	19.1
26. 6. 20. 오전 8:00:00	19.1
26. 6. 20. 오전 8:30:00	19.1
26. 6. 20. 오전 9:00:00	19.1
26. 6. 20. 오전 9:30:00	19.1
26. 6. 20. 오전 10:00:00	19.2
26. 6. 20. 오전 10:30:00	19.2
26. 6. 20. 오전 11:00:00	19.2
26. 6. 20. 오전 11:30:00	19.1
26. 6. 20. 오후 12:00:00	19.2
26. 6. 20. 오후 12:30:00	19.2
26. 6. 20. 오후 1:00:00	19.2
26. 6. 20. 오후 1:30:00	19.1
26. 6. 20. 오후 2:00:00	19.2
26. 6. 20. 오후 2:30:00	19.1
26. 6. 20. 오후 3:00:00	19.1
26. 6. 20. 오후 3:30:00	19.1
26. 6. 20. 오후 4:00:00	19.1
26. 6. 20. 오후 4:30:00	19.1
26. 6. 20. 오후 5:00:00	19.1
26. 6. 20. 오후 5:30:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 20. 오후 6:00:00	19.1
26. 6. 20. 오후 6:30:00	19.1
26. 6. 20. 오후 7:00:00	19.1
26. 6. 20. 오후 7:30:00	19.1
26. 6. 20. 오후 8:00:00	19.1
26. 6. 20. 오후 8:30:00	19.1
26. 6. 20. 오후 9:00:00	19.1
26. 6. 20. 오후 9:30:00	19.1
26. 6. 20. 오후 10:00:00	19.1
26. 6. 20. 오후 10:30:00	19.1
26. 6. 20. 오후 11:00:00	19.1
26. 6. 20. 오후 11:30:00	19.2
26. 6. 21. 오전 12:00:00	19.2
26. 6. 21. 오전 12:30:00	19.1
26. 6. 21. 오전 1:00:00	19.2
26. 6. 21. 오전 1:30:00	19.2
26. 6. 21. 오전 2:00:00	19.2
26. 6. 21. 오전 2:30:00	19.2
26. 6. 21. 오전 3:00:00	19.1
26. 6. 21. 오전 3:30:00	19.1
26. 6. 21. 오전 4:00:00	19.1
26. 6. 21. 오전 4:30:00	19.1
26. 6. 21. 오전 5:00:00	19.1
26. 6. 21. 오전 5:30:00	19.1
26. 6. 21. 오전 6:00:00	19.1
26. 6. 21. 오전 6:30:00	19.1
26. 6. 21. 오전 7:00:00	19.1
26. 6. 21. 오전 7:30:00	19.1
26. 6. 21. 오전 8:00:00	19.1
26. 6. 21. 오전 8:30:00	19.1
26. 6. 21. 오전 9:00:00	19.2
26. 6. 21. 오전 9:30:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 21. 오전 10:00:00	19.1
26. 6. 21. 오전 10:30:00	19.1
26. 6. 21. 오전 11:00:00	19.1
26. 6. 21. 오전 11:30:00	19.1
26. 6. 21. 오후 12:00:00	19.2
26. 6. 21. 오후 12:30:00	19.2
26. 6. 21. 오후 1:00:00	19.1
26. 6. 21. 오후 1:30:00	19.2
26. 6. 21. 오후 2:00:00	19.2
26. 6. 21. 오후 2:30:00	19.2
26. 6. 21. 오후 3:00:00	19.2
26. 6. 21. 오후 3:30:00	19.2
26. 6. 21. 오후 4:00:00	19.2
26. 6. 21. 오후 4:30:00	19.1
26. 6. 21. 오후 5:00:00	19.1
26. 6. 21. 오후 5:30:00	19.1
26. 6. 21. 오후 6:00:00	19.1
26. 6. 21. 오후 6:30:00	19.1
26. 6. 21. 오후 7:00:00	19.1
26. 6. 21. 오후 7:30:00	19.1
26. 6. 21. 오후 8:00:00	19.1
26. 6. 21. 오후 8:30:00	19.1
26. 6. 21. 오후 9:00:00	19.1
26. 6. 21. 오후 9:30:00	19.1
26. 6. 21. 오후 10:00:00	19.1
26. 6. 21. 오후 10:30:00	19.1
26. 6. 21. 오후 11:00:00	19.1
26. 6. 21. 오후 11:30:00	19.1
26. 6. 22. 오전 12:00:00	19.1
26. 6. 22. 오전 12:30:00	19.1
26. 6. 22. 오전 1:00:00	19.1
26. 6. 22. 오전 1:30:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 22. 오전 2:00:00	19.1
26. 6. 22. 오전 2:30:00	19.2
26. 6. 22. 오전 3:00:00	19.1
26. 6. 22. 오전 3:30:00	19.1
26. 6. 22. 오전 4:00:00	19.1
26. 6. 22. 오전 4:30:00	19.1
26. 6. 22. 오전 5:00:00	19.1
26. 6. 22. 오전 5:30:00	19.1
26. 6. 22. 오전 6:00:00	19.2
26. 6. 22. 오전 6:30:00	19.2
26. 6. 22. 오전 7:00:00	19.2
26. 6. 22. 오전 7:30:00	19.2
26. 6. 22. 오전 8:00:00	19.2
26. 6. 22. 오전 8:30:00	19.1
26. 6. 22. 오전 9:00:00	19.1
26. 6. 22. 오전 9:30:00	19.1
26. 6. 22. 오전 10:00:00	19.1
26. 6. 22. 오전 10:30:00	19.1
26. 6. 22. 오전 11:00:00	19.1
26. 6. 22. 오전 11:30:00	19.1
26. 6. 22. 오후 12:00:00	19.1
26. 6. 22. 오후 12:30:00	19.1
26. 6. 22. 오후 1:00:00	19.1
26. 6. 22. 오후 1:30:00	19.2
26. 6. 22. 오후 2:00:00	19.2
26. 6. 22. 오후 2:30:00	19.2
26. 6. 22. 오후 3:00:00	19.1
26. 6. 22. 오후 3:30:00	19.2
26. 6. 22. 오후 4:00:00	19.2
26. 6. 22. 오후 4:30:00	19.2
26. 6. 22. 오후 5:00:00	19.1
26. 6. 22. 오후 5:30:00	19.2



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 22. 오후 6:00:00	19.1
26. 6. 22. 오후 6:30:00	19.1
26. 6. 22. 오후 7:00:00	19.1
26. 6. 22. 오후 7:30:00	19.1
26. 6. 22. 오후 8:00:00	19.1
26. 6. 22. 오후 8:30:00	19.1
26. 6. 22. 오후 9:00:00	19.1
26. 6. 22. 오후 9:30:00	19.1
26. 6. 22. 오후 10:00:00	19.1
26. 6. 22. 오후 10:30:00	19.1
26. 6. 22. 오후 11:00:00	19.1
26. 6. 22. 오후 11:30:00	19.1
26. 6. 23. 오전 12:00:00	19.1
26. 6. 23. 오전 12:30:00	19.1
26. 6. 23. 오전 1:00:00	19.1
26. 6. 23. 오전 1:30:00	19.1
26. 6. 23. 오전 2:00:00	19.2
26. 6. 23. 오전 2:30:00	19.2
26. 6. 23. 오전 3:00:00	19.1
26. 6. 23. 오전 3:30:00	19.1
26. 6. 23. 오전 4:00:00	19.1
26. 6. 23. 오전 4:30:00	19.1
26. 6. 23. 오전 5:00:00	19.1
26. 6. 23. 오전 5:30:00	19.1
26. 6. 23. 오전 6:00:00	19.1
26. 6. 23. 오전 6:30:00	19.2
26. 6. 23. 오전 7:00:00	19.2
26. 6. 23. 오전 7:30:00	19.2
26. 6. 23. 오전 8:00:00	19.2
26. 6. 23. 오전 8:30:00	19.2
26. 6. 23. 오전 9:00:00	19.2
26. 6. 23. 오전 9:30:00	19.2



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 23. 오전 10:00:00	19.2
26. 6. 23. 오전 10:30:00	19.1
26. 6. 23. 오전 11:00:00	19.1
26. 6. 23. 오전 11:30:00	19.1
26. 6. 23. 오후 12:00:00	19.1
26. 6. 23. 오후 12:30:00	19.1
26. 6. 23. 오후 1:00:00	19.1
26. 6. 23. 오후 1:30:00	19.1
26. 6. 23. 오후 2:00:00	19.1
26. 6. 23. 오후 2:30:00	19.1
26. 6. 23. 오후 3:00:00	19.1
26. 6. 23. 오후 3:30:00	19.1
26. 6. 23. 오후 4:00:00	19.1
26. 6. 23. 오후 4:30:00	19.1
26. 6. 23. 오후 5:00:00	19.1
26. 6. 23. 오후 5:30:00	19.1
26. 6. 23. 오후 6:00:00	19.1
26. 6. 23. 오후 6:30:00	19.1
26. 6. 23. 오후 7:00:00	19.1
26. 6. 23. 오후 7:30:00	19.1
26. 6. 23. 오후 8:00:00	19.2
26. 6. 23. 오후 8:30:00	19.1
26. 6. 23. 오후 9:00:00	19.2
26. 6. 23. 오후 9:30:00	19.2
26. 6. 23. 오후 10:00:00	19.2
26. 6. 23. 오후 10:30:00	19.2
26. 6. 23. 오후 11:00:00	19.1
26. 6. 23. 오후 11:30:00	19.1
26. 6. 24. 오전 12:00:00	19.1
26. 6. 24. 오전 12:30:00	19.1
26. 6. 24. 오전 1:00:00	19.1
26. 6. 24. 오전 1:30:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 24. 오전 2:00:00	19.1
26. 6. 24. 오전 2:30:00	19.1
26. 6. 24. 오전 3:00:00	19.1
26. 6. 24. 오전 3:30:00	19.1
26. 6. 24. 오전 4:00:00	19.2
26. 6. 24. 오전 4:30:00	19.2
26. 6. 24. 오전 5:00:00	19.2
26. 6. 24. 오전 5:30:00	19.2
26. 6. 24. 오전 6:00:00	19.2
26. 6. 24. 오전 6:30:00	19.1
26. 6. 24. 오전 7:00:00	19.1
26. 6. 24. 오전 7:30:00	19.1
26. 6. 24. 오전 8:00:00	19.1
26. 6. 24. 오전 8:30:00	19.1
26. 6. 24. 오전 9:00:00	19.1
26. 6. 24. 오전 9:30:00	19.1
26. 6. 24. 오전 10:00:00	19.1
26. 6. 24. 오전 10:30:00	19.1
26. 6. 24. 오전 11:00:00	19.1
26. 6. 24. 오전 11:30:00	19.2
26. 6. 24. 오후 12:00:00	19.1
26. 6. 24. 오후 12:30:00	19.1
26. 6. 24. 오후 1:00:00	19.0
26. 6. 24. 오후 1:30:00	19.0
26. 6. 24. 오후 2:00:00	19.1
26. 6. 24. 오후 2:30:00	19.1
26. 6. 24. 오후 3:00:00	19.1
26. 6. 24. 오후 3:30:00	19.1
26. 6. 24. 오후 4:00:00	19.1
26. 6. 24. 오후 4:30:00	19.1
26. 6. 24. 오후 5:00:00	19.1
26. 6. 24. 오후 5:30:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 24. 오후 6:00:00	19.1
26. 6. 24. 오후 6:30:00	19.1
26. 6. 24. 오후 7:00:00	19.1
26. 6. 24. 오후 7:30:00	19.1
26. 6. 24. 오후 8:00:00	19.1
26. 6. 24. 오후 8:30:00	19.1
26. 6. 24. 오후 9:00:00	19.1
26. 6. 24. 오후 9:30:00	19.1
26. 6. 24. 오후 10:00:00	19.1
26. 6. 24. 오후 10:30:00	19.1
26. 6. 24. 오후 11:00:00	19.1
26. 6. 24. 오후 11:30:00	19.1
26. 6. 25. 오전 12:00:00	19.1
26. 6. 25. 오전 12:30:00	19.1
26. 6. 25. 오전 1:00:00	19.1
26. 6. 25. 오전 1:30:00	19.1
26. 6. 25. 오전 2:00:00	19.0
26. 6. 25. 오전 2:30:00	19.1
26. 6. 25. 오전 3:00:00	19.1
26. 6. 25. 오전 3:30:00	19.1
26. 6. 25. 오전 4:00:00	19.1
26. 6. 25. 오전 4:30:00	19.1
26. 6. 25. 오전 5:00:00	19.1
26. 6. 25. 오전 5:30:00	19.1
26. 6. 25. 오전 6:00:00	19.1
26. 6. 25. 오전 6:30:00	19.1
26. 6. 25. 오전 7:00:00	19.2
26. 6. 25. 오전 7:30:00	19.1
26. 6. 25. 오전 8:00:00	19.1
26. 6. 25. 오전 8:30:00	19.1
26. 6. 25. 오전 9:00:00	19.1
26. 6. 25. 오전 9:30:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 25. 오전 10:00:00	19.1
26. 6. 25. 오전 10:30:00	19.1
26. 6. 25. 오전 11:00:00	19.1
26. 6. 25. 오전 11:30:00	19.1
26. 6. 25. 오후 12:00:00	19.0
26. 6. 25. 오후 12:30:00	19.1
26. 6. 25. 오후 1:00:00	19.1
26. 6. 25. 오후 1:30:00	19.1
26. 6. 25. 오후 2:00:00	19.1
26. 6. 25. 오후 2:30:00	19.1
26. 6. 25. 오후 3:00:00	19.1
26. 6. 25. 오후 3:30:00	19.1
26. 6. 25. 오후 4:00:00	19.1
26. 6. 25. 오후 4:30:00	19.1
26. 6. 25. 오후 5:00:00	19.1
26. 6. 25. 오후 5:30:00	19.1
26. 6. 25. 오후 6:00:00	19.1
26. 6. 25. 오후 6:30:00	19.1
26. 6. 25. 오후 7:00:00	19.1
26. 6. 25. 오후 7:30:00	19.1
26. 6. 25. 오후 8:00:00	19.1
26. 6. 25. 오후 8:30:00	19.1
26. 6. 25. 오후 9:00:00	19.1
26. 6. 25. 오후 9:30:00	19.1
26. 6. 25. 오후 10:00:00	19.2
26. 6. 25. 오후 10:30:00	19.2
26. 6. 25. 오후 11:00:00	19.1
26. 6. 25. 오후 11:30:00	19.1
26. 6. 26. 오전 12:00:00	19.1
26. 6. 26. 오전 12:30:00	19.1
26. 6. 26. 오전 1:00:00	19.1
26. 6. 26. 오전 1:30:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 26. 오전 2:00:00	19.1
26. 6. 26. 오전 2:30:00	19.1
26. 6. 26. 오전 3:00:00	19.1
26. 6. 26. 오전 3:30:00	19.1
26. 6. 26. 오전 4:00:00	19.1
26. 6. 26. 오전 4:30:00	19.1
26. 6. 26. 오전 5:00:00	19.1
26. 6. 26. 오전 5:30:00	19.1
26. 6. 26. 오전 6:00:00	19.1
26. 6. 26. 오전 6:30:00	19.2
26. 6. 26. 오전 7:00:00	19.2
26. 6. 26. 오전 7:30:00	19.1
26. 6. 26. 오전 8:00:00	19.1
26. 6. 26. 오전 8:30:00	19.1
26. 6. 26. 오전 9:00:00	19.1
26. 6. 26. 오전 9:30:00	19.1
26. 6. 26. 오전 10:00:00	19.1
26. 6. 26. 오전 10:30:00	19.1
26. 6. 26. 오전 11:00:00	19.2
26. 6. 26. 오전 11:30:00	19.2
26. 6. 26. 오후 12:00:00	19.1
26. 6. 26. 오후 12:30:00	19.1
26. 6. 26. 오후 1:00:00	19.1
26. 6. 26. 오후 1:30:00	19.1
26. 6. 26. 오후 2:00:00	19.1
26. 6. 26. 오후 2:30:00	19.2
26. 6. 26. 오후 3:00:00	19.1
26. 6. 26. 오후 3:30:00	19.1
26. 6. 26. 오후 4:00:00	19.1
26. 6. 26. 오후 4:30:00	19.1
26. 6. 26. 오후 5:00:00	19.0
26. 6. 26. 오후 5:30:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25℃)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 26. 오후 6:00:00	19.1
26. 6. 26. 오후 6:30:00	19.1
26. 6. 26. 오후 7:00:00	19.1
26. 6. 26. 오후 7:30:00	19.1
26. 6. 26. 오후 8:00:00	19.1
26. 6. 26. 오후 8:30:00	19.1
26. 6. 26. 오후 9:00:00	19.1
26. 6. 26. 오후 9:30:00	19.1
26. 6. 26. 오후 10:00:00	19.1
26. 6. 26. 오후 10:30:00	19.1
26. 6. 26. 오후 11:00:00	19.1
26. 6. 26. 오후 11:30:00	19.1
26. 6. 27. 오전 12:00:00	19.1
26. 6. 27. 오전 12:30:00	19.1
26. 6. 27. 오전 1:00:00	19.1
26. 6. 27. 오전 1:30:00	19.1
26. 6. 27. 오전 2:00:00	19.1
26. 6. 27. 오전 2:30:00	19.1
26. 6. 27. 오전 3:00:00	19.1
26. 6. 27. 오전 3:30:00	19.2
26. 6. 27. 오전 4:00:00	19.2
26. 6. 27. 오전 4:30:00	19.2
26. 6. 27. 오전 5:00:00	19.1
26. 6. 27. 오전 5:30:00	19.1
26. 6. 27. 오전 6:00:00	19.1
26. 6. 27. 오전 6:30:00	19.1
26. 6. 27. 오전 7:00:00	19.2
26. 6. 27. 오전 7:30:00	19.2
26. 6. 27. 오전 8:00:00	19.1
26. 6. 27. 오전 8:30:00	19.1
26. 6. 27. 오전 9:00:00	19.1
26. 6. 27. 오전 9:30:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 27. 오전 10:00:00	19.1
26. 6. 27. 오전 10:30:00	19.1
26. 6. 27. 오전 11:00:00	19.1
26. 6. 27. 오전 11:30:00	19.1
26. 6. 27. 오후 12:00:00	19.2
26. 6. 27. 오후 12:30:00	19.1
26. 6. 27. 오후 1:00:00	19.1
26. 6. 27. 오후 1:30:00	19.1
26. 6. 27. 오후 2:00:00	19.1
26. 6. 27. 오후 2:30:00	19.1
26. 6. 27. 오후 3:00:00	19.1
26. 6. 27. 오후 3:30:00	19.2
26. 6. 27. 오후 4:00:00	19.2
26. 6. 27. 오후 4:30:00	19.1
26. 6. 27. 오후 5:00:00	19.1
26. 6. 27. 오후 5:30:00	19.1
26. 6. 27. 오후 6:00:00	19.1
26. 6. 27. 오후 6:30:00	19.1
26. 6. 27. 오후 7:00:00	19.1
26. 6. 27. 오후 7:30:00	19.2
26. 6. 27. 오후 8:00:00	19.1
26. 6. 27. 오후 8:30:00	19.1
26. 6. 27. 오후 9:00:00	19.1
26. 6. 27. 오후 9:30:00	19.1
26. 6. 27. 오후 10:00:00	19.1
26. 6. 27. 오후 10:30:00	19.1
26. 6. 27. 오후 11:00:00	19.2
26. 6. 27. 오후 11:30:00	19.1
26. 6. 28. 오전 12:00:00	19.1
26. 6. 28. 오전 12:30:00	19.1
26. 6. 28. 오전 1:00:00	19.1
26. 6. 28. 오전 1:30:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 28. 오전 2:00:00	19.0
26. 6. 28. 오전 2:30:00	19.1
26. 6. 28. 오전 3:00:00	19.1
26. 6. 28. 오전 3:30:00	19.1
26. 6. 28. 오전 4:00:00	19.1
26. 6. 28. 오전 4:30:00	19.1
26. 6. 28. 오전 5:00:00	19.1
26. 6. 28. 오전 5:30:00	19.1
26. 6. 28. 오전 6:00:00	19.1
26. 6. 28. 오전 6:30:00	19.2
26. 6. 28. 오전 7:00:00	19.2
26. 6. 28. 오전 7:30:00	19.2
26. 6. 28. 오전 8:00:00	19.1
26. 6. 28. 오전 8:30:00	19.1
26. 6. 28. 오전 9:00:00	19.1
26. 6. 28. 오전 9:30:00	19.2
26. 6. 28. 오전 10:00:00	19.1
26. 6. 28. 오전 10:30:00	19.1
26. 6. 28. 오전 11:00:00	19.1
26. 6. 28. 오전 11:30:00	19.1
26. 6. 28. 오후 12:00:00	19.1
26. 6. 28. 오후 12:30:00	19.1
26. 6. 28. 오후 1:00:00	19.1
26. 6. 28. 오후 1:30:00	19.1
26. 6. 28. 오후 2:00:00	19.2
26. 6. 28. 오후 2:30:00	19.2
26. 6. 28. 오후 3:00:00	19.1
26. 6. 28. 오후 3:30:00	19.1
26. 6. 28. 오후 4:00:00	19.1
26. 6. 28. 오후 4:30:00	19.1
26. 6. 28. 오후 5:00:00	19.1
26. 6. 28. 오후 5:30:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 28. 오후 6:00:00	19.2
26. 6. 28. 오후 6:30:00	19.2
26. 6. 28. 오후 7:00:00	19.2
26. 6. 28. 오후 7:30:00	19.1
26. 6. 28. 오후 8:00:00	19.1
26. 6. 28. 오후 8:30:00	19.1
26. 6. 28. 오후 9:00:00	19.1
26. 6. 28. 오후 9:30:00	19.1
26. 6. 28. 오후 10:00:00	19.2
26. 6. 28. 오후 10:30:00	19.2
26. 6. 28. 오후 11:00:00	19.1
26. 6. 28. 오후 11:30:00	19.1
26. 6. 29. 오전 12:00:00	19.1
26. 6. 29. 오전 12:30:00	19.0
26. 6. 29. 오전 1:00:00	19.1
26. 6. 29. 오전 1:30:00	19.1
26. 6. 29. 오전 2:00:00	19.1
26. 6. 29. 오전 2:30:00	19.2
26. 6. 29. 오전 3:00:00	19.1
26. 6. 29. 오전 3:30:00	19.1
26. 6. 29. 오전 4:00:00	19.1
26. 6. 29. 오전 4:30:00	19.1
26. 6. 29. 오전 5:00:00	19.1
26. 6. 29. 오전 5:30:00	19.2
26. 6. 29. 오전 6:00:00	19.1
26. 6. 29. 오전 6:30:00	19.1
26. 6. 29. 오전 7:00:00	19.1
26. 6. 29. 오전 7:30:00	19.1
26. 6. 29. 오전 8:00:00	19.2
26. 6. 29. 오전 8:30:00	19.2
26. 6. 29. 오전 9:00:00	19.2
26. 6. 29. 오전 9:30:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 29. 오전 10:00:00	19.1
26. 6. 29. 오전 10:30:00	19.1
26. 6. 29. 오전 11:00:00	19.0
26. 6. 29. 오전 11:30:00	19.1
26. 6. 29. 오후 12:00:00	19.1
26. 6. 29. 오후 12:30:00	19.1
26. 6. 29. 오후 1:00:00	19.1
26. 6. 29. 오후 1:30:00	19.1
26. 6. 29. 오후 2:00:00	19.1
26. 6. 29. 오후 2:30:00	19.1
26. 6. 29. 오후 3:00:00	19.1
26. 6. 29. 오후 3:30:00	19.1
26. 6. 29. 오후 4:00:00	19.1
26. 6. 29. 오후 4:30:00	19.1
26. 6. 29. 오후 5:00:00	19.1
26. 6. 29. 오후 5:30:00	19.1
26. 6. 29. 오후 6:00:00	19.1
26. 6. 29. 오후 6:30:00	19.1
26. 6. 29. 오후 7:00:00	19.2
26. 6. 29. 오후 7:30:00	19.1
26. 6. 29. 오후 8:00:00	19.1
26. 6. 29. 오후 8:30:00	19.1
26. 6. 29. 오후 9:00:00	19.1
26. 6. 29. 오후 9:30:00	19.1
26. 6. 29. 오후 10:00:00	19.2
26. 6. 29. 오후 10:30:00	19.1
26. 6. 29. 오후 11:00:00	19.1
26. 6. 29. 오후 11:30:00	19.1
26. 6. 30. 오전 12:00:00	19.1
26. 6. 30. 오전 12:30:00	19.1
26. 6. 30. 오전 1:00:00	19.1
26. 6. 30. 오전 1:30:00	19.2



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 30. 오전 2:00:00	19.1
26. 6. 30. 오전 2:30:00	19.1
26. 6. 30. 오전 3:00:00	19.1
26. 6. 30. 오전 3:30:00	19.1
26. 6. 30. 오전 4:00:00	19.1
26. 6. 30. 오전 4:30:00	19.2
26. 6. 30. 오전 5:00:00	19.1
26. 6. 30. 오전 5:30:00	19.1
26. 6. 30. 오전 6:00:00	19.1
26. 6. 30. 오전 6:30:00	19.1
26. 6. 30. 오전 7:00:00	19.2
26. 6. 30. 오전 7:30:00	19.2
26. 6. 30. 오전 8:00:00	19.2
26. 6. 30. 오전 8:30:00	19.1
26. 6. 30. 오전 9:00:00	19.0
26. 6. 30. 오전 9:30:00	19.1
26. 6. 30. 오전 10:00:00	19.1
26. 6. 30. 오전 10:30:00	19.1
26. 6. 30. 오전 11:00:00	19.1
26. 6. 30. 오전 11:30:00	19.1
26. 6. 30. 오후 12:00:00	19.1
26. 6. 30. 오후 12:30:00	19.1
26. 6. 30. 오후 1:00:00	19.1
26. 6. 30. 오후 1:30:00	19.2
26. 6. 30. 오후 2:00:00	19.2
26. 6. 30. 오후 2:30:00	19.2
26. 6. 30. 오후 3:00:00	19.2
26. 6. 30. 오후 3:30:00	19.2
26. 6. 30. 오후 4:00:00	19.2
26. 6. 30. 오후 4:30:00	19.1
26. 6. 30. 오후 5:00:00	19.1
26. 6. 30. 오후 5:30:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 6. 30. 오후 6:00:00	19.1
26. 6. 30. 오후 6:30:00	19.1
26. 6. 30. 오후 7:00:00	19.1
26. 6. 30. 오후 7:30:00	19.1
26. 6. 30. 오후 8:00:00	19.2
26. 6. 30. 오후 8:30:00	19.1
26. 6. 30. 오후 9:00:00	19.1
26. 6. 30. 오후 9:30:00	19.1
26. 6. 30. 오후 10:00:00	19.1
26. 6. 30. 오후 10:30:00	19.1
26. 6. 30. 오후 11:00:00	19.1
26. 6. 30. 오후 11:30:00	19.1



IP Storage Temperature Log_Monthly report

2026년 6월 1일 오전 12시 0분 0초 KST - 2026년 7월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 7월 1일 오전 6시 2분 53초 KST

알람 테이블

사용할 수 있는 값이 없습니다.