

Period of Record	Year. Month	2026.8
Storage	Location	동국대학교일산병원 임상시험센터 약국
	Description	전기냉장냉동기
	Manufacturer	SANYO
	Model Name	MPR-720R
	Serial Number	07120004
	Temperature Range	15~25°C
	Setting Point	19.0°C
Thermometer	Description	Digital Thermometer
	Manufacturer	SATO
	Model Name	PC-3500
	Serial Number	825289
	Min/Max Record	Reset after recording
Thermometer Calibration Certification	Date of Calibration	2026.1.27
	Date of Calibration Due	2027.1.27
	Certificate No.	SC2601-02474-4
Temperature deviation alarm system	Manufacturer	Daeil Tech
	Model Name	ACMS-DT01
	Serial Number	DT15-ACMS-000348-K
	Alarm threshold	Lower : 15°C, Upper : 25°C, alarm in 5 minutes
	Alarm Method	Call to Pharmacists
	Alarm Clear Method	000

Date	Time	Temperature(°C)			Verified by	Comment
		Current	Min	Max		
2026.08.01						
2026.08.02						
2026.08.03						
2026.08.04						
2026.08.05						
2026.08.06						
2026.08.07						
2026.08.08						
2026.08.09						
2026.08.10						
2026.08.11		Min	Current	Max		
2026.08.12	p:20	18	18	18	aml	원도메 115도
2026.08.13	p:20	18	20	19	aml	
2026.08.14	p:20	18	20	19	aml	
2026.08.15						광복절
2026.08.16						
2026.08.17						대체공휴일
2026.08.18	p:20	17	20	20	aml	
2026.08.19	p:20	18	20	20	aml	
2026.08.20	p:20	18	20	20	aml	
2026.08.21	p:20	18	18	20	aml	
2026.08.22						
2026.08.23						
2026.08.24	p:20	18	20	20	aml	
2026.08.25	p:20	18	18	20	aml	
2026.08.26	p:10	18	18	20	aml	
2026.08.27	p:20	18	18	20	aml	
2026.08.28	p:20	18	18	20	aml	
2026.08.29						
2026.08.30						
2026.08.31	p:20	18	20	20	aml	

IP STORAGE TEMPERATURE LOG

Period of Record	Year. Month	2026.8	
Storage	Location	동국대학교일산병원 임상시험센터 약국	9134Room
	Description	전기냉장냉동기	항온고 1
	Manufacturer	SANYO	
	Model Name	MPR-720R	
	Serial Number	07120004	
	Temperature Range	15~25℃	
	Setting Point	19.0℃	
Thermometer (including alarm system)	Description	Wi-Fi data logger	
	Manufacturer	Testo Korea	
	Model Name	Saveris2-T3	
	Serial Number	0054704396	
	Sensor Channel	Channel 1	
	Min/Max Record	Monthly report	This month Min (19.0)℃, Max(19.4)℃
Thermometer Calibration Certification	Date of Calibration	2025.11.25	
	Date of Calibration Due	2026.11.25	
	Certificate No.	A25-2356-001	
Alarm System	Alarm threshold	Lower : 15.5℃, Upper : 24.5℃, immediate alarm without delay	
	Alarm Method	Mail, SMS to Pharmacists	

Date	Time	Temperature(℃)	Deviation	Verified by	Comment
2026.08.01			☐ Y / ☐ N		
2026.08.02			☐ Y / ☐ N		
2026.08.03	P:20	19.2	☐ Y / ☒ N	ml	
2026.08.04	P:20	19.1	☐ Y / ☒ N	ml	
2026.08.05	P:05	19.1	☐ Y / ☒ N	ml	
2026.08.06	P:20	19.3	☐ Y / ☒ N	ml	
2026.08.07	P:30	19.3	☐ Y / ☒ N	정종민	
2026.08.08			☐ Y / ☐ N		
2026.08.09			☐ Y / ☐ N		
2026.08.10	P:20	19.1	☐ Y / ☒ N	ml	
2026.08.11	P:20	19.1	☐ Y / ☒ N	ml	
2026.08.12	P:20	19.1	☐ Y / ☒ N	ml	양조제 교체
2026.08.13			☐ Y / ☐ N		
2026.08.14			☐ Y / ☐ N		
2026.08.15			☐ Y / ☐ N		광복절
2026.08.16			☐ Y / ☐ N		
2026.08.17			☐ Y / ☐ N		대체공휴일
2026.08.18			☐ Y / ☐ N		
2026.08.19			☐ Y / ☐ N		
2026.08.20			☐ Y / ☐ N		
2026.08.21			☐ Y / ☐ N		
2026.08.22			☐ Y / ☐ N		
2026.08.23			☐ Y / ☐ N		
2026.08.24			☐ Y / ☐ N		
2026.08.25			☐ Y / ☐ N		
2026.08.26			☐ Y / ☐ N		
2026.08.27			☐ Y / ☐ N		
2026.08.28			☐ Y / ☐ N		
2026.08.29			☐ Y / ☐ N		
2026.08.30			☐ Y / ☐ N		
2026.08.31			☐ Y / ☐ N		

date of copying


sign



내용

동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1

임상시험용 의약품(항온고 15~25°C)_CH1





테스토 사베리스 보고서

kyon410@daum.net
생성됨: 2026년 8월 12일 오전 10시 16분 58초 KST

2026년 8월 1일 오전 12시 0분 0초 KST - 2026년 8월 12일 오후 11시 59분 59초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C) CH1
타입 스탬프 [Asia/Seoul]	온도 [°C]
Min	19.0
Max	19.4
Average	19.2
26. 8. 1. 오전 12:00:00	19.2
26. 8. 1. 오전 12:30:00	19.1
26. 8. 1. 오전 1:00:00	19.1
26. 8. 1. 오전 1:30:00	19.2
26. 8. 1. 오전 2:00:00	19.2
26. 8. 1. 오전 2:30:00	19.2
26. 8. 1. 오전 3:00:00	19.3
26. 8. 1. 오전 3:30:00	19.3
26. 8. 1. 오전 4:00:00	19.2
26. 8. 1. 오전 4:30:00	19.2
26. 8. 1. 오전 5:00:00	19.1
26. 8. 1. 오전 5:30:00	19.2
26. 8. 1. 오전 6:00:00	19.2
26. 8. 1. 오전 6:30:00	19.3
26. 8. 1. 오전 7:00:00	19.3
26. 8. 1. 오전 7:30:00	19.2
26. 8. 1. 오전 8:00:00	19.1
26. 8. 1. 오전 8:30:00	19.1
26. 8. 1. 오전 9:00:00	19.2
26. 8. 1. 오전 9:30:00	19.2
26. 8. 1. 오전 10:00:00	19.3
26. 8. 1. 오전 10:30:00	19.3
26. 8. 1. 오전 11:00:00	19.2
26. 8. 1. 오전 11:30:00	19.1
26. 8. 1. 오후 12:00:00	19.1
26. 8. 1. 오후 12:30:00	19.2
26. 8. 1. 오후 1:00:00	19.2
26. 8. 1. 오후 1:30:00	19.3
26. 8. 1. 오후 2:00:00	19.2



테스토 사베리스 보고서

kyon410@daum.net
생성됨: 2026년 8월 12일 오전 10시 16분 58초 KST

2026년 8월 1일 오전 12시 0분 0초 KST - 2026년 8월 12일 오후 11시 59분 59초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C) CH1
타입 스탬프 [Asia/Seoul]	온도 [°C]
26. 8. 1. 오후 2:30:00	19.1
26. 8. 1. 오후 3:00:00	19.2
26. 8. 1. 오후 3:30:00	19.2
26. 8. 1. 오후 4:00:00	19.3
26. 8. 1. 오후 4:30:00	19.3
26. 8. 1. 오후 5:00:00	19.1
26. 8. 1. 오후 5:30:00	19.2
26. 8. 1. 오후 6:00:00	19.2
26. 8. 1. 오후 6:30:00	19.3
26. 8. 1. 오후 7:00:00	19.3
26. 8. 1. 오후 7:30:00	19.3
26. 8. 1. 오후 8:00:00	19.1
26. 8. 1. 오후 8:30:00	19.2
26. 8. 1. 오후 9:00:00	19.2
26. 8. 1. 오후 9:30:00	19.3
26. 8. 1. 오후 10:00:00	19.3
26. 8. 1. 오후 10:30:00	19.2
26. 8. 1. 오후 11:00:00	19.1
26. 8. 1. 오후 11:30:00	19.3
26. 8. 2. 오전 12:00:00	19.3
26. 8. 2. 오전 12:30:00	19.3
26. 8. 2. 오전 1:00:00	19.2
26. 8. 2. 오전 1:30:00	19.1
26. 8. 2. 오전 2:00:00	19.1
26. 8. 2. 오전 2:30:00	19.2
26. 8. 2. 오전 3:00:00	19.3
26. 8. 2. 오전 3:30:00	19.3
26. 8. 2. 오전 4:00:00	19.2
26. 8. 2. 오전 4:30:00	19.1
26. 8. 2. 오전 5:00:00	19.2
26. 8. 2. 오전 5:30:00	19.2
26. 8. 2. 오전 6:00:00	19.3



테스토 사베리스 보고서

kyon410@daum.net
생성됨: 2026년 8월 12일 오전 10시 16분 58초 KST

2026년 8월 1일 오전 12시 0분 0초 KST - 2026년 8월 12일 오후 11시 59분 59초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C) CH1
타입 스탬프 [Asia/Seoul]	온도 [°C]
26. 8. 2. 오전 6:30:00	19.3
26. 8. 2. 오전 7:00:00	19.3
26. 8. 2. 오전 7:30:00	19.1
26. 8. 2. 오전 8:00:00	19.1
26. 8. 2. 오전 8:30:00	19.2
26. 8. 2. 오전 9:00:00	19.2
26. 8. 2. 오전 9:30:00	19.3
26. 8. 2. 오전 10:00:00	19.2
26. 8. 2. 오전 10:30:00	19.1
26. 8. 2. 오전 11:00:00	19.2
26. 8. 2. 오전 11:30:00	19.2
26. 8. 2. 오후 12:00:00	19.3
26. 8. 2. 오후 12:30:00	19.3
26. 8. 2. 오후 1:00:00	19.2
26. 8. 2. 오후 1:30:00	19.1
26. 8. 2. 오후 2:00:00	19.2
26. 8. 2. 오후 2:30:00	19.3
26. 8. 2. 오후 3:00:00	19.3
26. 8. 2. 오후 3:30:00	19.2
26. 8. 2. 오후 4:00:00	19.1
26. 8. 2. 오후 4:30:00	19.2
26. 8. 2. 오후 5:00:00	19.2
26. 8. 2. 오후 5:30:00	19.3
26. 8. 2. 오후 6:00:00	19.4
26. 8. 2. 오후 6:30:00	19.3
26. 8. 2. 오후 7:00:00	19.1
26. 8. 2. 오후 7:30:00	19.1
26. 8. 2. 오후 8:00:00	19.2
26. 8. 2. 오후 8:30:00	19.3
26. 8. 2. 오후 9:00:00	19.4
26. 8. 2. 오후 9:30:00	19.2
26. 8. 2. 오후 10:00:00	19.1



테스토 사베리스 보고서

kyon410@daum.net
생성됨: 2026년 8월 12일 오전 10시 16분 58초 KST

2026년 8월 1일 오전 12시 0분 0초 KST - 2026년 8월 12일 오후 11시 59분 59초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C) CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 8. 2. 오후 10:30:00	19.2
26. 8. 2. 오후 11:00:00	19.2
26. 8. 2. 오후 11:30:00	19.3
26. 8. 3. 오전 12:00:00	19.3
26. 8. 3. 오전 12:30:00	19.3
26. 8. 3. 오전 1:00:00	19.3
26. 8. 3. 오전 1:30:00	19.1
26. 8. 3. 오전 2:00:00	19.2
26. 8. 3. 오전 2:30:00	19.2
26. 8. 3. 오전 3:00:00	19.2
26. 8. 3. 오전 3:30:00	19.3
26. 8. 3. 오전 4:00:00	19.3
26. 8. 3. 오전 4:30:00	19.3
26. 8. 3. 오전 5:00:00	19.3
26. 8. 3. 오전 5:30:00	19.2
26. 8. 3. 오전 6:00:00	19.1
26. 8. 3. 오전 6:30:00	19.1
26. 8. 3. 오전 7:00:00	19.2
26. 8. 3. 오전 7:30:00	19.2
26. 8. 3. 오전 8:00:00	19.3
26. 8. 3. 오전 8:30:00	19.3
26. 8. 3. 오전 9:00:00	19.2
26. 8. 3. 오전 9:30:00	19.1
26. 8. 3. 오전 10:00:00	19.1
26. 8. 3. 오전 10:30:00	19.2
26. 8. 3. 오전 11:00:00	19.2
26. 8. 3. 오전 11:30:00	19.3
26. 8. 3. 오후 12:00:00	19.3
26. 8. 3. 오후 12:30:00	19.1
26. 8. 3. 오후 1:00:00	19.1
26. 8. 3. 오후 1:30:00	19.2
26. 8. 3. 오후 2:00:00	19.3



테스토 사베리스 보고서

kyon410@daum.net
생성됨: 2026년 8월 12일 오전 10시 16분 58초 KST

2026년 8월 1일 오전 12시 0분 0초 KST - 2026년 8월 12일 오후 11시 59분 59초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C) CH1
타입 스탬프 [Asia/Seoul]	온도 [°C]
26. 8. 3. 오후 2:30:00	19.3
26. 8. 3. 오후 3:00:00	19.4
26. 8. 3. 오후 3:30:00	19.2
26. 8. 3. 오후 4:00:00	19.2
26. 8. 3. 오후 4:30:00	19.2
26. 8. 3. 오후 5:00:00	19.3
26. 8. 3. 오후 5:30:00	19.3
26. 8. 3. 오후 6:00:00	19.1
26. 8. 3. 오후 6:30:00	19.2
26. 8. 3. 오후 7:00:00	19.3
26. 8. 3. 오후 7:30:00	19.3
26. 8. 3. 오후 8:00:00	19.3
26. 8. 3. 오후 8:30:00	19.2
26. 8. 3. 오후 9:00:00	19.2
26. 8. 3. 오후 9:30:00	19.2
26. 8. 3. 오후 10:00:00	19.3
26. 8. 3. 오후 10:30:00	19.3
26. 8. 3. 오후 11:00:00	19.1
26. 8. 3. 오후 11:30:00	19.1
26. 8. 4. 오전 12:00:00	19.2
26. 8. 4. 오전 12:30:00	19.3
26. 8. 4. 오전 1:00:00	19.3
26. 8. 4. 오전 1:30:00	19.1
26. 8. 4. 오전 2:00:00	19.1
26. 8. 4. 오전 2:30:00	19.2
26. 8. 4. 오전 3:00:00	19.3
26. 8. 4. 오전 3:30:00	19.3
26. 8. 4. 오전 4:00:00	19.2
26. 8. 4. 오전 4:30:00	19.2
26. 8. 4. 오전 5:00:00	19.2
26. 8. 4. 오전 5:30:00	19.4
26. 8. 4. 오전 6:00:00	19.3



테스토 사베리스 보고서

kyon410@daum.net
생성됨: 2026년 8월 12일 오전 10시 16분 58초 KST

2026년 8월 1일 오전 12시 0분 0초 KST - 2026년 8월 12일 오후 11시 59분 59초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C) CH1
타입 스탬프 [Asia/Seoul]	온도 [°C]
26. 8. 4. 오전 6:30:00	19.2
26. 8. 4. 오전 7:00:00	19.2
26. 8. 4. 오전 7:30:00	19.2
26. 8. 4. 오전 8:00:00	19.3
26. 8. 4. 오전 8:30:00	19.3
26. 8. 4. 오전 9:00:00	19.1
26. 8. 4. 오전 9:30:00	19.1
26. 8. 4. 오전 10:00:00	19.2
26. 8. 4. 오전 10:30:00	19.3
26. 8. 4. 오전 11:00:00	19.3
26. 8. 4. 오전 11:30:00	19.2
26. 8. 4. 오후 12:00:00	19.1
26. 8. 4. 오후 12:30:00	19.2
26. 8. 4. 오후 1:00:00	19.3
26. 8. 4. 오후 1:30:00	19.3
26. 8. 4. 오후 2:00:00	19.2
26. 8. 4. 오후 2:30:00	19.2
26. 8. 4. 오후 3:00:00	19.3
26. 8. 4. 오후 3:30:00	19.4
26. 8. 4. 오후 4:00:00	19.2
26. 8. 4. 오후 4:30:00	19.2
26. 8. 4. 오후 5:00:00	19.3
26. 8. 4. 오후 5:30:00	19.4
26. 8. 4. 오후 6:00:00	19.3
26. 8. 4. 오후 6:30:00	19.1
26. 8. 4. 오후 7:00:00	19.2
26. 8. 4. 오후 7:30:00	19.3
26. 8. 4. 오후 8:00:00	19.3
26. 8. 4. 오후 8:30:00	19.2
26. 8. 4. 오후 9:00:00	19.2
26. 8. 4. 오후 9:30:00	19.3
26. 8. 4. 오후 10:00:00	19.4



테스토 사베리스 보고서

kyon410@daum.net
생성됨: 2026년 8월 12일 오전 10시 16분 58초 KST

2026년 8월 1일 오전 12시 0분 0초 KST - 2026년 8월 12일 오후 11시 59분 59초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C) CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 8. 4. 오후 10:30:00	19.2
26. 8. 4. 오후 11:00:00	19.2
26. 8. 4. 오후 11:30:00	19.2
26. 8. 5. 오전 12:00:00	19.3
26. 8. 5. 오전 12:30:00	19.3
26. 8. 5. 오전 1:00:00	19.2
26. 8. 5. 오전 1:30:00	19.1
26. 8. 5. 오전 2:00:00	19.2
26. 8. 5. 오전 2:30:00	19.3
26. 8. 5. 오전 3:00:00	19.3
26. 8. 5. 오전 3:30:00	19.1
26. 8. 5. 오전 4:00:00	19.2
26. 8. 5. 오전 4:30:00	19.3
26. 8. 5. 오전 5:00:00	19.3
26. 8. 5. 오전 5:30:00	19.3
26. 8. 5. 오전 6:00:00	19.2
26. 8. 5. 오전 6:30:00	19.2
26. 8. 5. 오전 7:00:00	19.2
26. 8. 5. 오전 7:30:00	19.3
26. 8. 5. 오전 8:00:00	19.3
26. 8. 5. 오전 8:30:00	19.1
26. 8. 5. 오전 9:00:00	19.1
26. 8. 5. 오전 9:30:00	19.2
26. 8. 5. 오전 10:00:00	19.2
26. 8. 5. 오전 10:30:00	19.3
26. 8. 5. 오전 11:00:00	19.3
26. 8. 5. 오전 11:30:00	19.1
26. 8. 5. 오후 12:00:00	19.2
26. 8. 5. 오후 12:30:00	19.3
26. 8. 5. 오후 1:00:00	19.3
26. 8. 5. 오후 1:30:00	19.1
26. 8. 5. 오후 2:00:00	19.2



테스토 사베리스 보고서

kyon410@daum.net
생성됨: 2026년 8월 12일 오전 10시 16분 58초 KST

2026년 8월 1일 오전 12시 0분 0초 KST - 2026년 8월 12일 오후 11시 59분 59초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C) CH1
타입 스탬프 [Asia/Seoul]	온도 [°C]
26. 8. 5. 오후 2:30:00	19.2
26. 8. 5. 오후 3:00:00	19.3
26. 8. 5. 오후 3:30:00	19.2
26. 8. 5. 오후 4:00:00	19.2
26. 8. 5. 오후 4:30:00	19.3
26. 8. 5. 오후 5:00:00	19.4
26. 8. 5. 오후 5:30:00	19.2
26. 8. 5. 오후 6:00:00	19.2
26. 8. 5. 오후 6:30:00	19.3
26. 8. 5. 오후 7:00:00	19.4
26. 8. 5. 오후 7:30:00	19.2
26. 8. 5. 오후 8:00:00	19.2
26. 8. 5. 오후 8:30:00	19.3
26. 8. 5. 오후 9:00:00	19.3
26. 8. 5. 오후 9:30:00	19.4
26. 8. 5. 오후 10:00:00	19.2
26. 8. 5. 오후 10:30:00	19.1
26. 8. 5. 오후 11:00:00	19.2
26. 8. 5. 오후 11:30:00	19.3
26. 8. 6. 오전 12:00:00	19.3
26. 8. 6. 오전 12:30:00	19.3
26. 8. 6. 오전 1:00:00	19.2
26. 8. 6. 오전 1:30:00	19.2
26. 8. 6. 오전 2:00:00	19.2
26. 8. 6. 오전 2:30:00	19.4
26. 8. 6. 오전 3:00:00	19.2
26. 8. 6. 오전 3:30:00	19.1
26. 8. 6. 오전 4:00:00	19.1
26. 8. 6. 오전 4:30:00	19.2
26. 8. 6. 오전 5:00:00	19.3
26. 8. 6. 오전 5:30:00	19.3
26. 8. 6. 오전 6:00:00	19.3



테스토 사베리스 보고서

kyon410@daum.net
생성됨: 2026년 8월 12일 오전 10시 16분 58초 KST

2026년 8월 1일 오전 12시 0분 0초 KST - 2026년 8월 12일 오후 11시 59분 59초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C) CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 8. 6. 오전 6:30:00	19.1
26. 8. 6. 오전 7:00:00	19.1
26. 8. 6. 오전 7:30:00	19.2
26. 8. 6. 오전 8:00:00	19.3
26. 8. 6. 오전 8:30:00	19.3
26. 8. 6. 오전 9:00:00	19.3
26. 8. 6. 오전 9:30:00	19.1
26. 8. 6. 오전 10:00:00	19.1
26. 8. 6. 오전 10:30:00	19.2
26. 8. 6. 오전 11:00:00	19.3
26. 8. 6. 오전 11:30:00	19.4
26. 8. 6. 오후 12:00:00	19.2
26. 8. 6. 오후 12:30:00	19.1
26. 8. 6. 오후 1:00:00	19.2
26. 8. 6. 오후 1:30:00	19.2
26. 8. 6. 오후 2:00:00	19.3
26. 8. 6. 오후 2:30:00	19.3
26. 8. 6. 오후 3:00:00	19.2
26. 8. 6. 오후 3:30:00	19.2
26. 8. 6. 오후 4:00:00	19.3
26. 8. 6. 오후 4:30:00	19.3
26. 8. 6. 오후 5:00:00	19.3
26. 8. 6. 오후 5:30:00	19.2
26. 8. 6. 오후 6:00:00	19.1
26. 8. 6. 오후 6:30:00	19.2
26. 8. 6. 오후 7:00:00	19.3
26. 8. 6. 오후 7:30:00	19.3
26. 8. 6. 오후 8:00:00	19.3
26. 8. 6. 오후 8:30:00	19.2
26. 8. 6. 오후 9:00:00	19.1
26. 8. 6. 오후 9:30:00	19.2
26. 8. 6. 오후 10:00:00	19.2



테스토 사베리스 보고서

kyon410@daum.net
생성됨: 2026년 8월 12일 오전 10시 16분 58초 KST

2026년 8월 1일 오전 12시 0분 0초 KST - 2026년 8월 12일 오후 11시 59분 59초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C) CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 8. 6. 오후 10:30:00	19.3
26. 8. 6. 오후 11:00:00	19.3
26. 8. 6. 오후 11:30:00	19.3
26. 8. 7. 오전 12:00:00	19.2
26. 8. 7. 오전 12:30:00	19.1
26. 8. 7. 오전 1:00:00	19.1
26. 8. 7. 오전 1:30:00	19.2
26. 8. 7. 오전 2:00:00	19.3
26. 8. 7. 오전 2:30:00	19.4
26. 8. 7. 오전 3:00:00	19.3
26. 8. 7. 오전 3:30:00	19.1
26. 8. 7. 오전 4:00:00	19.1
26. 8. 7. 오전 4:30:00	19.2
26. 8. 7. 오전 5:00:00	19.2
26. 8. 7. 오전 5:30:00	19.3
26. 8. 7. 오전 6:00:00	19.3
26. 8. 7. 오전 6:30:00	19.3
26. 8. 7. 오전 7:00:00	19.2
26. 8. 7. 오전 7:30:00	19.1
26. 8. 7. 오전 8:00:00	19.1
26. 8. 7. 오전 8:30:00	19.2
26. 8. 7. 오전 9:00:00	19.2
26. 8. 7. 오전 9:30:00	19.3
26. 8. 7. 오전 10:00:00	19.3
26. 8. 7. 오전 10:30:00	19.3
26. 8. 7. 오전 11:00:00	19.1
26. 8. 7. 오전 11:30:00	19.1
26. 8. 7. 오후 12:00:00	19.1
26. 8. 7. 오후 12:30:00	19.2
26. 8. 7. 오후 1:00:00	19.2
26. 8. 7. 오후 1:30:00	19.3
26. 8. 7. 오후 2:00:00	19.4



테스토 사베리스 보고서

kyon410@daum.net
생성됨: 2026년 8월 12일 오전 10시 16분 58초 KST

2026년 8월 1일 오전 12시 0분 0초 KST - 2026년 8월 12일 오후 11시 59분 59초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C) CH1
타입 스탬프 [Asia/Seoul]	온도 [°C]
26. 8. 7. 오후 2:30:00	19.2
26. 8. 7. 오후 3:00:00	19.1
26. 8. 7. 오후 3:30:00	19.2
26. 8. 7. 오후 4:00:00	19.2
26. 8. 7. 오후 4:30:00	19.2
26. 8. 7. 오후 5:00:00	19.3
26. 8. 7. 오후 5:30:00	19.3
26. 8. 7. 오후 6:00:00	19.2
26. 8. 7. 오후 6:30:00	19.1
26. 8. 7. 오후 7:00:00	19.1
26. 8. 7. 오후 7:30:00	19.2
26. 8. 7. 오후 8:00:00	19.3
26. 8. 7. 오후 8:30:00	19.3
26. 8. 7. 오후 9:00:00	19.3
26. 8. 7. 오후 9:30:00	19.2
26. 8. 7. 오후 10:00:00	19.1
26. 8. 7. 오후 10:30:00	19.1
26. 8. 7. 오후 11:00:00	19.2
26. 8. 7. 오후 11:30:00	19.2
26. 8. 8. 오전 12:00:00	19.3
26. 8. 8. 오전 12:30:00	19.3
26. 8. 8. 오전 1:00:00	19.3
26. 8. 8. 오전 1:30:00	19.2
26. 8. 8. 오전 2:00:00	19.1
26. 8. 8. 오전 2:30:00	19.1
26. 8. 8. 오전 3:00:00	19.1
26. 8. 8. 오전 3:30:00	19.2
26. 8. 8. 오전 4:00:00	19.2
26. 8. 8. 오전 4:30:00	19.3
26. 8. 8. 오전 5:00:00	19.3
26. 8. 8. 오전 5:30:00	19.3
26. 8. 8. 오전 6:00:00	19.3



테스토 사베리스 보고서

kyon410@daum.net
생성됨: 2026년 8월 12일 오전 10시 16분 58초 KST

2026년 8월 1일 오전 12시 0분 0초 KST - 2026년 8월 12일 오후 11시 59분 59초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C) CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 8. 8. 오전 6:30:00	19.3
26. 8. 8. 오전 7:00:00	19.2
26. 8. 8. 오전 7:30:00	19.1
26. 8. 8. 오전 8:00:00	19.1
26. 8. 8. 오전 8:30:00	19.1
26. 8. 8. 오전 9:00:00	19.1
26. 8. 8. 오전 9:30:00	19.1
26. 8. 8. 오전 10:00:00	19.2
26. 8. 8. 오전 10:30:00	19.2
26. 8. 8. 오전 11:00:00	19.2
26. 8. 8. 오전 11:30:00	19.3
26. 8. 8. 오후 12:00:00	19.2
26. 8. 8. 오후 12:30:00	19.2
26. 8. 8. 오후 1:00:00	19.1
26. 8. 8. 오후 1:30:00	19.1
26. 8. 8. 오후 2:00:00	19.1
26. 8. 8. 오후 2:30:00	19.2
26. 8. 8. 오후 3:00:00	19.2
26. 8. 8. 오후 3:30:00	19.4
26. 8. 8. 오후 4:00:00	19.3
26. 8. 8. 오후 4:30:00	19.3
26. 8. 8. 오후 5:00:00	19.1
26. 8. 8. 오후 5:30:00	19.1
26. 8. 8. 오후 6:00:00	19.1
26. 8. 8. 오후 6:30:00	19.2
26. 8. 8. 오후 7:00:00	19.2
26. 8. 8. 오후 7:30:00	19.3
26. 8. 8. 오후 8:00:00	19.3
26. 8. 8. 오후 8:30:00	19.3
26. 8. 8. 오후 9:00:00	19.3
26. 8. 8. 오후 9:30:00	19.2
26. 8. 8. 오후 10:00:00	19.1



테스토 사베리스 보고서

kyon410@daum.net
생성됨: 2026년 8월 12일 오전 10시 16분 58초 KST

2026년 8월 1일 오전 12시 0분 0초 KST - 2026년 8월 12일 오후 11시 59분 59초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C) CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 8. 8. 오후 10:30:00	19.1
26. 8. 8. 오후 11:00:00	19.1
26. 8. 8. 오후 11:30:00	19.1
26. 8. 9. 오전 12:00:00	19.1
26. 8. 9. 오전 12:30:00	19.1
26. 8. 9. 오전 1:00:00	19.1
26. 8. 9. 오전 1:30:00	19.1
26. 8. 9. 오전 2:00:00	19.1
26. 8. 9. 오전 2:30:00	19.1
26. 8. 9. 오전 3:00:00	19.1
26. 8. 9. 오전 3:30:00	19.2
26. 8. 9. 오전 4:00:00	19.2
26. 8. 9. 오전 4:30:00	19.2
26. 8. 9. 오전 5:00:00	19.2
26. 8. 9. 오전 5:30:00	19.2
26. 8. 9. 오전 6:00:00	19.2
26. 8. 9. 오전 6:30:00	19.2
26. 8. 9. 오전 7:00:00	19.2
26. 8. 9. 오전 7:30:00	19.2
26. 8. 9. 오전 8:00:00	19.2
26. 8. 9. 오전 8:30:00	19.2
26. 8. 9. 오전 9:00:00	19.2
26. 8. 9. 오전 9:30:00	19.3
26. 8. 9. 오전 10:00:00	19.2
26. 8. 9. 오전 10:30:00	19.2
26. 8. 9. 오전 11:00:00	19.1
26. 8. 9. 오전 11:30:00	19.1
26. 8. 9. 오후 12:00:00	19.1
26. 8. 9. 오후 12:30:00	19.0
26. 8. 9. 오후 1:00:00	19.0
26. 8. 9. 오후 1:30:00	19.0
26. 8. 9. 오후 2:00:00	19.1



테스토 사베리스 보고서

kyon410@daum.net
생성됨: 2026년 8월 12일 오전 10시 16분 58초 KST

2026년 8월 1일 오전 12시 0분 0초 KST - 2026년 8월 12일 오후 11시 59분 59초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C) CH1
타입 스탬프 [Asia/Seoul]	온도 [°C]
26. 8. 9. 오후 2:30:00	19.1
26. 8. 9. 오후 3:00:00	19.1
26. 8. 9. 오후 3:30:00	19.1
26. 8. 9. 오후 4:00:00	19.1
26. 8. 9. 오후 4:30:00	19.2
26. 8. 9. 오후 5:00:00	19.2
26. 8. 9. 오후 5:30:00	19.2
26. 8. 9. 오후 6:00:00	19.3
26. 8. 9. 오후 6:30:00	19.3
26. 8. 9. 오후 7:00:00	19.3
26. 8. 9. 오후 7:30:00	19.3
26. 8. 9. 오후 8:00:00	19.3
26. 8. 9. 오후 8:30:00	19.3
26. 8. 9. 오후 9:00:00	19.2
26. 8. 9. 오후 9:30:00	19.2
26. 8. 9. 오후 10:00:00	19.2
26. 8. 9. 오후 10:30:00	19.2
26. 8. 9. 오후 11:00:00	19.1
26. 8. 9. 오후 11:30:00	19.1
26. 8. 10. 오전 12:00:00	19.0
26. 8. 10. 오전 12:30:00	19.0
26. 8. 10. 오전 1:00:00	19.0
26. 8. 10. 오전 1:30:00	19.1
26. 8. 10. 오전 2:00:00	19.0
26. 8. 10. 오전 2:30:00	19.1
26. 8. 10. 오전 3:00:00	19.0
26. 8. 10. 오전 3:30:00	19.1
26. 8. 10. 오전 4:00:00	19.1
26. 8. 10. 오전 4:30:00	19.1
26. 8. 10. 오전 5:00:00	19.1
26. 8. 10. 오전 5:30:00	19.1
26. 8. 10. 오전 6:00:00	19.2



테스토 사베리스 보고서

kyon410@daum.net
생성됨: 2026년 8월 12일 오전 10시 16분 58초 KST

2026년 8월 1일 오전 12시 0분 0초 KST - 2026년 8월 12일 오후 11시 59분 59초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C) CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 8. 10. 오전 6:30:00	19.2
26. 8. 10. 오전 7:00:00	19.2
26. 8. 10. 오전 7:30:00	19.2
26. 8. 10. 오전 8:00:00	19.2
26. 8. 10. 오전 8:30:00	19.2
26. 8. 10. 오전 9:00:00	19.1
26. 8. 10. 오전 9:30:00	19.1
26. 8. 10. 오전 10:00:00	19.1
26. 8. 10. 오전 10:30:00	19.1
26. 8. 10. 오전 11:00:00	19.1
26. 8. 10. 오전 11:30:00	19.1
26. 8. 10. 오후 12:00:00	19.1
26. 8. 10. 오후 12:30:00	19.1
26. 8. 10. 오후 1:00:00	19.1
26. 8. 10. 오후 1:30:00	19.2
26. 8. 10. 오후 2:00:00	19.2
26. 8. 10. 오후 2:30:00	19.2
26. 8. 10. 오후 3:00:00	19.3
26. 8. 10. 오후 3:30:00	19.3
26. 8. 10. 오후 4:00:00	19.3
26. 8. 10. 오후 4:30:00	19.3
26. 8. 10. 오후 5:00:00	19.3
26. 8. 10. 오후 5:30:00	19.3
26. 8. 10. 오후 6:00:00	19.2
26. 8. 10. 오후 6:30:00	19.1
26. 8. 10. 오후 7:00:00	19.1
26. 8. 10. 오후 7:30:00	19.1
26. 8. 10. 오후 8:00:00	19.1
26. 8. 10. 오후 8:30:00	19.1
26. 8. 10. 오후 9:00:00	19.1
26. 8. 10. 오후 9:30:00	19.1
26. 8. 10. 오후 10:00:00	19.1



테스토 사베리스 보고서

kyon410@daum.net
생성됨: 2026년 8월 12일 오전 10시 16분 58초 KST

2026년 8월 1일 오전 12시 0분 0초 KST - 2026년 8월 12일 오후 11시 59분 59초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C) CH1
타입 스탬프 [Asia/Seoul]	온도 [°C]
26. 8. 10. 오후 10:30:00	19.1
26. 8. 10. 오후 11:00:00	19.1
26. 8. 10. 오후 11:30:00	19.2
26. 8. 11. 오전 12:00:00	19.1
26. 8. 11. 오전 12:30:00	19.1
26. 8. 11. 오전 1:00:00	19.2
26. 8. 11. 오전 1:30:00	19.2
26. 8. 11. 오전 2:00:00	19.2
26. 8. 11. 오전 2:30:00	19.1
26. 8. 11. 오전 3:00:00	19.2
26. 8. 11. 오전 3:30:00	19.2
26. 8. 11. 오전 4:00:00	19.2
26. 8. 11. 오전 4:30:00	19.2
26. 8. 11. 오전 5:00:00	19.2
26. 8. 11. 오전 5:30:00	19.2
26. 8. 11. 오전 6:00:00	19.2
26. 8. 11. 오전 6:30:00	19.3
26. 8. 11. 오전 7:00:00	19.3
26. 8. 11. 오전 7:30:00	19.2
26. 8. 11. 오전 8:00:00	19.2
26. 8. 11. 오전 8:30:00	19.2
26. 8. 11. 오전 9:00:00	19.1
26. 8. 11. 오전 9:30:00	19.1
26. 8. 11. 오전 10:00:00	19.2
26. 8. 11. 오전 10:30:00	19.1
26. 8. 11. 오전 11:00:00	19.1
26. 8. 11. 오전 11:30:00	19.1
26. 8. 11. 오후 12:00:00	19.1
26. 8. 11. 오후 12:30:00	19.1
26. 8. 11. 오후 1:00:00	19.1
26. 8. 11. 오후 1:30:00	19.1
26. 8. 11. 오후 2:00:00	19.2



테스토 사베리스 보고서

kyon410@daum.net
생성됨: 2026년 8월 12일 오전 10시 16분 58초 KST

2026년 8월 1일 오전 12시 0분 0초 KST - 2026년 8월 12일 오후 11시 59분 59초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C) CH1
타입 스탬프 [Asia/Seoul]	온도 [°C]
26. 8. 11. 오후 2:30:00	19.2
26. 8. 11. 오후 3:00:00	19.2
26. 8. 11. 오후 3:30:00	19.3
26. 8. 11. 오후 4:00:00	19.4
26. 8. 11. 오후 4:30:00	19.3
26. 8. 11. 오후 5:00:00	19.2
26. 8. 11. 오후 5:30:00	19.2
26. 8. 11. 오후 6:00:00	19.2
26. 8. 11. 오후 6:30:00	19.1
26. 8. 11. 오후 7:00:00	19.1
26. 8. 11. 오후 7:30:00	19.1
26. 8. 11. 오후 8:00:00	19.1
26. 8. 11. 오후 8:30:00	19.1
26. 8. 11. 오후 9:00:00	19.1
26. 8. 11. 오후 9:30:00	19.1
26. 8. 11. 오후 10:00:00	19.1
26. 8. 11. 오후 10:30:00	19.1
26. 8. 11. 오후 11:00:00	19.1
26. 8. 11. 오후 11:30:00	19.1
26. 8. 12. 오전 12:00:00	19.1
26. 8. 12. 오전 12:30:00	19.1
26. 8. 12. 오전 1:00:00	19.1
26. 8. 12. 오전 1:30:00	19.1
26. 8. 12. 오전 2:00:00	19.2
26. 8. 12. 오전 2:30:00	19.1
26. 8. 12. 오전 3:00:00	19.2
26. 8. 12. 오전 3:30:00	19.2
26. 8. 12. 오전 4:00:00	19.2
26. 8. 12. 오전 4:30:00	19.3
26. 8. 12. 오전 5:00:00	19.2
26. 8. 12. 오전 5:30:00	19.1
26. 8. 12. 오전 6:00:00	19.1



테스토 사베리스 보고서

kyon410@daum.net
생성됨: 2026년 8월 12일 오전 10시 16분 58초 KST

2026년 8월 1일 오전 12시 0분 0초 KST - 2026년 8월 12일 오후 11시 59분 59초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(항온고 15~25°C)_CH1
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 8. 12. 오전 6:30:00	19.1
26. 8. 12. 오전 7:00:00	19.1
26. 8. 12. 오전 7:30:00	19.0
26. 8. 12. 오전 8:00:00	19.1
26. 8. 12. 오전 8:30:00	19.1
26. 8. 12. 오전 9:00:00	19.1
26. 8. 12. 오전 9:30:00	19.1



테스토 사베리스 보고서

kyon410@daum.net
생성됨: 2026년 8월 12일 오전 10시 16분 58초 KST

2026년 8월 1일 오전 12시 0분 0초 KST - 2026년 8월 12일 오후 11시 59분 59초 KST

알람 테이블

사용할 수 있는 값이 없습니다.