

IP STORAGE TEMPERATURE LOG

Period of Record	Year. Month	2026.9	
Storage	Location	동국대학교일산병원 임상시험센터 약국	9134Room
	Description		실온보관장
	Manufacturer	NA	
	Model Name	NA	
	Serial Number	NA	
	Temperature Range	1~30℃	
	Setting Point	20~23℃	
Thermometer (including alarm system)	Description	Wi-Fi data logger	
	Manufacturer	Testo Korea	
	Model Name	Saveris2-T3	
	Serial Number	0054704396	
	Sensor Channel	Channel 2	
Thermometer Calibration Certification	Min/Max Record	Monthly report	This month Min (23.0)℃, Max(25.4)℃
	Date of Calibration	2025.11.25	
	Date of Calibration Due	2026.11.25	
Alarm System	Certificate No.	A25-2356-001	
	Alarm threshold	Lower : 10.0℃, Upper : 29.0℃, immediate alarm without delay	
	Alarm Method	Mail, SMS to Pharmacists	


Date	Time	Temperature(℃)	Deviation	Verified by	Comment
2026.09.01	p:20	24.0	☐ Y / ☒ N	ml	
2026.09.02	p:10	24.0	☐ Y / ☒ N	ml	
2026.09.03	p:20	24.0	☐ Y / ☒ N	ml	
2026.09.04	p:20	23.9	☐ Y / ☒ N	ml	
2026.09.05			☐ Y / ☐ N		
2026.09.06			☐ Y / ☐ N		
2026.09.07	p:30	23.9	☐ Y / ☒ N	김지현	(HT-006) TP (HD-CARS) return
2026.09.08	p:20	23.8	☐ Y / ☒ N	ml	
2026.09.09	p:20	23.8	☐ Y / ☒ N	ml	
2026.09.10	p:20	24.3	☐ Y / ☒ N	ml	
2026.09.11	p:20	24.2	☐ Y / ☒ N	ml	
2026.09.12			☐ Y / ☐ N		
2026.09.13			☐ Y / ☐ N		
2026.09.14	p:20	23.9	☐ Y / ☒ N	ml	
2026.09.15	p:20	23.7	☐ Y / ☒ N	ml	
2026.09.16	p:20	23.7	☐ Y / ☒ N	ml	
2026.09.17	p:20	23.8	☐ Y / ☒ N	ml	
2026.09.18	p:20	23.7	☐ Y / ☒ N	ml	
2026.09.19			☐ Y / ☐ N		
2026.09.20			☐ Y / ☐ N		
2026.09.21	p:20	23.7	☐ Y / ☒ N	ml	
2026.09.22	p:25	23.8	☐ Y / ☒ N	ml	
2026.09.23	p:20	23.8	☐ Y / ☒ N	ml	
2026.09.24			☐ Y / ☐ N		추석연휴 추석
2026.09.25			☐ Y / ☐ N		
2026.09.26			☐ Y / ☐ N		
2026.09.27			☐ Y / ☐ N		
2026.09.28	p:20	24.0	☐ Y / ☒ N	ml	
2026.09.29	p:20	23.7	☐ Y / ☒ N	ml	
2026.09.30	p:10	23.7	☐ Y / ☒ N	ml	
			☐ Y / ☐ N		



측정 포인트 요약

동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2

임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2





IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
최소값	23.0
최대값	25.4
평균	23.9
26. 9. 1. 오전 12:00:00	23.8
26. 9. 1. 오전 12:30:00	23.9
26. 9. 1. 오전 1:00:00	23.9
26. 9. 1. 오전 1:30:00	23.9
26. 9. 1. 오전 2:00:00	23.9
26. 9. 1. 오전 2:30:00	23.9
26. 9. 1. 오전 3:00:00	23.9
26. 9. 1. 오전 3:30:00	23.9
26. 9. 1. 오전 4:00:00	23.9
26. 9. 1. 오전 4:30:00	23.9
26. 9. 1. 오전 5:00:00	23.9
26. 9. 1. 오전 5:30:00	23.9
26. 9. 1. 오전 6:00:00	24.0
26. 9. 1. 오전 6:30:00	24.0
26. 9. 1. 오전 7:00:00	23.9
26. 9. 1. 오전 7:30:00	24.0
26. 9. 1. 오전 8:00:00	24.0
26. 9. 1. 오전 8:30:00	24.0
26. 9. 1. 오전 9:00:00	24.0
26. 9. 1. 오전 9:30:00	24.1
26. 9. 1. 오전 10:00:00	24.0
26. 9. 1. 오전 10:30:00	23.9
26. 9. 1. 오전 11:00:00	23.8
26. 9. 1. 오전 11:30:00	23.9
26. 9. 1. 오후 12:00:00	23.9
26. 9. 1. 오후 12:30:00	23.7
26. 9. 1. 오후 1:00:00	23.6
26. 9. 1. 오후 1:30:00	23.5



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 1. 오후 2:00:00	23.5
26. 9. 1. 오후 2:30:00	23.6
26. 9. 1. 오후 3:00:00	23.8
26. 9. 1. 오후 3:30:00	23.8
26. 9. 1. 오후 4:00:00	23.5
26. 9. 1. 오후 4:30:00	23.5
26. 9. 1. 오후 5:00:00	23.5
26. 9. 1. 오후 5:30:00	23.5
26. 9. 1. 오후 6:00:00	23.3
26. 9. 1. 오후 6:30:00	23.2
26. 9. 1. 오후 7:00:00	23.4
26. 9. 1. 오후 7:30:00	23.3
26. 9. 1. 오후 8:00:00	23.3
26. 9. 1. 오후 8:30:00	23.5
26. 9. 1. 오후 9:00:00	23.5
26. 9. 1. 오후 9:30:00	23.6
26. 9. 1. 오후 10:00:00	23.6
26. 9. 1. 오후 10:30:00	23.7
26. 9. 1. 오후 11:00:00	23.7
26. 9. 1. 오후 11:30:00	23.7
26. 9. 2. 오전 12:00:00	23.7
26. 9. 2. 오전 12:30:00	23.7
26. 9. 2. 오전 1:00:00	23.8
26. 9. 2. 오전 1:30:00	23.8
26. 9. 2. 오전 2:00:00	23.8
26. 9. 2. 오전 2:30:00	23.7
26. 9. 2. 오전 3:00:00	23.8
26. 9. 2. 오전 3:30:00	23.8
26. 9. 2. 오전 4:00:00	23.8
26. 9. 2. 오전 4:30:00	23.8
26. 9. 2. 오전 5:00:00	23.8



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 2. 오전 5:30:00	23.8
26. 9. 2. 오전 6:00:00	23.8
26. 9. 2. 오전 6:30:00	23.8
26. 9. 2. 오전 7:00:00	23.8
26. 9. 2. 오전 7:30:00	23.9
26. 9. 2. 오전 8:00:00	23.9
26. 9. 2. 오전 8:30:00	23.9
26. 9. 2. 오전 9:00:00	24.0
26. 9. 2. 오전 9:30:00	24.1
26. 9. 2. 오전 10:00:00	24.0
26. 9. 2. 오전 10:30:00	23.8
26. 9. 2. 오전 11:00:00	23.7
26. 9. 2. 오전 11:30:00	23.6
26. 9. 2. 오후 12:00:00	23.7
26. 9. 2. 오후 12:30:00	23.7
26. 9. 2. 오후 1:00:00	23.6
26. 9. 2. 오후 1:30:00	23.7
26. 9. 2. 오후 2:00:00	23.8
26. 9. 2. 오후 2:30:00	23.9
26. 9. 2. 오후 3:00:00	23.9
26. 9. 2. 오후 3:30:00	23.9
26. 9. 2. 오후 4:00:00	23.8
26. 9. 2. 오후 4:30:00	23.8
26. 9. 2. 오후 5:00:00	23.8
26. 9. 2. 오후 5:30:00	23.8
26. 9. 2. 오후 6:00:00	23.7
26. 9. 2. 오후 6:30:00	23.6
26. 9. 2. 오후 7:00:00	23.6
26. 9. 2. 오후 7:30:00	23.5
26. 9. 2. 오후 8:00:00	23.5
26. 9. 2. 오후 8:30:00	23.4



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 2. 오후 9:00:00	23.3
26. 9. 2. 오후 9:30:00	23.4
26. 9. 2. 오후 10:00:00	23.4
26. 9. 2. 오후 10:30:00	23.4
26. 9. 2. 오후 11:00:00	23.4
26. 9. 2. 오후 11:30:00	23.3
26. 9. 3. 오전 12:00:00	23.2
26. 9. 3. 오전 12:30:00	23.2
26. 9. 3. 오전 1:00:00	23.3
26. 9. 3. 오전 1:30:00	23.4
26. 9. 3. 오전 2:00:00	23.5
26. 9. 3. 오전 2:30:00	23.6
26. 9. 3. 오전 3:00:00	23.6
26. 9. 3. 오전 3:30:00	23.7
26. 9. 3. 오전 4:00:00	23.7
26. 9. 3. 오전 4:30:00	23.7
26. 9. 3. 오전 5:00:00	23.7
26. 9. 3. 오전 5:30:00	23.7
26. 9. 3. 오전 6:00:00	23.7
26. 9. 3. 오전 6:30:00	23.7
26. 9. 3. 오전 7:00:00	23.8
26. 9. 3. 오전 7:30:00	23.8
26. 9. 3. 오전 8:00:00	23.8
26. 9. 3. 오전 8:30:00	23.9
26. 9. 3. 오전 9:00:00	24.0
26. 9. 3. 오전 9:30:00	24.2
26. 9. 3. 오전 10:00:00	24.1
26. 9. 3. 오전 10:30:00	24.0
26. 9. 3. 오전 11:00:00	23.9
26. 9. 3. 오전 11:30:00	23.9
26. 9. 3. 오후 12:00:00	24.0



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 3. 오후 12:30:00	24.1
26. 9. 3. 오후 1:00:00	24.1
26. 9. 3. 오후 1:30:00	24.1
26. 9. 3. 오후 2:00:00	24.1
26. 9. 3. 오후 2:30:00	24.1
26. 9. 3. 오후 3:00:00	24.2
26. 9. 3. 오후 3:30:00	24.1
26. 9. 3. 오후 4:00:00	24.0
26. 9. 3. 오후 4:30:00	23.9
26. 9. 3. 오후 5:00:00	23.8
26. 9. 3. 오후 5:30:00	23.8
26. 9. 3. 오후 6:00:00	23.8
26. 9. 3. 오후 6:30:00	23.7
26. 9. 3. 오후 7:00:00	23.5
26. 9. 3. 오후 7:30:00	23.6
26. 9. 3. 오후 8:00:00	23.4
26. 9. 3. 오후 8:30:00	23.4
26. 9. 3. 오후 9:00:00	23.4
26. 9. 3. 오후 9:30:00	23.4
26. 9. 3. 오후 10:00:00	23.3
26. 9. 3. 오후 10:30:00	23.2
26. 9. 3. 오후 11:00:00	23.3
26. 9. 3. 오후 11:30:00	23.3
26. 9. 4. 오전 12:00:00	23.3
26. 9. 4. 오전 12:30:00	23.3
26. 9. 4. 오전 1:00:00	23.3
26. 9. 4. 오전 1:30:00	23.4
26. 9. 4. 오전 2:00:00	23.5
26. 9. 4. 오전 2:30:00	23.6
26. 9. 4. 오전 3:00:00	23.6
26. 9. 4. 오전 3:30:00	23.6



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 4. 오전 4:00:00	23.6
26. 9. 4. 오전 4:30:00	23.6
26. 9. 4. 오전 5:00:00	23.6
26. 9. 4. 오전 5:30:00	23.6
26. 9. 4. 오전 6:00:00	23.6
26. 9. 4. 오전 6:30:00	23.5
26. 9. 4. 오전 7:00:00	23.6
26. 9. 4. 오전 7:30:00	23.6
26. 9. 4. 오전 8:00:00	23.6
26. 9. 4. 오전 8:30:00	23.6
26. 9. 4. 오전 9:00:00	23.8
26. 9. 4. 오전 9:30:00	23.9
26. 9. 4. 오전 10:00:00	24.1
26. 9. 4. 오전 10:30:00	24.1
26. 9. 4. 오전 11:00:00	24.0
26. 9. 4. 오전 11:30:00	23.9
26. 9. 4. 오후 12:00:00	24.0
26. 9. 4. 오후 12:30:00	23.9
26. 9. 4. 오후 1:00:00	23.9
26. 9. 4. 오후 1:30:00	23.9
26. 9. 4. 오후 2:00:00	23.9
26. 9. 4. 오후 2:30:00	24.0
26. 9. 4. 오후 3:00:00	24.0
26. 9. 4. 오후 3:30:00	23.9
26. 9. 4. 오후 4:00:00	23.9
26. 9. 4. 오후 4:30:00	23.9
26. 9. 4. 오후 5:00:00	23.8
26. 9. 4. 오후 5:30:00	23.7
26. 9. 4. 오후 6:00:00	23.6
26. 9. 4. 오후 6:30:00	23.6
26. 9. 4. 오후 7:00:00	23.4



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 4. 오후 7:30:00	23.5
26. 9. 4. 오후 8:00:00	23.4
26. 9. 4. 오후 8:30:00	23.4
26. 9. 4. 오후 9:00:00	23.3
26. 9. 4. 오후 9:30:00	23.3
26. 9. 4. 오후 10:00:00	23.3
26. 9. 4. 오후 10:30:00	23.4
26. 9. 4. 오후 11:00:00	23.5
26. 9. 4. 오후 11:30:00	23.6
26. 9. 5. 오전 12:00:00	23.6
26. 9. 5. 오전 12:30:00	23.6
26. 9. 5. 오전 1:00:00	23.6
26. 9. 5. 오전 1:30:00	23.7
26. 9. 5. 오전 2:00:00	23.7
26. 9. 5. 오전 2:30:00	23.7
26. 9. 5. 오전 3:00:00	23.6
26. 9. 5. 오전 3:30:00	23.7
26. 9. 5. 오전 4:00:00	23.7
26. 9. 5. 오전 4:30:00	23.6
26. 9. 5. 오전 5:00:00	23.6
26. 9. 5. 오전 5:30:00	23.6
26. 9. 5. 오전 6:00:00	23.6
26. 9. 5. 오전 6:30:00	23.6
26. 9. 5. 오전 7:00:00	23.6
26. 9. 5. 오전 7:30:00	23.6
26. 9. 5. 오전 8:00:00	23.6
26. 9. 5. 오전 8:30:00	23.6
26. 9. 5. 오전 9:00:00	23.8
26. 9. 5. 오전 9:30:00	24.0
26. 9. 5. 오전 10:00:00	24.1
26. 9. 5. 오전 10:30:00	24.0



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 5. 오전 11:00:00	23.9
26. 9. 5. 오전 11:30:00	23.9
26. 9. 5. 오후 12:00:00	23.9
26. 9. 5. 오후 12:30:00	23.9
26. 9. 5. 오후 1:00:00	24.0
26. 9. 5. 오후 1:30:00	24.0
26. 9. 5. 오후 2:00:00	24.0
26. 9. 5. 오후 2:30:00	24.0
26. 9. 5. 오후 3:00:00	23.9
26. 9. 5. 오후 3:30:00	23.9
26. 9. 5. 오후 4:00:00	23.8
26. 9. 5. 오후 4:30:00	23.8
26. 9. 5. 오후 5:00:00	23.7
26. 9. 5. 오후 5:30:00	23.5
26. 9. 5. 오후 6:00:00	23.4
26. 9. 5. 오후 6:30:00	23.5
26. 9. 5. 오후 7:00:00	23.3
26. 9. 5. 오후 7:30:00	23.3
26. 9. 5. 오후 8:00:00	23.3
26. 9. 5. 오후 8:30:00	23.4
26. 9. 5. 오후 9:00:00	23.4
26. 9. 5. 오후 9:30:00	23.4
26. 9. 5. 오후 10:00:00	23.4
26. 9. 5. 오후 10:30:00	23.4
26. 9. 5. 오후 11:00:00	23.5
26. 9. 5. 오후 11:30:00	23.6
26. 9. 6. 오전 12:00:00	23.7
26. 9. 6. 오전 12:30:00	23.7
26. 9. 6. 오전 1:00:00	23.8
26. 9. 6. 오전 1:30:00	23.7
26. 9. 6. 오전 2:00:00	23.8



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 6. 오전 2:30:00	23.8
26. 9. 6. 오전 3:00:00	23.8
26. 9. 6. 오전 3:30:00	23.7
26. 9. 6. 오전 4:00:00	23.7
26. 9. 6. 오전 4:30:00	23.7
26. 9. 6. 오전 5:00:00	23.7
26. 9. 6. 오전 5:30:00	23.8
26. 9. 6. 오전 6:00:00	23.7
26. 9. 6. 오전 6:30:00	23.8
26. 9. 6. 오전 7:00:00	23.7
26. 9. 6. 오전 7:30:00	23.8
26. 9. 6. 오전 8:00:00	23.8
26. 9. 6. 오전 8:30:00	23.8
26. 9. 6. 오전 9:00:00	24.0
26. 9. 6. 오전 9:30:00	24.2
26. 9. 6. 오전 10:00:00	24.1
26. 9. 6. 오전 10:30:00	24.0
26. 9. 6. 오전 11:00:00	23.8
26. 9. 6. 오전 11:30:00	23.9
26. 9. 6. 오후 12:00:00	23.9
26. 9. 6. 오후 12:30:00	24.1
26. 9. 6. 오후 1:00:00	24.0
26. 9. 6. 오후 1:30:00	24.1
26. 9. 6. 오후 2:00:00	24.1
26. 9. 6. 오후 2:30:00	24.1
26. 9. 6. 오후 3:00:00	24.0
26. 9. 6. 오후 3:30:00	23.9
26. 9. 6. 오후 4:00:00	23.8
26. 9. 6. 오후 4:30:00	23.8
26. 9. 6. 오후 5:00:00	23.7
26. 9. 6. 오후 5:30:00	23.6



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 6. 오후 6:00:00	23.6
26. 9. 6. 오후 6:30:00	23.4
26. 9. 6. 오후 7:00:00	23.4
26. 9. 6. 오후 7:30:00	23.4
26. 9. 6. 오후 8:00:00	23.3
26. 9. 6. 오후 8:30:00	23.3
26. 9. 6. 오후 9:00:00	23.3
26. 9. 6. 오후 9:30:00	23.4
26. 9. 6. 오후 10:00:00	23.3
26. 9. 6. 오후 10:30:00	23.3
26. 9. 6. 오후 11:00:00	23.3
26. 9. 6. 오후 11:30:00	23.4
26. 9. 7. 오전 12:00:00	23.4
26. 9. 7. 오전 12:30:00	23.4
26. 9. 7. 오전 1:00:00	23.4
26. 9. 7. 오전 1:30:00	23.4
26. 9. 7. 오전 2:00:00	23.2
26. 9. 7. 오전 2:30:00	23.2
26. 9. 7. 오전 3:00:00	23.4
26. 9. 7. 오전 3:30:00	23.4
26. 9. 7. 오전 4:00:00	23.4
26. 9. 7. 오전 4:30:00	23.5
26. 9. 7. 오전 5:00:00	23.5
26. 9. 7. 오전 5:30:00	23.7
26. 9. 7. 오전 6:00:00	23.7
26. 9. 7. 오전 6:30:00	23.7
26. 9. 7. 오전 7:00:00	23.7
26. 9. 7. 오전 7:30:00	23.8
26. 9. 7. 오전 8:00:00	23.9
26. 9. 7. 오전 8:30:00	24.0
26. 9. 7. 오전 9:00:00	24.0



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 7. 오전 9:30:00	23.9
26. 9. 7. 오전 10:00:00	23.9
26. 9. 7. 오전 10:30:00	23.9
26. 9. 7. 오전 11:00:00	23.9
26. 9. 7. 오전 11:30:00	24.0
26. 9. 7. 오후 12:00:00	24.0
26. 9. 7. 오후 12:30:00	24.0
26. 9. 7. 오후 1:00:00	24.0
26. 9. 7. 오후 1:30:00	24.1
26. 9. 7. 오후 2:00:00	24.1
26. 9. 7. 오후 2:30:00	24.1
26. 9. 7. 오후 3:00:00	24.1
26. 9. 7. 오후 3:30:00	24.0
26. 9. 7. 오후 4:00:00	23.9
26. 9. 7. 오후 4:30:00	23.9
26. 9. 7. 오후 5:00:00	23.8
26. 9. 7. 오후 5:30:00	23.7
26. 9. 7. 오후 6:00:00	23.6
26. 9. 7. 오후 6:30:00	23.6
26. 9. 7. 오후 7:00:00	23.4
26. 9. 7. 오후 7:30:00	23.4
26. 9. 7. 오후 8:00:00	23.4
26. 9. 7. 오후 8:30:00	23.3
26. 9. 7. 오후 9:00:00	23.3
26. 9. 7. 오후 9:30:00	23.3
26. 9. 7. 오후 10:00:00	23.3
26. 9. 7. 오후 10:30:00	23.4
26. 9. 7. 오후 11:00:00	23.4
26. 9. 7. 오후 11:30:00	23.4
26. 9. 8. 오전 12:00:00	23.4
26. 9. 8. 오전 12:30:00	23.4



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 8. 오전 1:00:00	23.3
26. 9. 8. 오전 1:30:00	23.5
26. 9. 8. 오전 2:00:00	23.5
26. 9. 8. 오전 2:30:00	23.6
26. 9. 8. 오전 3:00:00	23.7
26. 9. 8. 오전 3:30:00	23.7
26. 9. 8. 오전 4:00:00	23.7
26. 9. 8. 오전 4:30:00	23.6
26. 9. 8. 오전 5:00:00	23.7
26. 9. 8. 오전 5:30:00	23.7
26. 9. 8. 오전 6:00:00	23.7
26. 9. 8. 오전 6:30:00	23.6
26. 9. 8. 오전 7:00:00	23.6
26. 9. 8. 오전 7:30:00	23.7
26. 9. 8. 오전 8:00:00	23.7
26. 9. 8. 오전 8:30:00	23.7
26. 9. 8. 오전 9:00:00	23.8
26. 9. 8. 오전 9:30:00	24.0
26. 9. 8. 오전 10:00:00	24.0
26. 9. 8. 오전 10:30:00	24.0
26. 9. 8. 오전 11:00:00	23.9
26. 9. 8. 오전 11:30:00	23.9
26. 9. 8. 오후 12:00:00	24.1
26. 9. 8. 오후 12:30:00	24.2
26. 9. 8. 오후 1:00:00	24.2
26. 9. 8. 오후 1:30:00	24.1
26. 9. 8. 오후 2:00:00	24.2
26. 9. 8. 오후 2:30:00	24.2
26. 9. 8. 오후 3:00:00	24.0
26. 9. 8. 오후 3:30:00	24.0
26. 9. 8. 오후 4:00:00	23.9



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 8. 오후 4:30:00	23.7
26. 9. 8. 오후 5:00:00	23.7
26. 9. 8. 오후 5:30:00	23.7
26. 9. 8. 오후 6:00:00	23.6
26. 9. 8. 오후 6:30:00	23.5
26. 9. 8. 오후 7:00:00	23.5
26. 9. 8. 오후 7:30:00	23.3
26. 9. 8. 오후 8:00:00	23.4
26. 9. 8. 오후 8:30:00	23.4
26. 9. 8. 오후 9:00:00	23.4
26. 9. 8. 오후 9:30:00	23.4
26. 9. 8. 오후 10:00:00	23.3
26. 9. 8. 오후 10:30:00	23.3
26. 9. 8. 오후 11:00:00	23.3
26. 9. 8. 오후 11:30:00	23.3
26. 9. 9. 오전 12:00:00	23.3
26. 9. 9. 오전 12:30:00	23.4
26. 9. 9. 오전 1:00:00	23.4
26. 9. 9. 오전 1:30:00	23.3
26. 9. 9. 오전 2:00:00	23.3
26. 9. 9. 오전 2:30:00	23.3
26. 9. 9. 오전 3:00:00	23.5
26. 9. 9. 오전 3:30:00	23.5
26. 9. 9. 오전 4:00:00	23.6
26. 9. 9. 오전 4:30:00	23.6
26. 9. 9. 오전 5:00:00	23.6
26. 9. 9. 오전 5:30:00	23.6
26. 9. 9. 오전 6:00:00	23.6
26. 9. 9. 오전 6:30:00	23.7
26. 9. 9. 오전 7:00:00	23.6
26. 9. 9. 오전 7:30:00	23.7



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 9. 오전 8:00:00	23.7
26. 9. 9. 오전 8:30:00	23.7
26. 9. 9. 오전 9:00:00	23.8
26. 9. 9. 오전 9:30:00	23.9
26. 9. 9. 오전 10:00:00	24.1
26. 9. 9. 오전 10:30:00	24.2
26. 9. 9. 오전 11:00:00	24.3
26. 9. 9. 오전 11:30:00	24.5
26. 9. 9. 오후 12:00:00	24.6
26. 9. 9. 오후 12:30:00	24.5
26. 9. 9. 오후 1:00:00	24.5
26. 9. 9. 오후 1:30:00	24.3
26. 9. 9. 오후 2:00:00	24.2
26. 9. 9. 오후 2:30:00	24.1
26. 9. 9. 오후 3:00:00	24.1
26. 9. 9. 오후 3:30:00	24.0
26. 9. 9. 오후 4:00:00	24.0
26. 9. 9. 오후 4:30:00	23.9
26. 9. 9. 오후 5:00:00	23.9
26. 9. 9. 오후 5:30:00	23.7
26. 9. 9. 오후 6:00:00	23.8
26. 9. 9. 오후 6:30:00	23.6
26. 9. 9. 오후 7:00:00	23.5
26. 9. 9. 오후 7:30:00	23.4
26. 9. 9. 오후 8:00:00	23.4
26. 9. 9. 오후 8:30:00	23.5
26. 9. 9. 오후 9:00:00	23.6
26. 9. 9. 오후 9:30:00	23.6
26. 9. 9. 오후 10:00:00	23.6
26. 9. 9. 오후 10:30:00	23.5
26. 9. 9. 오후 11:00:00	23.4



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 9. 오후 11:30:00	23.6
26. 9. 10. 오전 12:00:00	23.6
26. 9. 10. 오전 12:30:00	23.5
26. 9. 10. 오전 1:00:00	23.4
26. 9. 10. 오전 1:30:00	23.4
26. 9. 10. 오전 2:00:00	23.5
26. 9. 10. 오전 2:30:00	23.5
26. 9. 10. 오전 3:00:00	23.7
26. 9. 10. 오전 3:30:00	23.8
26. 9. 10. 오전 4:00:00	23.8
26. 9. 10. 오전 4:30:00	23.9
26. 9. 10. 오전 5:00:00	23.9
26. 9. 10. 오전 5:30:00	23.9
26. 9. 10. 오전 6:00:00	24.0
26. 9. 10. 오전 6:30:00	24.1
26. 9. 10. 오전 7:00:00	24.1
26. 9. 10. 오전 7:30:00	24.2
26. 9. 10. 오전 8:00:00	24.3
26. 9. 10. 오전 8:30:00	24.3
26. 9. 10. 오전 9:00:00	24.3
26. 9. 10. 오전 9:30:00	24.2
26. 9. 10. 오전 10:00:00	24.0
26. 9. 10. 오전 10:30:00	23.8
26. 9. 10. 오전 11:00:00	23.8
26. 9. 10. 오전 11:30:00	23.7
26. 9. 10. 오후 12:00:00	23.6
26. 9. 10. 오후 12:30:00	23.6
26. 9. 10. 오후 1:00:00	23.7
26. 9. 10. 오후 1:30:00	23.8
26. 9. 10. 오후 2:00:00	23.9
26. 9. 10. 오후 2:30:00	23.9



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 10. 오후 3:00:00	23.9
26. 9. 10. 오후 3:30:00	23.7
26. 9. 10. 오후 4:00:00	23.9
26. 9. 10. 오후 4:30:00	23.8
26. 9. 10. 오후 5:00:00	23.6
26. 9. 10. 오후 5:30:00	23.7
26. 9. 10. 오후 6:00:00	23.5
26. 9. 10. 오후 6:30:00	23.6
26. 9. 10. 오후 7:00:00	23.5
26. 9. 10. 오후 7:30:00	23.6
26. 9. 10. 오후 8:00:00	23.4
26. 9. 10. 오후 8:30:00	23.5
26. 9. 10. 오후 9:00:00	23.6
26. 9. 10. 오후 9:30:00	23.7
26. 9. 10. 오후 10:00:00	23.6
26. 9. 10. 오후 10:30:00	23.6
26. 9. 10. 오후 11:00:00	23.6
26. 9. 10. 오후 11:30:00	23.7
26. 9. 11. 오전 12:00:00	23.8
26. 9. 11. 오전 12:30:00	23.8
26. 9. 11. 오전 1:00:00	23.8
26. 9. 11. 오전 1:30:00	23.9
26. 9. 11. 오전 2:00:00	24.0
26. 9. 11. 오전 2:30:00	24.0
26. 9. 11. 오전 3:00:00	24.0
26. 9. 11. 오전 3:30:00	23.9
26. 9. 11. 오전 4:00:00	23.9
26. 9. 11. 오전 4:30:00	23.9
26. 9. 11. 오전 5:00:00	23.8
26. 9. 11. 오전 5:30:00	23.8
26. 9. 11. 오전 6:00:00	23.9



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 11. 오전 6:30:00	23.9
26. 9. 11. 오전 7:00:00	23.9
26. 9. 11. 오전 7:30:00	23.8
26. 9. 11. 오전 8:00:00	23.9
26. 9. 11. 오전 8:30:00	24.1
26. 9. 11. 오전 9:00:00	24.2
26. 9. 11. 오전 9:30:00	24.3
26. 9. 11. 오전 10:00:00	24.2
26. 9. 11. 오전 10:30:00	24.3
26. 9. 11. 오전 11:00:00	24.6
26. 9. 11. 오전 11:30:00	24.6
26. 9. 11. 오후 12:00:00	24.8
26. 9. 11. 오후 12:30:00	25.0
26. 9. 11. 오후 1:00:00	25.0
26. 9. 11. 오후 1:30:00	25.1
26. 9. 11. 오후 2:00:00	25.0
26. 9. 11. 오후 2:30:00	24.9
26. 9. 11. 오후 3:00:00	24.9
26. 9. 11. 오후 3:30:00	24.7
26. 9. 11. 오후 4:00:00	24.5
26. 9. 11. 오후 4:30:00	24.4
26. 9. 11. 오후 5:00:00	24.4
26. 9. 11. 오후 5:30:00	24.3
26. 9. 11. 오후 6:00:00	24.2
26. 9. 11. 오후 6:30:00	24.1
26. 9. 11. 오후 7:00:00	23.9
26. 9. 11. 오후 7:30:00	23.7
26. 9. 11. 오후 8:00:00	23.8
26. 9. 11. 오후 8:30:00	23.7
26. 9. 11. 오후 9:00:00	23.7
26. 9. 11. 오후 9:30:00	23.6



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 11. 오후 10:00:00	23.7
26. 9. 11. 오후 10:30:00	23.6
26. 9. 11. 오후 11:00:00	23.7
26. 9. 11. 오후 11:30:00	23.5
26. 9. 12. 오전 12:00:00	23.6
26. 9. 12. 오전 12:30:00	23.7
26. 9. 12. 오전 1:00:00	23.6
26. 9. 12. 오전 1:30:00	23.4
26. 9. 12. 오전 2:00:00	23.5
26. 9. 12. 오전 2:30:00	23.4
26. 9. 12. 오전 3:00:00	23.5
26. 9. 12. 오전 3:30:00	23.6
26. 9. 12. 오전 4:00:00	23.6
26. 9. 12. 오전 4:30:00	23.7
26. 9. 12. 오전 5:00:00	23.6
26. 9. 12. 오전 5:30:00	23.4
26. 9. 12. 오전 6:00:00	23.4
26. 9. 12. 오전 6:30:00	23.4
26. 9. 12. 오전 7:00:00	23.4
26. 9. 12. 오전 7:30:00	23.4
26. 9. 12. 오전 8:00:00	23.6
26. 9. 12. 오전 8:30:00	23.6
26. 9. 12. 오전 9:00:00	23.5
26. 9. 12. 오전 9:30:00	23.7
26. 9. 12. 오전 10:00:00	23.7
26. 9. 12. 오전 10:30:00	23.7
26. 9. 12. 오전 11:00:00	23.8
26. 9. 12. 오전 11:30:00	23.8
26. 9. 12. 오후 12:00:00	23.8
26. 9. 12. 오후 12:30:00	24.0
26. 9. 12. 오후 1:00:00	24.0



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 12. 오후 1:30:00	23.9
26. 9. 12. 오후 2:00:00	23.8
26. 9. 12. 오후 2:30:00	23.8
26. 9. 12. 오후 3:00:00	23.8
26. 9. 12. 오후 3:30:00	23.8
26. 9. 12. 오후 4:00:00	23.9
26. 9. 12. 오후 4:30:00	23.9
26. 9. 12. 오후 5:00:00	23.8
26. 9. 12. 오후 5:30:00	24.0
26. 9. 12. 오후 6:00:00	23.9
26. 9. 12. 오후 6:30:00	23.8
26. 9. 12. 오후 7:00:00	23.8
26. 9. 12. 오후 7:30:00	23.7
26. 9. 12. 오후 8:00:00	23.6
26. 9. 12. 오후 8:30:00	23.6
26. 9. 12. 오후 9:00:00	23.6
26. 9. 12. 오후 9:30:00	23.7
26. 9. 12. 오후 10:00:00	23.5
26. 9. 12. 오후 10:30:00	23.6
26. 9. 12. 오후 11:00:00	23.7
26. 9. 12. 오후 11:30:00	23.6
26. 9. 13. 오전 12:00:00	23.6
26. 9. 13. 오전 12:30:00	23.7
26. 9. 13. 오전 1:00:00	23.6
26. 9. 13. 오전 1:30:00	23.5
26. 9. 13. 오전 2:00:00	23.6
26. 9. 13. 오전 2:30:00	23.6
26. 9. 13. 오전 3:00:00	23.5
26. 9. 13. 오전 3:30:00	23.5
26. 9. 13. 오전 4:00:00	23.5
26. 9. 13. 오전 4:30:00	23.6



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 13. 오전 5:00:00	23.5
26. 9. 13. 오전 5:30:00	23.4
26. 9. 13. 오전 6:00:00	23.4
26. 9. 13. 오전 6:30:00	23.6
26. 9. 13. 오전 7:00:00	23.6
26. 9. 13. 오전 7:30:00	23.5
26. 9. 13. 오전 8:00:00	23.4
26. 9. 13. 오전 8:30:00	23.6
26. 9. 13. 오전 9:00:00	23.8
26. 9. 13. 오전 9:30:00	23.8
26. 9. 13. 오전 10:00:00	23.7
26. 9. 13. 오전 10:30:00	23.8
26. 9. 13. 오전 11:00:00	23.9
26. 9. 13. 오전 11:30:00	23.9
26. 9. 13. 오후 12:00:00	24.0
26. 9. 13. 오후 12:30:00	23.9
26. 9. 13. 오후 1:00:00	23.8
26. 9. 13. 오후 1:30:00	23.9
26. 9. 13. 오후 2:00:00	23.8
26. 9. 13. 오후 2:30:00	23.8
26. 9. 13. 오후 3:00:00	23.9
26. 9. 13. 오후 3:30:00	23.9
26. 9. 13. 오후 4:00:00	24.0
26. 9. 13. 오후 4:30:00	24.0
26. 9. 13. 오후 5:00:00	23.9
26. 9. 13. 오후 5:30:00	24.0
26. 9. 13. 오후 6:00:00	23.9
26. 9. 13. 오후 6:30:00	23.8
26. 9. 13. 오후 7:00:00	23.8
26. 9. 13. 오후 7:30:00	23.8
26. 9. 13. 오후 8:00:00	23.9



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 13. 오후 8:30:00	23.8
26. 9. 13. 오후 9:00:00	23.8
26. 9. 13. 오후 9:30:00	23.8
26. 9. 13. 오후 10:00:00	23.7
26. 9. 13. 오후 10:30:00	23.8
26. 9. 13. 오후 11:00:00	23.7
26. 9. 13. 오후 11:30:00	23.8
26. 9. 14. 오전 12:00:00	23.7
26. 9. 14. 오전 12:30:00	23.8
26. 9. 14. 오전 1:00:00	23.8
26. 9. 14. 오전 1:30:00	23.7
26. 9. 14. 오전 2:00:00	23.7
26. 9. 14. 오전 2:30:00	23.7
26. 9. 14. 오전 3:00:00	23.5
26. 9. 14. 오전 3:30:00	23.7
26. 9. 14. 오전 4:00:00	23.8
26. 9. 14. 오전 4:30:00	23.6
26. 9. 14. 오전 5:00:00	23.6
26. 9. 14. 오전 5:30:00	23.8
26. 9. 14. 오전 6:00:00	23.7
26. 9. 14. 오전 6:30:00	23.7
26. 9. 14. 오전 7:00:00	23.9
26. 9. 14. 오전 7:30:00	23.8
26. 9. 14. 오전 8:00:00	23.8
26. 9. 14. 오전 8:30:00	23.7
26. 9. 14. 오전 9:00:00	23.9
26. 9. 14. 오전 9:30:00	23.9
26. 9. 14. 오전 10:00:00	24.1
26. 9. 14. 오전 10:30:00	24.2
26. 9. 14. 오전 11:00:00	24.4
26. 9. 14. 오전 11:30:00	24.6



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 14. 오후 12:00:00	24.7
26. 9. 14. 오후 12:30:00	24.9
26. 9. 14. 오후 1:00:00	24.9
26. 9. 14. 오후 1:30:00	25.2
26. 9. 14. 오후 2:00:00	25.3
26. 9. 14. 오후 2:30:00	25.3
26. 9. 14. 오후 3:00:00	25.2
26. 9. 14. 오후 3:30:00	25.0
26. 9. 14. 오후 4:00:00	24.9
26. 9. 14. 오후 4:30:00	24.8
26. 9. 14. 오후 5:00:00	24.8
26. 9. 14. 오후 5:30:00	24.6
26. 9. 14. 오후 6:00:00	24.4
26. 9. 14. 오후 6:30:00	24.2
26. 9. 14. 오후 7:00:00	24.2
26. 9. 14. 오후 7:30:00	24.2
26. 9. 14. 오후 8:00:00	24.1
26. 9. 14. 오후 8:30:00	24.0
26. 9. 14. 오후 9:00:00	23.8
26. 9. 14. 오후 9:30:00	23.8
26. 9. 14. 오후 10:00:00	23.8
26. 9. 14. 오후 10:30:00	23.9
26. 9. 14. 오후 11:00:00	23.7
26. 9. 14. 오후 11:30:00	23.7
26. 9. 15. 오전 12:00:00	23.6
26. 9. 15. 오전 12:30:00	23.5
26. 9. 15. 오전 1:00:00	23.5
26. 9. 15. 오전 1:30:00	23.6
26. 9. 15. 오전 2:00:00	23.5
26. 9. 15. 오전 2:30:00	23.6
26. 9. 15. 오전 3:00:00	23.6



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 15. 오전 3:30:00	23.5
26. 9. 15. 오전 4:00:00	23.5
26. 9. 15. 오전 4:30:00	23.5
26. 9. 15. 오전 5:00:00	23.6
26. 9. 15. 오전 5:30:00	23.6
26. 9. 15. 오전 6:00:00	23.6
26. 9. 15. 오전 6:30:00	23.6
26. 9. 15. 오전 7:00:00	23.7
26. 9. 15. 오전 7:30:00	23.7
26. 9. 15. 오전 8:00:00	23.6
26. 9. 15. 오전 8:30:00	23.7
26. 9. 15. 오전 9:00:00	23.7
26. 9. 15. 오전 9:30:00	23.8
26. 9. 15. 오전 10:00:00	24.1
26. 9. 15. 오전 10:30:00	24.5
26. 9. 15. 오전 11:00:00	24.7
26. 9. 15. 오전 11:30:00	24.9
26. 9. 15. 오후 12:00:00	25.1
26. 9. 15. 오후 12:30:00	25.2
26. 9. 15. 오후 1:00:00	25.2
26. 9. 15. 오후 1:30:00	25.2
26. 9. 15. 오후 2:00:00	25.2
26. 9. 15. 오후 2:30:00	25.3
26. 9. 15. 오후 3:00:00	25.2
26. 9. 15. 오후 3:30:00	25.1
26. 9. 15. 오후 4:00:00	25.0
26. 9. 15. 오후 4:30:00	24.8
26. 9. 15. 오후 5:00:00	24.7
26. 9. 15. 오후 5:30:00	24.6
26. 9. 15. 오후 6:00:00	24.4
26. 9. 15. 오후 6:30:00	24.4



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 15. 오후 7:00:00	24.2
26. 9. 15. 오후 7:30:00	24.1
26. 9. 15. 오후 8:00:00	24.1
26. 9. 15. 오후 8:30:00	24.0
26. 9. 15. 오후 9:00:00	23.9
26. 9. 15. 오후 9:30:00	23.9
26. 9. 15. 오후 10:00:00	23.9
26. 9. 15. 오후 10:30:00	23.7
26. 9. 15. 오후 11:00:00	23.7
26. 9. 15. 오후 11:30:00	23.9
26. 9. 16. 오전 12:00:00	23.8
26. 9. 16. 오전 12:30:00	23.8
26. 9. 16. 오전 1:00:00	23.8
26. 9. 16. 오전 1:30:00	23.7
26. 9. 16. 오전 2:00:00	23.7
26. 9. 16. 오전 2:30:00	23.5
26. 9. 16. 오전 3:00:00	23.6
26. 9. 16. 오전 3:30:00	23.4
26. 9. 16. 오전 4:00:00	23.6
26. 9. 16. 오전 4:30:00	23.7
26. 9. 16. 오전 5:00:00	23.5
26. 9. 16. 오전 5:30:00	23.7
26. 9. 16. 오전 6:00:00	23.6
26. 9. 16. 오전 6:30:00	23.5
26. 9. 16. 오전 7:00:00	23.6
26. 9. 16. 오전 7:30:00	23.5
26. 9. 16. 오전 8:00:00	23.5
26. 9. 16. 오전 8:30:00	23.6
26. 9. 16. 오전 9:00:00	23.7
26. 9. 16. 오전 9:30:00	23.9
26. 9. 16. 오전 10:00:00	24.1



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 16. 오전 10:30:00	24.2
26. 9. 16. 오전 11:00:00	24.4
26. 9. 16. 오전 11:30:00	24.5
26. 9. 16. 오후 12:00:00	24.8
26. 9. 16. 오후 12:30:00	24.9
26. 9. 16. 오후 1:00:00	25.1
26. 9. 16. 오후 1:30:00	25.3
26. 9. 16. 오후 2:00:00	25.3
26. 9. 16. 오후 2:30:00	25.1
26. 9. 16. 오후 3:00:00	25.0
26. 9. 16. 오후 3:30:00	25.1
26. 9. 16. 오후 4:00:00	25.1
26. 9. 16. 오후 4:30:00	24.9
26. 9. 16. 오후 5:00:00	24.6
26. 9. 16. 오후 5:30:00	24.5
26. 9. 16. 오후 6:00:00	24.4
26. 9. 16. 오후 6:30:00	24.2
26. 9. 16. 오후 7:00:00	24.0
26. 9. 16. 오후 7:30:00	23.9
26. 9. 16. 오후 8:00:00	23.8
26. 9. 16. 오후 8:30:00	23.8
26. 9. 16. 오후 9:00:00	23.8
26. 9. 16. 오후 9:30:00	23.7
26. 9. 16. 오후 10:00:00	23.7
26. 9. 16. 오후 10:30:00	23.8
26. 9. 16. 오후 11:00:00	23.6
26. 9. 16. 오후 11:30:00	23.7
26. 9. 17. 오전 12:00:00	23.7
26. 9. 17. 오전 12:30:00	23.7
26. 9. 17. 오전 1:00:00	23.6
26. 9. 17. 오전 1:30:00	23.6



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 17. 오전 2:00:00	23.6
26. 9. 17. 오전 2:30:00	23.5
26. 9. 17. 오전 3:00:00	23.5
26. 9. 17. 오전 3:30:00	23.5
26. 9. 17. 오전 4:00:00	23.5
26. 9. 17. 오전 4:30:00	23.6
26. 9. 17. 오전 5:00:00	23.7
26. 9. 17. 오전 5:30:00	23.6
26. 9. 17. 오전 6:00:00	23.7
26. 9. 17. 오전 6:30:00	23.6
26. 9. 17. 오전 7:00:00	23.7
26. 9. 17. 오전 7:30:00	23.5
26. 9. 17. 오전 8:00:00	23.4
26. 9. 17. 오전 8:30:00	23.6
26. 9. 17. 오전 9:00:00	23.8
26. 9. 17. 오전 9:30:00	24.0
26. 9. 17. 오전 10:00:00	24.0
26. 9. 17. 오전 10:30:00	24.0
26. 9. 17. 오전 11:00:00	24.3
26. 9. 17. 오전 11:30:00	24.7
26. 9. 17. 오후 12:00:00	24.9
26. 9. 17. 오후 12:30:00	24.9
26. 9. 17. 오후 1:00:00	25.3
26. 9. 17. 오후 1:30:00	25.3
26. 9. 17. 오후 2:00:00	25.3
26. 9. 17. 오후 2:30:00	25.2
26. 9. 17. 오후 3:00:00	24.9
26. 9. 17. 오후 3:30:00	24.6
26. 9. 17. 오후 4:00:00	24.6
26. 9. 17. 오후 4:30:00	24.6
26. 9. 17. 오후 5:00:00	24.4



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 17. 오후 5:30:00	24.3
26. 9. 17. 오후 6:00:00	24.1
26. 9. 17. 오후 6:30:00	24.0
26. 9. 17. 오후 7:00:00	24.0
26. 9. 17. 오후 7:30:00	23.8
26. 9. 17. 오후 8:00:00	23.8
26. 9. 17. 오후 8:30:00	23.8
26. 9. 17. 오후 9:00:00	23.8
26. 9. 17. 오후 9:30:00	23.8
26. 9. 17. 오후 10:00:00	23.7
26. 9. 17. 오후 10:30:00	23.7
26. 9. 17. 오후 11:00:00	23.8
26. 9. 17. 오후 11:30:00	23.8
26. 9. 18. 오전 12:00:00	23.7
26. 9. 18. 오전 12:30:00	23.6
26. 9. 18. 오전 1:00:00	23.6
26. 9. 18. 오전 1:30:00	23.6
26. 9. 18. 오전 2:00:00	23.6
26. 9. 18. 오전 2:30:00	23.5
26. 9. 18. 오전 3:00:00	23.4
26. 9. 18. 오전 3:30:00	23.5
26. 9. 18. 오전 4:00:00	23.5
26. 9. 18. 오전 4:30:00	23.7
26. 9. 18. 오전 5:00:00	23.6
26. 9. 18. 오전 5:30:00	23.5
26. 9. 18. 오전 6:00:00	23.6
26. 9. 18. 오전 6:30:00	23.4
26. 9. 18. 오전 7:00:00	23.5
26. 9. 18. 오전 7:30:00	23.5
26. 9. 18. 오전 8:00:00	23.4
26. 9. 18. 오전 8:30:00	23.6



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 18. 오전 9:00:00	23.7
26. 9. 18. 오전 9:30:00	23.8
26. 9. 18. 오전 10:00:00	24.0
26. 9. 18. 오전 10:30:00	24.3
26. 9. 18. 오전 11:00:00	24.5
26. 9. 18. 오전 11:30:00	24.8
26. 9. 18. 오후 12:00:00	25.1
26. 9. 18. 오후 12:30:00	25.2
26. 9. 18. 오후 1:00:00	25.3
26. 9. 18. 오후 1:30:00	25.2
26. 9. 18. 오후 2:00:00	25.3
26. 9. 18. 오후 2:30:00	25.4
26. 9. 18. 오후 3:00:00	25.3
26. 9. 18. 오후 3:30:00	25.4
26. 9. 18. 오후 4:00:00	25.3
26. 9. 18. 오후 4:30:00	25.0
26. 9. 18. 오후 5:00:00	24.7
26. 9. 18. 오후 5:30:00	24.6
26. 9. 18. 오후 6:00:00	24.4
26. 9. 18. 오후 6:30:00	24.3
26. 9. 18. 오후 7:00:00	24.2
26. 9. 18. 오후 7:30:00	24.1
26. 9. 18. 오후 8:00:00	24.1
26. 9. 18. 오후 8:30:00	24.0
26. 9. 18. 오후 9:00:00	23.9
26. 9. 18. 오후 9:30:00	23.9
26. 9. 18. 오후 10:00:00	23.8
26. 9. 18. 오후 10:30:00	23.8
26. 9. 18. 오후 11:00:00	23.9
26. 9. 18. 오후 11:30:00	23.9
26. 9. 19. 오전 12:00:00	23.7



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 19. 오전 12:30:00	23.6
26. 9. 19. 오전 1:00:00	23.7
26. 9. 19. 오전 1:30:00	23.7
26. 9. 19. 오전 2:00:00	23.7
26. 9. 19. 오전 2:30:00	23.6
26. 9. 19. 오전 3:00:00	23.5
26. 9. 19. 오전 3:30:00	23.5
26. 9. 19. 오전 4:00:00	23.6
26. 9. 19. 오전 4:30:00	23.5
26. 9. 19. 오전 5:00:00	23.5
26. 9. 19. 오전 5:30:00	23.6
26. 9. 19. 오전 6:00:00	23.6
26. 9. 19. 오전 6:30:00	23.5
26. 9. 19. 오전 7:00:00	23.6
26. 9. 19. 오전 7:30:00	23.5
26. 9. 19. 오전 8:00:00	23.5
26. 9. 19. 오전 8:30:00	23.7
26. 9. 19. 오전 9:00:00	23.8
26. 9. 19. 오전 9:30:00	23.7
26. 9. 19. 오전 10:00:00	23.8
26. 9. 19. 오전 10:30:00	23.9
26. 9. 19. 오전 11:00:00	24.2
26. 9. 19. 오전 11:30:00	24.6
26. 9. 19. 오후 12:00:00	25.0
26. 9. 19. 오후 12:30:00	25.0
26. 9. 19. 오후 1:00:00	25.2
26. 9. 19. 오후 1:30:00	25.2
26. 9. 19. 오후 2:00:00	25.2
26. 9. 19. 오후 2:30:00	25.1
26. 9. 19. 오후 3:00:00	25.1
26. 9. 19. 오후 3:30:00	24.9



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 19. 오후 4:00:00	24.9
26. 9. 19. 오후 4:30:00	24.8
26. 9. 19. 오후 5:00:00	24.6
26. 9. 19. 오후 5:30:00	24.5
26. 9. 19. 오후 6:00:00	24.3
26. 9. 19. 오후 6:30:00	24.3
26. 9. 19. 오후 7:00:00	24.3
26. 9. 19. 오후 7:30:00	24.1
26. 9. 19. 오후 8:00:00	24.1
26. 9. 19. 오후 8:30:00	24.0
26. 9. 19. 오후 9:00:00	24.0
26. 9. 19. 오후 9:30:00	24.0
26. 9. 19. 오후 10:00:00	24.0
26. 9. 19. 오후 10:30:00	24.0
26. 9. 19. 오후 11:00:00	23.9
26. 9. 19. 오후 11:30:00	23.9
26. 9. 20. 오전 12:00:00	23.9
26. 9. 20. 오전 12:30:00	23.8
26. 9. 20. 오전 1:00:00	23.8
26. 9. 20. 오전 1:30:00	23.8
26. 9. 20. 오전 2:00:00	23.8
26. 9. 20. 오전 2:30:00	23.7
26. 9. 20. 오전 3:00:00	23.7
26. 9. 20. 오전 3:30:00	23.7
26. 9. 20. 오전 4:00:00	23.6
26. 9. 20. 오전 4:30:00	23.6
26. 9. 20. 오전 5:00:00	23.6
26. 9. 20. 오전 5:30:00	23.5
26. 9. 20. 오전 6:00:00	23.6
26. 9. 20. 오전 6:30:00	23.5
26. 9. 20. 오전 7:00:00	23.6



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 20. 오전 7:30:00	23.5
26. 9. 20. 오전 8:00:00	23.6
26. 9. 20. 오전 8:30:00	23.6
26. 9. 20. 오전 9:00:00	23.6
26. 9. 20. 오전 9:30:00	23.7
26. 9. 20. 오전 10:00:00	23.9
26. 9. 20. 오전 10:30:00	24.1
26. 9. 20. 오전 11:00:00	24.4
26. 9. 20. 오전 11:30:00	24.7
26. 9. 20. 오후 12:00:00	24.8
26. 9. 20. 오후 12:30:00	25.1
26. 9. 20. 오후 1:00:00	25.3
26. 9. 20. 오후 1:30:00	25.3
26. 9. 20. 오후 2:00:00	25.3
26. 9. 20. 오후 2:30:00	25.3
26. 9. 20. 오후 3:00:00	25.3
26. 9. 20. 오후 3:30:00	25.2
26. 9. 20. 오후 4:00:00	25.0
26. 9. 20. 오후 4:30:00	24.9
26. 9. 20. 오후 5:00:00	24.8
26. 9. 20. 오후 5:30:00	24.7
26. 9. 20. 오후 6:00:00	24.4
26. 9. 20. 오후 6:30:00	24.3
26. 9. 20. 오후 7:00:00	24.3
26. 9. 20. 오후 7:30:00	24.1
26. 9. 20. 오후 8:00:00	24.2
26. 9. 20. 오후 8:30:00	24.2
26. 9. 20. 오후 9:00:00	24.1
26. 9. 20. 오후 9:30:00	24.1
26. 9. 20. 오후 10:00:00	23.9
26. 9. 20. 오후 10:30:00	23.9



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 20. 오후 11:00:00	24.0
26. 9. 20. 오후 11:30:00	23.9
26. 9. 21. 오전 12:00:00	24.0
26. 9. 21. 오전 12:30:00	23.9
26. 9. 21. 오전 1:00:00	23.8
26. 9. 21. 오전 1:30:00	23.7
26. 9. 21. 오전 2:00:00	23.8
26. 9. 21. 오전 2:30:00	23.6
26. 9. 21. 오전 3:00:00	23.7
26. 9. 21. 오전 3:30:00	23.7
26. 9. 21. 오전 4:00:00	23.7
26. 9. 21. 오전 4:30:00	23.6
26. 9. 21. 오전 5:00:00	23.5
26. 9. 21. 오전 5:30:00	23.6
26. 9. 21. 오전 6:00:00	23.5
26. 9. 21. 오전 6:30:00	23.5
26. 9. 21. 오전 7:00:00	23.4
26. 9. 21. 오전 7:30:00	23.5
26. 9. 21. 오전 8:00:00	23.6
26. 9. 21. 오전 8:30:00	23.5
26. 9. 21. 오전 9:00:00	23.7
26. 9. 21. 오전 9:30:00	23.8
26. 9. 21. 오전 10:00:00	23.9
26. 9. 21. 오전 10:30:00	24.3
26. 9. 21. 오전 11:00:00	24.6
26. 9. 21. 오전 11:30:00	24.6
26. 9. 21. 오후 12:00:00	24.8
26. 9. 21. 오후 12:30:00	24.9
26. 9. 21. 오후 1:00:00	25.0
26. 9. 21. 오후 1:30:00	25.2
26. 9. 21. 오후 2:00:00	25.2



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 21. 오후 2:30:00	25.0
26. 9. 21. 오후 3:00:00	24.9
26. 9. 21. 오후 3:30:00	24.9
26. 9. 21. 오후 4:00:00	24.9
26. 9. 21. 오후 4:30:00	24.6
26. 9. 21. 오후 5:00:00	24.4
26. 9. 21. 오후 5:30:00	24.3
26. 9. 21. 오후 6:00:00	24.3
26. 9. 21. 오후 6:30:00	24.2
26. 9. 21. 오후 7:00:00	24.2
26. 9. 21. 오후 7:30:00	24.2
26. 9. 21. 오후 8:00:00	24.1
26. 9. 21. 오후 8:30:00	23.9
26. 9. 21. 오후 9:00:00	23.9
26. 9. 21. 오후 9:30:00	23.9
26. 9. 21. 오후 10:00:00	23.9
26. 9. 21. 오후 10:30:00	23.9
26. 9. 21. 오후 11:00:00	23.9
26. 9. 21. 오후 11:30:00	23.9
26. 9. 22. 오전 12:00:00	24.0
26. 9. 22. 오전 12:30:00	23.8
26. 9. 22. 오전 1:00:00	23.7
26. 9. 22. 오전 1:30:00	23.8
26. 9. 22. 오전 2:00:00	23.6
26. 9. 22. 오전 2:30:00	23.7
26. 9. 22. 오전 3:00:00	23.5
26. 9. 22. 오전 3:30:00	23.7
26. 9. 22. 오전 4:00:00	23.7
26. 9. 22. 오전 4:30:00	23.6
26. 9. 22. 오전 5:00:00	23.6
26. 9. 22. 오전 5:30:00	23.5



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 22. 오전 6:00:00	23.6
26. 9. 22. 오전 6:30:00	23.6
26. 9. 22. 오전 7:00:00	23.6
26. 9. 22. 오전 7:30:00	23.6
26. 9. 22. 오전 8:00:00	23.6
26. 9. 22. 오전 8:30:00	23.8
26. 9. 22. 오전 9:00:00	23.7
26. 9. 22. 오전 9:30:00	23.8
26. 9. 22. 오전 10:00:00	24.0
26. 9. 22. 오전 10:30:00	24.1
26. 9. 22. 오전 11:00:00	24.3
26. 9. 22. 오전 11:30:00	24.5
26. 9. 22. 오후 12:00:00	24.7
26. 9. 22. 오후 12:30:00	25.0
26. 9. 22. 오후 1:00:00	25.1
26. 9. 22. 오후 1:30:00	25.1
26. 9. 22. 오후 2:00:00	25.3
26. 9. 22. 오후 2:30:00	25.2
26. 9. 22. 오후 3:00:00	25.2
26. 9. 22. 오후 3:30:00	25.4
26. 9. 22. 오후 4:00:00	25.2
26. 9. 22. 오후 4:30:00	25.0
26. 9. 22. 오후 5:00:00	24.8
26. 9. 22. 오후 5:30:00	24.6
26. 9. 22. 오후 6:00:00	24.5
26. 9. 22. 오후 6:30:00	24.5
26. 9. 22. 오후 7:00:00	24.2
26. 9. 22. 오후 7:30:00	24.1
26. 9. 22. 오후 8:00:00	24.0
26. 9. 22. 오후 8:30:00	24.0
26. 9. 22. 오후 9:00:00	24.0



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 22. 오후 9:30:00	24.0
26. 9. 22. 오후 10:00:00	23.9
26. 9. 22. 오후 10:30:00	23.9
26. 9. 22. 오후 11:00:00	23.8
26. 9. 22. 오후 11:30:00	23.7
26. 9. 23. 오전 12:00:00	23.8
26. 9. 23. 오전 12:30:00	23.8
26. 9. 23. 오전 1:00:00	23.7
26. 9. 23. 오전 1:30:00	23.6
26. 9. 23. 오전 2:00:00	23.7
26. 9. 23. 오전 2:30:00	23.6
26. 9. 23. 오전 3:00:00	23.7
26. 9. 23. 오전 3:30:00	23.6
26. 9. 23. 오전 4:00:00	23.5
26. 9. 23. 오전 4:30:00	23.6
26. 9. 23. 오전 5:00:00	23.5
26. 9. 23. 오전 5:30:00	23.5
26. 9. 23. 오전 6:00:00	23.5
26. 9. 23. 오전 6:30:00	23.6
26. 9. 23. 오전 7:00:00	23.5
26. 9. 23. 오전 7:30:00	23.6
26. 9. 23. 오전 8:00:00	23.5
26. 9. 23. 오전 8:30:00	23.6
26. 9. 23. 오전 9:00:00	23.8
26. 9. 23. 오전 9:30:00	24.0
26. 9. 23. 오전 10:00:00	24.0
26. 9. 23. 오전 10:30:00	24.2
26. 9. 23. 오전 11:00:00	24.5
26. 9. 23. 오전 11:30:00	24.8
26. 9. 23. 오후 12:00:00	24.9
26. 9. 23. 오후 12:30:00	25.1



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 23. 오후 1:00:00	25.1
26. 9. 23. 오후 1:30:00	25.0
26. 9. 23. 오후 2:00:00	25.2
26. 9. 23. 오후 2:30:00	25.2
26. 9. 23. 오후 3:00:00	25.2
26. 9. 23. 오후 3:30:00	25.2
26. 9. 23. 오후 4:00:00	24.9
26. 9. 23. 오후 4:30:00	24.6
26. 9. 23. 오후 5:00:00	24.3
26. 9. 23. 오후 5:30:00	24.2
26. 9. 23. 오후 6:00:00	24.1
26. 9. 23. 오후 6:30:00	24.0
26. 9. 23. 오후 7:00:00	24.0
26. 9. 23. 오후 7:30:00	24.0
26. 9. 23. 오후 8:00:00	24.0
26. 9. 23. 오후 8:30:00	23.9
26. 9. 23. 오후 9:00:00	23.8
26. 9. 23. 오후 9:30:00	23.8
26. 9. 23. 오후 10:00:00	23.8
26. 9. 23. 오후 10:30:00	23.9
26. 9. 23. 오후 11:00:00	23.8
26. 9. 23. 오후 11:30:00	23.7
26. 9. 24. 오전 12:00:00	23.7
26. 9. 24. 오전 12:30:00	23.7
26. 9. 24. 오전 1:00:00	23.6
26. 9. 24. 오전 1:30:00	23.6
26. 9. 24. 오전 2:00:00	23.5
26. 9. 24. 오전 2:30:00	23.7
26. 9. 24. 오전 3:00:00	23.5
26. 9. 24. 오전 3:30:00	23.5
26. 9. 24. 오전 4:00:00	23.4



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 24. 오전 4:30:00	23.5
26. 9. 24. 오전 5:00:00	23.5
26. 9. 24. 오전 5:30:00	23.4
26. 9. 24. 오전 6:00:00	23.5
26. 9. 24. 오전 6:30:00	23.4
26. 9. 24. 오전 7:00:00	23.4
26. 9. 24. 오전 7:30:00	23.5
26. 9. 24. 오전 8:00:00	23.4
26. 9. 24. 오전 8:30:00	23.5
26. 9. 24. 오전 9:00:00	23.5
26. 9. 24. 오전 9:30:00	23.6
26. 9. 24. 오전 10:00:00	23.7
26. 9. 24. 오전 10:30:00	23.8
26. 9. 24. 오전 11:00:00	23.9
26. 9. 24. 오전 11:30:00	24.0
26. 9. 24. 오후 12:00:00	23.9
26. 9. 24. 오후 12:30:00	24.0
26. 9. 24. 오후 1:00:00	24.0
26. 9. 24. 오후 1:30:00	24.2
26. 9. 24. 오후 2:00:00	24.1
26. 9. 24. 오후 2:30:00	23.9
26. 9. 24. 오후 3:00:00	24.0
26. 9. 24. 오후 3:30:00	23.9
26. 9. 24. 오후 4:00:00	23.9
26. 9. 24. 오후 4:30:00	23.8
26. 9. 24. 오후 5:00:00	23.7
26. 9. 24. 오후 5:30:00	23.7
26. 9. 24. 오후 6:00:00	23.7
26. 9. 24. 오후 6:30:00	23.6
26. 9. 24. 오후 7:00:00	23.5
26. 9. 24. 오후 7:30:00	23.6



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 24. 오후 8:00:00	23.5
26. 9. 24. 오후 8:30:00	23.5
26. 9. 24. 오후 9:00:00	23.4
26. 9. 24. 오후 9:30:00	23.5
26. 9. 24. 오후 10:00:00	23.4
26. 9. 24. 오후 10:30:00	23.4
26. 9. 24. 오후 11:00:00	23.4
26. 9. 24. 오후 11:30:00	23.3
26. 9. 25. 오전 12:00:00	23.3
26. 9. 25. 오전 12:30:00	23.3
26. 9. 25. 오전 1:00:00	23.3
26. 9. 25. 오전 1:30:00	23.3
26. 9. 25. 오전 2:00:00	23.3
26. 9. 25. 오전 2:30:00	23.3
26. 9. 25. 오전 3:00:00	23.2
26. 9. 25. 오전 3:30:00	23.2
26. 9. 25. 오전 4:00:00	23.2
26. 9. 25. 오전 4:30:00	23.4
26. 9. 25. 오전 5:00:00	23.5
26. 9. 25. 오전 5:30:00	23.5
26. 9. 25. 오전 6:00:00	23.6
26. 9. 25. 오전 6:30:00	23.6
26. 9. 25. 오전 7:00:00	23.6
26. 9. 25. 오전 7:30:00	23.7
26. 9. 25. 오전 8:00:00	23.8
26. 9. 25. 오전 8:30:00	23.8
26. 9. 25. 오전 9:00:00	23.9
26. 9. 25. 오전 9:30:00	23.9
26. 9. 25. 오전 10:00:00	23.7
26. 9. 25. 오전 10:30:00	23.6
26. 9. 25. 오전 11:00:00	23.5



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 25. 오전 11:30:00	23.4
26. 9. 25. 오후 12:00:00	23.4
26. 9. 25. 오후 12:30:00	23.4
26. 9. 25. 오후 1:00:00	23.3
26. 9. 25. 오후 1:30:00	23.3
26. 9. 25. 오후 2:00:00	23.2
26. 9. 25. 오후 2:30:00	23.3
26. 9. 25. 오후 3:00:00	23.3
26. 9. 25. 오후 3:30:00	23.2
26. 9. 25. 오후 4:00:00	23.3
26. 9. 25. 오후 4:30:00	23.2
26. 9. 25. 오후 5:00:00	23.2
26. 9. 25. 오후 5:30:00	23.2
26. 9. 25. 오후 6:00:00	23.2
26. 9. 25. 오후 6:30:00	23.2
26. 9. 25. 오후 7:00:00	23.1
26. 9. 25. 오후 7:30:00	23.1
26. 9. 25. 오후 8:00:00	23.1
26. 9. 25. 오후 8:30:00	23.1
26. 9. 25. 오후 9:00:00	23.0
26. 9. 25. 오후 9:30:00	23.0
26. 9. 25. 오후 10:00:00	23.1
26. 9. 25. 오후 10:30:00	23.2
26. 9. 25. 오후 11:00:00	23.3
26. 9. 25. 오후 11:30:00	23.4
26. 9. 26. 오전 12:00:00	23.4
26. 9. 26. 오전 12:30:00	23.5
26. 9. 26. 오전 1:00:00	23.6
26. 9. 26. 오전 1:30:00	23.7
26. 9. 26. 오전 2:00:00	23.8
26. 9. 26. 오전 2:30:00	23.7



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net
생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 26. 오전 3:00:00	23.8
26. 9. 26. 오전 3:30:00	23.8
26. 9. 26. 오전 4:00:00	23.8
26. 9. 26. 오전 4:30:00	23.7
26. 9. 26. 오전 5:00:00	23.7
26. 9. 26. 오전 5:30:00	23.9
26. 9. 26. 오전 6:00:00	23.9
26. 9. 26. 오전 6:30:00	23.9
26. 9. 26. 오전 7:00:00	23.9
26. 9. 26. 오전 7:30:00	23.9
26. 9. 26. 오전 8:00:00	23.9
26. 9. 26. 오전 8:30:00	24.0
26. 9. 26. 오전 9:00:00	24.1
26. 9. 26. 오전 9:30:00	24.2
26. 9. 26. 오전 10:00:00	24.3
26. 9. 26. 오전 10:30:00	24.5
26. 9. 26. 오전 11:00:00	24.6
26. 9. 26. 오전 11:30:00	24.7
26. 9. 26. 오후 12:00:00	24.8
26. 9. 26. 오후 12:30:00	24.9
26. 9. 26. 오후 1:00:00	24.9
26. 9. 26. 오후 1:30:00	24.9
26. 9. 26. 오후 2:00:00	24.8
26. 9. 26. 오후 2:30:00	24.7
26. 9. 26. 오후 3:00:00	24.6
26. 9. 26. 오후 3:30:00	24.6
26. 9. 26. 오후 4:00:00	24.5
26. 9. 26. 오후 4:30:00	24.4
26. 9. 26. 오후 5:00:00	24.3
26. 9. 26. 오후 5:30:00	24.1
26. 9. 26. 오후 6:00:00	23.9



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 26. 오후 6:30:00	23.9
26. 9. 26. 오후 7:00:00	23.8
26. 9. 26. 오후 7:30:00	23.6
26. 9. 26. 오후 8:00:00	23.6
26. 9. 26. 오후 8:30:00	23.5
26. 9. 26. 오후 9:00:00	23.5
26. 9. 26. 오후 9:30:00	23.6
26. 9. 26. 오후 10:00:00	23.9
26. 9. 26. 오후 10:30:00	24.1
26. 9. 26. 오후 11:00:00	24.1
26. 9. 26. 오후 11:30:00	24.1
26. 9. 27. 오전 12:00:00	24.2
26. 9. 27. 오전 12:30:00	24.3
26. 9. 27. 오전 1:00:00	24.2
26. 9. 27. 오전 1:30:00	24.3
26. 9. 27. 오전 2:00:00	24.2
26. 9. 27. 오전 2:30:00	24.2
26. 9. 27. 오전 3:00:00	24.4
26. 9. 27. 오전 3:30:00	24.3
26. 9. 27. 오전 4:00:00	24.3
26. 9. 27. 오전 4:30:00	24.3
26. 9. 27. 오전 5:00:00	24.3
26. 9. 27. 오전 5:30:00	24.3
26. 9. 27. 오전 6:00:00	24.3
26. 9. 27. 오전 6:30:00	24.3
26. 9. 27. 오전 7:00:00	24.3
26. 9. 27. 오전 7:30:00	24.3
26. 9. 27. 오전 8:00:00	24.4
26. 9. 27. 오전 8:30:00	24.5
26. 9. 27. 오전 9:00:00	24.5
26. 9. 27. 오전 9:30:00	24.7



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 27. 오전 10:00:00	24.9
26. 9. 27. 오전 10:30:00	25.0
26. 9. 27. 오전 11:00:00	25.0
26. 9. 27. 오전 11:30:00	24.9
26. 9. 27. 오후 12:00:00	25.0
26. 9. 27. 오후 12:30:00	25.0
26. 9. 27. 오후 1:00:00	25.0
26. 9. 27. 오후 1:30:00	25.1
26. 9. 27. 오후 2:00:00	25.1
26. 9. 27. 오후 2:30:00	25.1
26. 9. 27. 오후 3:00:00	25.0
26. 9. 27. 오후 3:30:00	24.9
26. 9. 27. 오후 4:00:00	24.9
26. 9. 27. 오후 4:30:00	24.7
26. 9. 27. 오후 5:00:00	24.6
26. 9. 27. 오후 5:30:00	24.4
26. 9. 27. 오후 6:00:00	24.3
26. 9. 27. 오후 6:30:00	24.2
26. 9. 27. 오후 7:00:00	24.1
26. 9. 27. 오후 7:30:00	24.0
26. 9. 27. 오후 8:00:00	23.9
26. 9. 27. 오후 8:30:00	23.9
26. 9. 27. 오후 9:00:00	23.8
26. 9. 27. 오후 9:30:00	23.8
26. 9. 27. 오후 10:00:00	23.8
26. 9. 27. 오후 10:30:00	23.7
26. 9. 27. 오후 11:00:00	23.6
26. 9. 27. 오후 11:30:00	23.6
26. 9. 28. 오전 12:00:00	23.8
26. 9. 28. 오전 12:30:00	24.0
26. 9. 28. 오전 1:00:00	24.2



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 28. 오전 1:30:00	24.3
26. 9. 28. 오전 2:00:00	24.4
26. 9. 28. 오전 2:30:00	24.4
26. 9. 28. 오전 3:00:00	24.4
26. 9. 28. 오전 3:30:00	24.5
26. 9. 28. 오전 4:00:00	24.5
26. 9. 28. 오전 4:30:00	24.5
26. 9. 28. 오전 5:00:00	24.5
26. 9. 28. 오전 5:30:00	24.5
26. 9. 28. 오전 6:00:00	24.5
26. 9. 28. 오전 6:30:00	24.4
26. 9. 28. 오전 7:00:00	24.2
26. 9. 28. 오전 7:30:00	23.9
26. 9. 28. 오전 8:00:00	23.8
26. 9. 28. 오전 8:30:00	23.8
26. 9. 28. 오전 9:00:00	24.0
26. 9. 28. 오전 9:30:00	24.0
26. 9. 28. 오전 10:00:00	24.1
26. 9. 28. 오전 10:30:00	24.2
26. 9. 28. 오전 11:00:00	24.4
26. 9. 28. 오전 11:30:00	24.5
26. 9. 28. 오후 12:00:00	24.6
26. 9. 28. 오후 12:30:00	24.7
26. 9. 28. 오후 1:00:00	24.8
26. 9. 28. 오후 1:30:00	24.8
26. 9. 28. 오후 2:00:00	24.9
26. 9. 28. 오후 2:30:00	24.9
26. 9. 28. 오후 3:00:00	24.8
26. 9. 28. 오후 3:30:00	24.7
26. 9. 28. 오후 4:00:00	24.6
26. 9. 28. 오후 4:30:00	24.5



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 28. 오후 5:00:00	24.3
26. 9. 28. 오후 5:30:00	24.2
26. 9. 28. 오후 6:00:00	24.0
26. 9. 28. 오후 6:30:00	23.9
26. 9. 28. 오후 7:00:00	23.7
26. 9. 28. 오후 7:30:00	23.6
26. 9. 28. 오후 8:00:00	23.6
26. 9. 28. 오후 8:30:00	23.6
26. 9. 28. 오후 9:00:00	23.5
26. 9. 28. 오후 9:30:00	23.4
26. 9. 28. 오후 10:00:00	23.4
26. 9. 28. 오후 10:30:00	23.3
26. 9. 28. 오후 11:00:00	23.4
26. 9. 28. 오후 11:30:00	23.2
26. 9. 29. 오전 12:00:00	23.2
26. 9. 29. 오전 12:30:00	23.3
26. 9. 29. 오전 1:00:00	23.1
26. 9. 29. 오전 1:30:00	23.1
26. 9. 29. 오전 2:00:00	23.2
26. 9. 29. 오전 2:30:00	23.2
26. 9. 29. 오전 3:00:00	23.1
26. 9. 29. 오전 3:30:00	23.1
26. 9. 29. 오전 4:00:00	23.1
26. 9. 29. 오전 4:30:00	23.1
26. 9. 29. 오전 5:00:00	23.1
26. 9. 29. 오전 5:30:00	23.2
26. 9. 29. 오전 6:00:00	23.3
26. 9. 29. 오전 6:30:00	23.3
26. 9. 29. 오전 7:00:00	23.4
26. 9. 29. 오전 7:30:00	23.5
26. 9. 29. 오전 8:00:00	23.6



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 29. 오전 8:30:00	23.7
26. 9. 29. 오전 9:00:00	23.7
26. 9. 29. 오전 9:30:00	23.7
26. 9. 29. 오전 10:00:00	23.7
26. 9. 29. 오전 10:30:00	23.8
26. 9. 29. 오전 11:00:00	24.0
26. 9. 29. 오전 11:30:00	24.2
26. 9. 29. 오후 12:00:00	24.3
26. 9. 29. 오후 12:30:00	24.5
26. 9. 29. 오후 1:00:00	24.5
26. 9. 29. 오후 1:30:00	24.6
26. 9. 29. 오후 2:00:00	24.6
26. 9. 29. 오후 2:30:00	24.7
26. 9. 29. 오후 3:00:00	24.6
26. 9. 29. 오후 3:30:00	24.6
26. 9. 29. 오후 4:00:00	24.5
26. 9. 29. 오후 4:30:00	24.3
26. 9. 29. 오후 5:00:00	24.2
26. 9. 29. 오후 5:30:00	24.1
26. 9. 29. 오후 6:00:00	23.9
26. 9. 29. 오후 6:30:00	23.8
26. 9. 29. 오후 7:00:00	23.7
26. 9. 29. 오후 7:30:00	23.6
26. 9. 29. 오후 8:00:00	23.5
26. 9. 29. 오후 8:30:00	23.3
26. 9. 29. 오후 9:00:00	23.3
26. 9. 29. 오후 9:30:00	23.5
26. 9. 29. 오후 10:00:00	23.6
26. 9. 29. 오후 10:30:00	23.7
26. 9. 29. 오후 11:00:00	23.7
26. 9. 29. 오후 11:30:00	23.7



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 30. 오전 12:00:00	23.7
26. 9. 30. 오전 12:30:00	23.7
26. 9. 30. 오전 1:00:00	23.7
26. 9. 30. 오전 1:30:00	23.7
26. 9. 30. 오전 2:00:00	23.7
26. 9. 30. 오전 2:30:00	23.7
26. 9. 30. 오전 3:00:00	23.7
26. 9. 30. 오전 3:30:00	23.7
26. 9. 30. 오전 4:00:00	23.7
26. 9. 30. 오전 4:30:00	23.6
26. 9. 30. 오전 5:00:00	23.6
26. 9. 30. 오전 5:30:00	23.6
26. 9. 30. 오전 6:00:00	23.7
26. 9. 30. 오전 6:30:00	23.7
26. 9. 30. 오전 7:00:00	23.7
26. 9. 30. 오전 7:30:00	23.7
26. 9. 30. 오전 8:00:00	23.7
26. 9. 30. 오전 8:30:00	23.7
26. 9. 30. 오전 9:00:00	23.7
26. 9. 30. 오전 9:30:00	23.8
26. 9. 30. 오전 10:00:00	23.8
26. 9. 30. 오전 10:30:00	24.0
26. 9. 30. 오전 11:00:00	24.1
26. 9. 30. 오전 11:30:00	24.2
26. 9. 30. 오후 12:00:00	24.4
26. 9. 30. 오후 12:30:00	24.5
26. 9. 30. 오후 1:00:00	24.5
26. 9. 30. 오후 1:30:00	24.6
26. 9. 30. 오후 2:00:00	24.6
26. 9. 30. 오후 2:30:00	24.7
26. 9. 30. 오후 3:00:00	24.6



IP Storage Temperature Log_Monthly report

2026년 9월 1일 오전 12시 0분 0초 KST - 2026년 10월 1일 오전 12시 0분 0초 KST

kyon410@daum.net

생성됨: 2026년 10월 1일 오전 6시 3분 22초 KST

측정 포인트	동국대학교 일산병원 임상시험센터 > 임상시험용의약품 보관실 > 임상시험용 의약품(실온보관장 1~30°C)_0054704396_CH2
타임 스탬프 [Asia/Seoul]	온도 [°C]
26. 9. 30. 오후 3:30:00	24.6
26. 9. 30. 오후 4:00:00	24.5
26. 9. 30. 오후 4:30:00	24.5
26. 9. 30. 오후 5:00:00	24.4
26. 9. 30. 오후 5:30:00	24.3
26. 9. 30. 오후 6:00:00	24.1
26. 9. 30. 오후 6:30:00	24.1
26. 9. 30. 오후 7:00:00	24.0
26. 9. 30. 오후 7:30:00	24.0
26. 9. 30. 오후 8:00:00	23.9
26. 9. 30. 오후 8:30:00	24.1
26. 9. 30. 오후 9:00:00	24.2
26. 9. 30. 오후 9:30:00	24.3
26. 9. 30. 오후 10:00:00	24.3
26. 9. 30. 오후 10:30:00	24.4
26. 9. 30. 오후 11:00:00	24.3
26. 9. 30. 오후 11:30:00	24.3