

Period of Record	Year. Month	2026.9	
Storage	Location	동국대학교일산병원 임상시험센터 약국	9134Room
	Description	전기냉장냉동기	항온고 1
	Manufacturer	SANYO	
	Model Name	MPR-720R	
	Serial Number	07120004	
	Temperature Range	15~25℃	
	Setting Point	19.0℃	
Thermometer	Description	Digital Thermometer	
	Manufacturer	SATO	
	Model Name	PC-3500	
	Serial Number	825289	
	Min/Max Record	Reset after recording	
Thermometer Calibration Certification	Date of Calibration	2026.1.27	
	Date of Calibration Due	2027.1.27	
	Certificate No.	SC2601-02474-4	
Temperature deviation alarm system	Manufacturer	Daeil Tech	
	Model Name	ACMS-DT01	
	Serial Number	DT15-ACMS-000348-K	
	Alarm threshold	Lower : 15℃, Upper : 25℃, alarm in 5 minutes	
	Alarm Method	Call to Pharmacists	
	Alarm Clear Method	000	

Date	Time	Temperature(℃)			Verified by	Comment
		Current	Min	Max		
2026.09.01	p:20	19	18	20	ml	
2026.09.02	p:10	19	18	20	ml	
2026.09.03	p:20	19	17	20	ml	
2026.09.04	p:20	19	18	20	ml	
2026.09.05						
2026.09.06						
2026.09.07	p:20	19	18	20	7658	
2026.09.08	p:20	19	17	20	ml	
2026.09.09	p:20	20	17	20	ml	
2026.09.10	p:10	18	17	20	ml	
2026.09.11	p:20	18	17	20	ml	
2026.09.12						
2026.09.13						
2026.09.14	p:20	19	17	20	ml	
2026.09.15	p:20	18	17	20	ml	
2026.09.16	p:20	18	17	20	ml	
2026.09.17	p:20	18	17	20	ml	
2026.09.18	p:20	18	17	20	ml	
2026.09.19						
2026.09.20						
2026.09.21	p:20	19	17	20	ml	
2026.09.22	p:15	20	17	20	ml	
2026.09.23	p:20	19	17	20	ml	
2026.09.24						추석연휴
2026.09.25						추석
2026.09.26						
2026.09.27						
2026.09.28	p:20	19	17	20	ml	
2026.09.29	p:10	18	17	20	ml	
2026.09.30	p:10	18	17	20	ml	
1900-01-00						